

**Parking brake**

For secure sitting.

Padded seat

can easily be folded back.

Aluminium frame

can easily be folded up.

EVA rear wheels

with drum brake,

Optional backwards stop

prevents inadvertent rolling back.

Tyres

special dimensions for outdoors.

I am your outdoor posterior walking aid.

Just like my brother, I am always with you when you really want to get moving, and will offer you the best-possible support to do so. Thanks to my outdoor outfit, I am happiest outdoors, and am the perfect companion for a walk: my large wheels will roll along any surface and my drum and parking brakes provide additional security.



Optional anti-tip mechanism protects you from tipping back.

Frame colours size 1**Frame colours size 2****Frame colours size 3**



Data	Size 1	Size 2	Size 3
Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm
Height of forearm rests	46 - 64 cm	63 - 87 cm	79 - 103 cm
Handle spacing	22 - 36 cm	30 - 44 cm	38 - 52 cm
Seat height	33 cm	45 cm	60 cm
Total width	56 cm	66 cm	73 cm
Total length	66 cm	75 cm	82 cm
Frame height (folded)	25 cm	32 cm	34 cm
Transport dimensions (L x W x H)	66 x 56 x 25 cm	75 x 66 x 32 cm	82 x 73 x 34 cm
max. load	30 kg	50 kg	70 kg
Weight	7,1 kg	8,3 kg	9,8 kg
HMV No.	10.46.02.3052		

Did you know that I am one of the few walking trainers which has an outdoor approval? With larger and wider wheels and a drum and parking brake, I am well-equipped for big adventures outdoors. Brilliant!



You can find more information on accessories in the price list.



Scan this QR code and watch videos about the product directly on our homepage.

