

GMFCS  
Levels III-V



### I am the car seat which grows with you.

I promise to stay loyal to you for many years. It doesn't matter how big you are now; I will adapt to fit your measurements – right up to a height of 1.50 m! My back height and seat width can be adjusted to suit you in one easy step. When the back height is adjusted, the straps will automatically adjust without sliding out of the strap guidance. So that you sit really securely, I have plenty of features for you: side pads stabilise your upper body, a positioning vest keeps you firmly in the saddle, and a torso-holding pad offers you a lot of security. But even safer than all this is the innovative side impact rotation if I am installed with the Isofix. In case of a side impact, I pivot away and thus protect you from an impact and catapulted parts.



15 bis  
36 kg



# siggi 3 fix.



**Side impact rotation (SIR)**  
If the car seat is fastened using the Isofix, it will pivot away from the impact point in case of a side impact.



**Fixation vest**  
provides you with even more security.





| Data                | siggi 3 fix.  |
|---------------------|---------------|
| Seat width at front | 32 cm         |
| Seat width at rear  | 27 cm         |
| Seat surface depth  | 25 cm         |
| Height of headrest  | 55 - 75 cm    |
| Total width         | 53 cm         |
| Total height        | 68 - 87 cm    |
| Back angle          | 90° / 110 °   |
| Group               | 15 - 36 kg    |
| Weight              | 8 kg          |
| HMV No.             | 26.11.06.0023 |

**Side guide and torso-holding pad**  
is placed on your lap as additional security.

**Optional foot support**  
is adjustable in height and angle.



You can find more information on accessories in the Price List.

