



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The universal toilet seat system.

- Usable as an attachment on the toilet, with subframe freestanding or as a shower chair
- Incl. pelvic harness and removable armrests
- Adjustable seat depth and back angle



HMV-No.: 33.40.01.4005
(use as a toilet attachment)

HMV-No.: 33.40.04.1008
(use as a toilet chair)

Now with
height-
adjustable
backrest for size 2

neu



Basic model

	hts. size 1	hts. size 2	hts. size 3
	Price / piece	Price / piece	Price / piece
	57 01000 1.134 €	57 02 100 1.134 €	57 03 100 1.238 €

Subframes

Mobile subframe hts. – incl. footplate and 75 mm castors (with locking device and directional lock), adjustable in height, with seat angle adjustment (-)5° - 25°



Mobile subframe hts. – incl. footplate and 75 mm castors (with locking device and directional lock), adjustable in height, with seat angle adjustment (-)15° - 15°



Fixed subframe hts.
incl. footplate, adjustable in height



Price / piece		Price / piece		Price / piece	
for size 1 (57 01025)	1.519 €	for size 1 (57 01002)	1.519 €	for size 1 (57 01026)	1.144 €
for size 2 (57 02025)	1.519 €	for size 2 (57 02002)	1.519 €	for size 2 (57 02026)	1.144 €
for size 3 (57 03025)	1.519 €	for size 3 (57 03002)	1.519 €	for size 3 (57 03026)	1.144 €

Untergestelle

Transfer and care subframe hts. – made of stainless steel, powder-coated, incl. footplate and 100-mm double castors, height-adjustable. With the transfer and care subframe **hts.**, we provide a further subframe for the Rifton **hts.** shower and toilet chair in order to better fulfil the transfer and care requirements of the user and their assistant. With the help of a pedal on both sides of the subframe the height can be adjusted. The seat can simply be tucked forward and backwards.

	for size 1+2	for size 3
Width	67 cm	71 cm
Length	84 cm	87 cm
Tilt-in-space	(-)5° - 25°	(-)5° - 25°
Seat height	51 - 78 cm	51 - 78 cm
Lower leg length	27 - 40 cm	40 - 53 cm
Price / piece		
for size 1+2 (57 01 030)	3.012 €	
for size 3 (57 03 030)	3.257 €	



Subframe

Bath subframe hts. – not possible for size 3, only for use without footplate

Width (bottom / top)	41 / 49 cm
Length	58 cm
Height (front)	25 cm
Price / piece	
for size 1+2 (57 01 003)	414 €



Accessories

Attachment bar – to mount on the toilet (standard configuration, if a subframe has not been selected)



for size 1+2 **(57 01 004)**
for size 3 **(57 03 004)**


Price / piece

197 €
197 €

Dimensions

			Size 1	Size 2	Size 3
A	Width	Mobile subframes	53 cm		57 cm
		Fixed subframes	51 cm		53 cm
		Transfer and care subframe	67 cm		71 cm
B	Seat height	Mobile subframes	39 - 55 cm		44 - 58 cm
		Fixed subframes	33 - 43 cm	36 - 44 cm	38 - 48 cm
		Transfer and care subframe	51 - 78 cm		
C	Seat width		20* / 25* / 30 cm	20* / 25* / 30* / 36 cm	28* / 33* / 38* / 43 cm
D	Seat depth		20 / 23 / 25 / 28 cm	29 / 32 / 34 / 37 cm	37 / 40 / 42 / 44 cm
E	Armrest height		13,5 cm	17 cm	20 cm
F	Back width		29 cm	34 cm	41 cm
G	Lower leg length toilet seat		22 - 34 cm	27 - 39 cm	34 - 44 cm
H	Back height		42 cm	40 / 45 / 50 cm	50 / 56 / 62 cm
I	Max. position above toilet height	Mobile subframes	52 cm		52,5 cm
		Fixed subframes	38 cm		42 cm
J	Length	Mobile subframes	74 cm		81 cm
		Fixed subframes	71 cm		74 cm
		Transfer and care sub-frame	84 cm		87 cm
K	Front edge of hygiene opening to back edge of subframe		48 cm		
	Adjustable height lateral pads size1		24 - 39,5 cm		33 - 52 cm
	Adjustable height lateral pads size2		25,5 - 41 cm		34,5 - 53,5 cm
	Back angle		(-)6°, 6°, 18°		
	Tilt-in-space		(-)5° - 25°		
			(-)15° - 15°		
	Max. load		34 kg	68 kg	113 kg

* adjustable with optional thigh guide



Accessories

Arm rests for size 1
detachable (standard configuration for size 1, if nothing else stated)



Height 13,5cm

Piece / pair

57 01 001 incl.

Arm rests for size 1-3
detachable (standard configuration for size 2, if nothing else stated)



Height 17cm

Piece / pair

57 02 001 incl.

Arm rests for size 2+3
detachable (standard configuration for size 3, if nothing else stated)



Height 20cm

Piece / pair

57 03 001 incl.



Pelvic support

new Available from April!

Price / piece

for size 2 + 3 (**57 02 037**) 198 €

Padding

Zusatzpolster

Back pad and seat pad
(closed at back)



Hole size in seat pad

for size 1 8x20cm

for size 2 9x25cm

for size 3 11x30cm

Piece / set

for size 1 (**57 01 013**) 394 €

for size 2 (**57 02 113**) 394 €

for size 3 (**57 03 113**) 394 €

Back pad and seat pad
(open at back)



Hole size in seat pad

for size 1 10x23cm

for size 2 11x30cm

for size 3 15x36cm

Piece / set

for size 1 (**57 01 012**) 394 €

for size 2 (**57 02 112**) 394 €

for size 3 (**57 03 112**) 394 €

Back pad with soft seat pad
(closed at back)



Hole size in seat pad

for size 1 11x22cm

for size 2 11x31cm

for size 3 12x37cm

Piece / set

for size 1 (**57 01 035**) 505 €

for size 2 (**57 02 135**) 505 €

for size 3 (**57 03 135**) 505 €

Back pad and seat pad
(completely closed)



Piece / set

for size 1 (**57 01 028**) 472 €

for size 2 (**57 02 128**) 472 €

for size 3 (**57 03 128**) 472 €

Seat pad
(completely closed)



Price / piece

for size 1 (**57 01 029**) 266 €

for size 2 (**57 02 029**) 266 €

for size 3 (**57 03 029**) 266 €

Accessories

Lateral pads – width and height adjustable

	Size 1	Size 2
Height	10cm	13cm
Depth	13cm	17cm
to suit	Size 1+2	Size 2+3

Chest width for size 1 15 – 28cm

Chest width for size 2 19 – 32cm

Chest width for size 3 28 – 42cm



Piece / pair

Size 1 (**57 01 007**) 351 €

Size 2 (**57 02 007**) 351 €

Size 1 with additional chest strap (**57 01 008**) 415 €

Size 2 with additional chest strap (**57 02 008**) 415 €

Lateral pads flip away – width and height adjustable

	Size 1	Size 2
Height	10cm	13cm
Depth	13cm	15cm
to suit	Size 1+2	Size 2+3

Chest width for size 1 12 – 30cm

Chest width for size 2 17 – 37cm

Chest width for size 3 27 – 46cm



Piece / pair

Size 1 (**57 01 032**) 488 €

Size 2 (**57 02 032**) 488 €

Size 1 with additional chest strap (**57 01 033**) 552 €

Size 2 with additional chest strap (**57 02 033**) 552 €

Accessories

Thigh guides
width adjustable



	Pad height	Pad width	Seat width range		
			Size1	Size2	Size3
57 00 010	7,5cm	18cm	20 / 25cm	25 / 30cm	33 / 38cm
57 01 010	7,5cm	28cm	-	25 / 30cm	33 / 38cm
57 02 010	7,5cm	28cm	-	20 / 25cm	28 / 33cm

Piece / pair

for size 1 – 3 (**57 00 010**) 189 €

for size 2 + 3 (**57 01 010**) 211 €

for size 2 + 3 extra thick thigh guide (**57 02 010**) 233 €

Tray / arm rest



Price / piece

for size 1 (**57 01 011**) 431 €

for size 2 (**57 02 011**) 431 €

for size 3 (**57 03 011**) 431 €

Shaped headrest – height, angle and depth adjustable



Price / piece

57 00 006 415 €

Accessories

Standard deflector – not usable when the seat pad is completely closed and in conjunction with an abduction wedge



Price / piece

for size 1+2 **(57 01014)**
for size 3 **(57 03014)**

83 €
83 €

Deflector with splash guard – not usable when the seat pad is completely closed



Price / piece

for size 1+2 **(57 01015)**
for size 3 **(57 03015)**

168 €
168 €

Large dynamic deflector and abductor
not usable when the seat pad is completely closed



Price / piece

for size 2+3 **(57 02034)**

394 €

Accessories

Abductor wedge



Price / piece

for size 1+2 **(57 01016)**
for size 3 **(57 03016)**

193 €
193 €

Toilet bucket



Price / piece

for size 1 **(57 01017)**
for size 2 **(57 02017)**
for size 3 **(57 03017)**

138 €
138 €
138 €

Bowl adapter – can only be used with open or partially open seat pad



Price / piece

for size 1+2 **(57 01036)**
for size 3 **(57 03036)**

109 €
109 €

Accessories

Footplate – when used as a toilet attachment



Price / piece

for size 1 **(57 01018)**
for size 2 **(57 02018)**
for size 3 **(57 03018)**

525 €
525 €
525 €

Calf guide – only for use with the footplate



Price / piece

for size 1+2 **(57 01019)**
for size 3 **(57 03019)**

109 €
109 €

Ankle cuffs – only for use with the footplate



Price / pair

for size 1+2 **(57 01020)**
for size 3 **(57 03020)**

124 €
124 €

Chest-shoulder strap



for size 1

for size 2+3

Price / piece

for size 1 **(57 01009)**
for size 2 **(57 02009)**
for size 3 **(57 03009)**

275 €
275 €
275 €

Accessories

Push handles – for mobile subframes



Price / pair

57 00030

237 €

Castors 100 mm – for mobile subframes



Upcharge / Set

57 00031

157 €

Travel base hts. – with bag, max. toilet height = 47 cm, not possible for size 3, only for use without footplate (only for temporary use when travelling and only in combination with the attachment bar)



Price / set

for size 1+2 **(57 01022)**

776 €