



Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The superior therapy chair.


- Incl. padding and cover
- Incl. back
- Incl. armrests
- Incl. 4 castors (100 mm) with directional locks at the front and brakes at the rear.
- Adjustment options:
 - Height, seat angle and backrest angle adjustment via electric motor with manual control
 - Seat depth and seat width
 - Height-, width- and depth-adjustable armrests
- Cover: Slate grey (No.124)
- Frame colour: Black grey

HMV-No.: 26.11.05.0051




Dimensions

		Size1	Size1b	Size 2	Size 3
A	Seat depth	30 - 42 cm	38 - 48 cm		45 - 55 cm
B	Seat width	28 - 38 cm		36 - 48 cm	
C	Short back height	42 - 60 cm		48 - 60 cm	
C	Long back height	52 - 70 cm		58 - 70 cm	
D	Seat height		40 - 55 cm		
E	Lower leg length	27 - 45 cm		35 - 52 cm	
F	Total width		60,5 cm		
G	Base frame length		84,5 cm		
	Seat angle		(-)20° - 20°		
	Back angle		(-)5° - 25°		
	max. load		125 kg		
	Weight		45 kg		



Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:

I Axilla height	_____	IV Pelvic width	_____
II Chest width	_____	V Seat depth	_____
III Back height	_____	VI Lower leg length	_____

	Price
62 00 080	incl.

Basic model

maxi. size 1	maxi. size 1b	maxi. size 2	maxi. size 3
Price		Price	Price
62 01 000	3.502 €	62 02 000	3.502 €
			62 03 000
			3.708 €

HD-Kit

maxi. HD-Kit – with this option, the maxi. is reinforced at key points (optional accessories are only affected if selected).

- Back plate made of steel instead of aluminium
- Backrest angle adjusted via a gas spring, using a release lever on the push handle, instead of an electric motor
- The foot angle of the flip-up footrest is adjusted using screws instead of pull pins



- Thoracic support brackets reinforced with an additional metal plate
- Armrests only available including the therapy table mounting – made of wood; armrest angle adjusted using a screw instead of a manni handle
- Push handle must be selected

Upcharge

for size 1+1b (62 01 047)	650 €
for size 2+3 (62 02 047)	650 €

Covers



Minor position

Slate grey (No. 124)		Soft grey (No. 127)		Flexi light grey (No. 50)		Without seat and back upholstery and cover – for the installation of custom cushions
Price		Price		Upcharge		Price reduction
91 00 124	incl.	91 00 127	incl.	62 00 050	272 €	62 00 001
						Netto -100 €

Headrests

Headrest (standard) incl. standard holder, padded, incl. cover, height and angle adjustable		Headrest (shell-shaped) incl. standard holder, padded, incl. cover, height and angle adjustable		Headrest (lateral guide) incl. standard holder, padded, incl. cover, height and angle adjustable	
Price / piece		Price / piece		Price / piece	
for size 1+1b (62 01 005)	335 €	for size 1+1b (62 01 006)	453 €	for size 1+1b (62 01 007)	391 €
for size 2+3 (62 02 005)	335 €	for size 2+3 (62 02 006)	453 €	for size 2+3 (62 02 007)	391 €

Headrests

Headrest (Lateral guide, angle adjustable) – incl. standard holder, padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable		Adapter for headrests for lateral adjustment		Double joint for headrest mounts to extend the depth adjustment range		Universal adapter for tailored headrests in the version 15 x 15 mm	
							
A = Width contact surface B = Height headrest C = Width lateral guide							
		Price / piece		Price / piece		Price / set	
Size 1		13,5 cm	10,5 cm	10 cm			
Size 2		15 cm	10,5 cm	13,5 cm			
Size 3		17 cm	12,5 cm	16,5 cm			
		Price / piece		Price / piece		Price / set	
Size 1 (62 01 008)		651 €		62 02 009	154 €	62 00 009	98 €
Size 2 (62 02 008)		651 €					
Size 3 (62 03 008)		651 €					

Accessories




Back – incl. pad and cover (short back is standard unless otherwise specified)

Back extension unit – incl. pad and cover (selectable for short back and for sizes 2 and 3)

Price / piece			Price / piece		
short for size 1+1b	(62 01 002)	incl.	for size 2+3	62 02 023	132 €
short for size 2+3	(62 02 002)	incl.			
long for size 1+1b	(62 01 003)	incl.			
long for size 2+3	(62 02 003)	incl.			

Zubehör

Lateral pads – padded, incl. cover, angle-, depth and width adjustabler



Lateral pads, long (trunk guidance)
padded, incl. cover, angle-, depth and width adjustabler



Lateral pads	Height	Depth
Size 1-3	12cm	16cm

Lateral pads	Height	Depth
Size 1-3	20cm	16cm

	Width 1	Width 2
Chest width for maxi. size 1+1b	17 - 25 cm	25 - 33 cm
Chest width for maxi. size 2+3	25 - 33 cm	33 - 41 cm
Price / pair		
Width 1 (62 01 022)	185 €	
Width 2 (62 02 022)	185 €	

	Width 1	Width 2
Chest width for maxi. size 1+1b	17 - 25 cm	25 - 33 cm
Chest width for maxi. size 2+3	25 - 33 cm	33 - 41 cm
Price / pair		
Width 1 (62 01 010)	185 €	
Width 2 (62 02 010)	185 €	

Accessories


Thigh guides – padded, incl. cover, adjustable in angle, depth and width



Removable thigh guide – padded, incl. cover, adjustable in angle, depth and width



Knee abduction wedge – pre-pluggable, padded, incl. cover




Price / pair			Price / pair			Price / piece		
for size 1+1b (Length 22cm)	(62 02 011)	218 €	for size 1+1b (Length 22cm)	(62 01 043)	320 €	62 00 014		249 €
for size 2+3 (Length 25cm)	(62 03 011)	218 €	for size 2+3 (Length 25cm)	(62 02 043)	320 €			

Accessories

Adjustable pelvic guidance – padded, incl. cover, angle adjustable

	Seat width A
for size 1	28 - 38cm
for size 2+3	36 - 48cm



Therapy tray – made of wood with rail, prefixable, fixed on both sides (can only be used with arm-rests incl. therapy table insertion)



	Height	Depth
Size 1	12cm	12cm
Size 2	15cm	16cm

	Front width	Rear width	Depth
for size 1-3	49,5cm	60cm	50cm
for size 2+3	59,5cm	62cm	50cm

Price / piece			Price / piece		
for size 1+1b	62 01 012	452 €	for size 1-3	(62 00 013)	294 €
for size 2+3	62 02 012	452 €	for size 2+3	(62 02 024)	294 €

Accessories

Armrests standard – width and height adjustable (not compatible with HD kit)



Price / pair

62 00 004 incl.

Accessories

Armrests incl. therapy table insertion – width, height and angle adjustable



Upcharge / pair

fixed **(62 01 004)** 247 €
fold-down **(62 02 004)** 371 €

Armrest padding, size 2 – for armrests incl. therapy table insertion.
(Standard configuration for the Basic model, size 1–2, unless otherwise specified.)

Length: 27 cm Width: 7 cm

62 02 046

Armrest padding, size 3 – for armrests incl. therapy table insertion.
(Standard configuration for the Basic model, size 3, unless otherwise specified.)

Length: 34 cm Width: 8 cm

62 03 046

Accessories

Plug-on footplate – height approx. 13 cm above the floor



Price / piece

62 00 015 402 €

Wooden footplate – height and knee angle adjustable, can be folded up, permanently attached to the seat unit and slots for attaching harnesses



Width

Depth

for size 1+1b

41 cm

26 cm


for size 2+3

45 cm

Price / piece

for size 1+1b **(62 01 016)** 525 €
for size 2+3 **(62 02 016)** 525 €

Calf plate for footplate – covered with black Flexi fabric (not possible with plug-on footplate)



Width

Height

for size 1+1b

32 cm

6,5 cm

for size 2+3

38 cm

Price / piece

for size 1+1b **(62 01 017)** 228 €
for size 2+3 **(62 02 017)** 228 €

Padding

Footrest padding for the attachable footrest

upholstered in black Flexi fabric, height: 1 cm (not compatible with the foot harness)

Price / piece

62 00 038 165 €

Footrest padding for height-adjustable footrest.

upholstered in black Flexi fabric, height: 1 cm (not compatible with the foot harness)

Price / piece

for size 1+1b **(62 01 038)** 133 €
for size 2+3 **(62 02 038)** 133 €

Accessories

Castors with locking devices – (standard with 2 locking devices at the rear, unless otherwise specified – with 4 locking devices, a directional lock is not possible)



Upcharge / set

100 mm castors 2 with locking devices **(62 00 018)** incl.
100 mm castors 4 with locking devices **(62 01 018)** 101 €
125 mm castors 2 with locking devices **(62 00 041)** 42 €
125 mm castors 4 with locking devices **(62 01 041)** 121 €

Central parking brake for rear castors

parking lever optionally mounted on the right or left



Upcharge / set

mounted left (100 mm castors) **(62 00 019)** 474 €
mounted right (100 mm castors) **(62 01 019)** 474 €
mounted left (100 mm front castors with locking device) **(62 02 019)** 567 €
mounted right (100 mm front castors with locking device) **(62 03 019)** 567 €
mounted left (125 mm castors) **(62 00 042)** 516 €
mounted right (125 mm castors) **(62 01 042)** 516 €
mounted left (125 mm front castors with locking device) **(62 02 042)** 609 €
mounted right (125 mm front castors with locking device) **(62 03 042)** 609 €

Accessories

Push bar – with rubber cover



Price / piece

for size 1+1b **(62 01 021)** 216 €
for size 2+3 **(62 02 021)** 216 €

Set for tool-free adjustment

Headrest, thoracic supports, back height, knee angle, seat depth, thigh guide, and lower leg length are adjustable without tools.

Price / set

62 00 045 78 €

Accessories

Harness assembly kit for the seat

consisting of 2 clip buckles and 2 clip buckles (drilled)

Harness assembly kit for the back

consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

Price / set

Price / set

25mm **(62 01 030)**

48 €


25mm **(62 01 031)**

48 €

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 1	Size 2	Size 3
A	13 cm	18 cm	18 cm
B	19 cm	20 cm	21 cm
C	15 cm	19 cm	28 cm
D	5 cm	5 cm	5 cm
E	28 cm	34 cm	44 cm
Strap width	2 cm	2,5 cm	2,5 cm



Price / piece

Size 1 **(62 01 032)**

197 €

Size 2 **(62 02 032)**


197 €

Size 3 **(62 03 032)**

197 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 1	Size 2
A	11 cm	18 cm
B	21 cm	22 cm
C	12 cm	15 cm
D	3,5 cm	3,5 cm
E	28 cm	34 cm
Strap width	2 cm	2,5 cm



Price / piece

Size 1 **(62 01 033)**


228 €

Size 2 **(62 02 033)**

228 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece

Size 1 **(62 01 034)**

99 €

Size 2 **(62 02 034)**

99 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece

Size 1 **(62 01 035)**

172 €

Size 2 **(62 02 035)**

172 €

Harness

Pelvic harness – by „Pelvi.Loc“, colour: black



	Size 2 (S)	Size 3 (M)	Size 4 (L)
for pelvic width	25 – 34 cm	32 – 40 cm	38 – 46 cm
Padding width	7 cm	7,5 cm	8,5 cm
Padding Length per side	19 cm	22 cm	24 cm
Strap Length per side	25 cm	45 cm	45 cm
Strap width	2,5 cm	2,5 cm	2,5 cm

Price / piece

Size 2 (S) **(62 02 039)**

325 €


Size 3 (M) **(62 03 039)**

325 €

Size 4 (L) **(62 04 039)**

325 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black



	Size 2	Size 3	Size 4
A	29 cm	33 cm	36 cm
B	17 cm	22 cm	26 cm
C	43 cm	49 cm	55 cm
Strap width	2,5 cm	2,5 cm	2,5 cm

Size 2 **(62 02 040)**

165 €

Size 3 **(62 03 040)**


165 €

Size 4 **(62 04 040)**

165 €


Harness

Ankle cuffs – with velcro fastener and idler pulley, colour: black (not possible in combination with padded cover for footplate/split footplate)



	Ankle circumference	Strap width
Size 1	12 – 19 cm	2,5 cm
Size 2	15 – 23 cm	2,5 cm
Size 3	19 – 27 cm	4 cm
		Price / pair
Size 1 (62 01 036)		236 €
Size 2 (62 02 036)		236 €
Size 3 (62 03 036)		236 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



	A	Strap width
Size 1	6,5 cm	2,5 cm
Size 2	8 cm	2,5 cm
Size 3	9,5 cm	2,5 cm
		Price / pair
Size 1 (62 01 037)		217 €
Size 2 (62 02 037)		217 €
Size 3 (62 03 037)		217 €