



rifton stander. prone

117

Invoice address:

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Dispatch address:

See invoice address

Commission:

The versatile standing trainer.

- Open frame with height adjustment and push handle
- Angle adjustment of the main column from 0° - 90° using a gas spring
- Swivelling trunk laterals for easy transfer
- Leg support with hip abduction and individual leg length adjustment
- Incl. trunk strap and pelvic harness for secure positioning
- Use in the supine position thanks to the mountable upper body support



Dimensions

	Size 1	Size 2	Size 3
Total length	81,5cm	99 cm	109 cm
Total width	53,5 cm	68,5cm	73,5cm
Horizontal body support height above floor	56 - 71cm	56 - 91,5 cm	60 - 96,5cm
Sandal to knee height	15 - 28cm	25,5 - 38 cm	31 - 48 cm
Knee to hip height	15 - 25,5cm	20,5 - 33 cm	30,5 - 44 cm
Sandal to top of headrest	78,5 - 109cm	104 - 139,5 cm	124,5 - 169 cm
Distance between trunk lateral	12,5 - 25,5cm	18 - 28 cm	23 - 35,5 cm
Distance between thigh medial and lateral support	5 - 10 cm	7,5 - 12,5 cm	10 - 14 cm
Basic item weight	21kg	32,5kg	52 kg
Body size	71 - 106cm	96 - 134 cm	122 - 165 cm
Max load	34 kg	45kg	84 kg



Extension for use in the supine position



Upper trunk support with a headrest – absolutely essential for extended use in the supine position!

Head laterals



	Price / piece		Price / pair
for size 1 (45 01 206)	435 €	for size 1+2 45 01 203	258 €
for size 2 (45 02 206)	447 €	for size 2+3 45 02 203	267 €
for size 3 (45 03 206)	459 €		

Extension for use in the supine position

Hip strap – Standard, if not stated otherwise



Hip strap with rotation control



One of these options is mandatory for extended use in the supine position!

	Price / piece		Price / piece
for size 1 (45 01 220)	383 €	for size 1 (45 01 221)	452 €
for size 2 (45 02 220)	395 €	for size 2 (45 02 221)	474 €
for size 3 (45 03 220)	405 €	for size 3 (45 03 221)	495 €

Extension for use in the supine position

Butterfly harness



	Price / piece
for size 1 (45 01 219)	452 €
for size 2 (45 02 219)	474 €
for size 3 (45 03 219)	495 €

Extension for use in the supine position

Knee supports



Swivelling knee supports
(only for size 2 available)



One of these options is mandatory for extended use in the supine position!

	Price / pair		Price / pair
for size 1 (45 01 205)	254 €	for size 2 (45 02 212)	532 €
for size 2 (45 02 205)	266 €	for size 3 (45 03 212)	562 €
for size 3 (45 03 205)	279 €		