



114 **rifton stander. supine**

Invoice address:

Dispatch address:

See invoice address

Order

Cost estimate

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The versatile standing trainer.

- Open frame with height adjustment and push handle
- Angle adjustment of the main column from 0° - 90° using a gas spring
- Swivelling trunk laterals for easy transfer
- Leg support with hip abduction and individual leg length adjustment
- Incl. hip strap and trunk strap for secure positioning
- Incl. soft knee straps with a cutout for the patella
- Use in prone position thanks to removable upper body support



Dimensions

	Size 1	Size 2	Size 3
Total length	81,5 cm	99 cm	109 cm
Total width	53,5 cm	68,5 cm	73,5 cm
Horizontal body support height above floor	56 - 71 cm	56 - 91,5 cm	60 - 96,5 cm
Sandal to knee height	15 - 28 cm	25,5 - 38 cm	31 - 48 cm
Knee to hip height	15 - 25,5 cm	20,5 - 33 cm	30,5 - 44 cm
Sandal to top of headrest	78,5 - 109 cm	104 - 139,5 cm	124,5 - 169 cm
Distance between trunk lateral	12,5 - 25,5 cm	18 - 28 cm	23 - 35,5 cm
Distance between thigh medial and lateral support	5 - 10 cm	7,5 - 12,5 cm	10 - 14 cm
Basic item weight	23,5 kg	35,5 kg	52 kg
Body size	71 - 106 cm	96 - 134 cm	122 - 165 cm
Max load	34 kg	45 kg	84 kg



Basic model

rifton stander. supine size 1 For users with a height of 71 – 106 cm Frame colour: red		rifton stander. supine size 2 For users with a height of 96 – 134 cm Frame colour: blue		rifton stander. supine size 3 For users with a height of 122 – 165 cm Frame colour: blue	
				 	
Price / piece		Price / piece		Price / piece	
45 01 100	4.800 €	45 02 100	5.400 €	45 03 100	6.200 €

Accessories

Hip strap – Standard, if not stated otherwise



Hip strap with rotation control



Price / piece initial equipment		Price / piece retrofitted		Upcharge / piece initial equipment		Price / piece retrofitted	
for size 1	45 01 120	incl.	45 01 120ET	383 €	for size 1	45 01 121	68 €
for size 2	45 02 120	incl.	45 02 120ET	395 €	for size 2	45 02 121	78 €
for size 3	45 03 120	incl.	45 03 120ET	405 €	for size 3	45 03 121	88 €

Accessories

Trunk strap – Standard, if not stated otherwise



Trunk strap with rotation control



Price / piece initial equipment		Price / piece retrofitted		Upcharge / piece initial equipment		Price / piece retrofitted	
for size 1	45 01 117	incl.	45 01 117ET	383 €	for size 1	45 01 118	68 €
for size 2	45 02 117	incl.	45 02 117ET	395 €	for size 2	45 02 118	78 €
for size 3	45 03 117	incl.	45 03 117ET	405 €	for size 3	45 03 118	88 €

Accessories

Butterfly harness



Upcharge / piece initial equipment		Price / piece retrofitted	
for size 1	45 01 119	68 €	45 01 119ET
for size 2	45 02 119	78 €	45 02 119ET
for size 3	45 03 119	88 €	45 03 119ET

Accessories

Knee supports

Standard, if not stated otherwise



Swivelling knee supports



Price / piece initial equipment		Price / piece retrofitted		Upcharge / piece initial equipment		Price / piece retrofitted	
for size 1	45 01 105	incl.	45 01 105ET	254 €	for size 2	45 02 112	273 €
for size 2	45 02 105	incl.	45 02 105ET	266 €	for size 3	45 03 112	293 €
for size 3	45 03 105	incl.	45 03 105ET	279 €			

Accessories

Head laterals	Additional trunk laterals (can be used as laterals at chest or hip level)	Medial thigh supports	Lateral thigh supports
			
Piece / pair	Piece / pair	Piece / pair	Piece / pair
for size 1+2 (45 01 103) 258 € for size 2+3 (45 02 103) 267 €	for size 1+2 (45 01 104) 258 € for size 2+3 (45 02 104) 267 €	for size 1+2 (45 01 108) 258 € for size 2+3 (45 02 108) 267 €	for size 1+2 (45 01 109) 258 € for size 2+3 (45 02 109) 267 €

Accessories

Legrest incl. abduction adjustment Standard, if not stated otherwise				Legrest incl. abduction adjustment, knee contracture compensation and hip & knee flexion					
									
Price / pair initial equipment		Price / pair retrofitted		Upcharge / pair initial equipment		Price / pair retrofitted			
for size 1	45 01 101	incl.	45 01 101ET	826 €	for size 2	45 02 102	387 €	45 02 102ET	1.233 €
for size 2	45 02 101	incl.	45 02 101ET	851 €	for size 3	45 03 102	397 €	45 03 102ET	1.257 €
for size 3	45 03 101	incl.	45 03 101ET	879 €					

Accessories

Sandal lift kit (only for size 1 available)	Tray – angle and height adjustable, removable	Elbow blocks for tray
		
Piece / pair	Price / piece	Piece / pair
for size 1 (45 01 111) 126 €	for size 1+2 (45 01 110) 711 € for size 3 (45 03 110) 733 €	for size 1+2 (45 01 113) 290 €

Extension for use in the prone position

Pelvic harness – absolutely essential for extended use in the prone position	Knee cuffs with rotation control	Upper body support
		
Price / piece	Piece / pair	Price / piece
for size 1 (45 01 122) 381 € for size 2 (45 02 122) 393 € for size 3 (45 03 122) 405 €	for size 1+2 (45 01 107) 258 € for size 3 (45 02 107) 279 €	for size 3 (45 03 222) 312 €