



**Order
Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Invoice address:

Dispatch address:

See invoice address

Commission:

The subframe with cradle system.

- Electric adjustment of the height, seat angle, backrest angle, knee angle and seat tilt by hand control (incl. locking function)
- Universal mounting points for all types of cradles
- Universal pick-up points for attachments
- Locking brake
- Depth-adjustable seat section
- Adjustable push bar in height and length
- Frame colour: white-blue




new

Dimensions - cradle subframe


		tina.
A	Height of cradle support from floor	53 - 79 cm
B	Adjustable back angle* (if C is set to 10°, then also 120° - 180°)	130° - 180°
C	Seat angle	0° - 40°
D	Adjustable knee angle	110° - 180°
	Base frame width	64 cm
	Max. load	100 kg
	Weight	70 kg

*if **C** is set to 10°, then **B** can also be set from 120° - 180°



Dimensions - cradle system

		Size 1	Size 2	Size 3
E	Width	59 cm (with fall out protection + 12 cm)		
F	Length of head section	60 cm	75 cm	90 cm
G	Length seat section	35 - 50 cm		
H	Length leg section	45 cm	55 cm	65 cm
I	Total length	140 - 155 cm	165 - 180 cm	190 - 205 cm
J	Height	6 cm (with fall out protection + 9 cm)		



Basic model cradle subframe

tina.

with **16" PU rear wheels** with drum brake & **8" PU steering wheels**



tina.

with **4 swivel castors** (125 mm) incl. rear directional adjustment & total lock (only for indoor use)



Price / piece

40 00 200

6.277 €


Price / piece

40 00 300

6.488 €

Cradle system

Cradle system – The system consists of a head section, seat section and leg section



Price / set

Size 1 (**40 01 201**)

791 €

Size 2 (**40 02 201**)

791 €

Size 3 (**40 03 201**) **Note:** This size includes an extended wheelbase!

1.240 €

Mattress

Vacuum mattress – for **tina.** incl. valve and electric pump with two-way switch (The vacuum mattress must be used with the reversible cover provided or with a suitable cover of your own.)



Support mattress (contourable with lateral guide) ncl. washable cover made of flexi-fabric & positioning cushion (set of 2) (The support mattress must be used with the reversible cover provided or with a suitable cover of your own.)

Filling made of polystyrene beads and soft foam



	Length	Width		Length	Width
for size 1	ca. 150 cm	ca. 75 cm	for size 1	ca. 150 cm	ca. 75 cm
for size 2	ca. 175 cm	ca. 75 cm	for size 2	ca. 175 cm	ca. 75 cm
for size 3	ca. 200 cm	ca. 80 cm	for size 3	ca. 200 cm	ca. 80 cm
Price / piece			Price / piece		
for size 1 (41 01 300)		3.588 €	for size 1 (41 01 302)		978 €
for size 2 (41 02 300)		3.588 €	for size 2 (41 02 302)		978 €
for size 3 (41 03 300)		3.588 €	for size 3 (41 03 302)		1.102 €

Mattress cover

Reversible cover for vacuum mattress. – made of knitted fabric in Flexi grey or deep sea blue



Price / piece

for size 1 (41 01 311)	420 €
for size 2 (41 02 311)	420 €
for size 3 (41 03 311)	420 €

Mattress cover

Reversible cover support mattress – made of knitted fabric in Flexi grey or deep sea blue



Price / piece

for size 1 (41 01 312)	420 €
for size 2 (41 02 312)	420 €
for size 3 (41 03 312)	420 €

Side protection

Side protection for size 1 – for head and seat section, fold-down


	Length head section	Length seat section	Length leg section	Height all
for size 1	20 cm	20 cm	–	15 cm
				Price / set
for size 1 (40 01 202)				343 €

Side protection for size 2 & 3 – for head, seat and leg sections, fold-down


	Length head section	Length seat section	Length leg section	Height all
for size 2	30 cm	20 cm	20 cm	15 cm
for size 3	50 cm		30 cm	
				Price / set
for size 2	(40 02 202)			512 €
for size 3	(40 03 202)			512 €

Accessories

Positioning cushion (set of 2)

(can only be ordered for the support mattress)



Width	Height
28 cm	22 cm

Price / set

4100304 195 €

Large support cushion


Width	Height
100 cm	39 cm

Price / piece

4100305 195 €

Half-moon cushion


Width	Height
90 cm	62 cm

Price / piece

4100306 119 €

T-cushion


Width	Height
60 cm	45 cm

Price / piece

4100307 119 €

Accessories

Wide seat belt – strap with sliding unit to fix on the tube frame


Price / piece

4100308 275 €

Preparing the base frame for fitting a push aid

Price / piece

for **Enjo** **(4000211)** 237 €
for **via Go** **(4001211)** 237 €

Accessories

Wheel base extension (10 cm) for size 1+2 (is already included in size 3)


for 16" rear wheels



for 125 mm castors

Price / pair

for 16" rear wheels **(4000206)** 195 €
for 125 mm rear castors **(4001206)** 195 €

Infusion holder

telescopic (is mounted on the right)




Price / piece

4000209 284 €

Accessories

Tray for medical devices on the push bar


Height and depth adjustable, with slotted holes for webbing (25 mm) to secure the medical devices



Width	56 cm
Depth	22 - 37,5 cm
Device height	max. 50 cm

Tray for medical devices on the base frame

with slotted holes for webbing (25 mm) for securing the medical devices



	tina. with PU wheels	tina. with swivel castors
Width	40 cm	
Depth	22 cm / 32 cm*	30 cm / 40 cm*
Device height	25 cm	
	*with wheel base extension	

Price / piece		Price / piece	
40 00 207	443 €	40 00 208	443 €

Accessories

Attachment loops for transporting the base frame in a car – mounted on the base frame, universally compatible, e.g., with systems from Q'Straint (approved according to DIN EN 7176-19; use in a 50° upright position of the headrest with a pelvic and diagonal shoulder belt or the vehicle's three-point seatbelt; if the required upright position cannot be achieved, a medical certificate must be carried in accordance with the Third Regulation on Exceptions to Road Traffic Regulations).



Here is the video for „installing the tina. into the vehicle.“



	Price / set
40 00 210	939 €