



**Invoice address:**

**Dispatch address:**

See invoice address

**Order**

**Cost estimate**

Customer No.: \_\_\_\_\_

Orderer: \_\_\_\_\_

Date: \_\_\_\_\_

Signature: \_\_\_\_\_

**Commission:**


**The dynamic standing trainer.**

- Open central column with height adjustment via a hand-tightened screw
- Supplied with a flexible pelvic pad and a protective pad to cover the manual adjustment screw
- Dynamic centre column:
  - Tilt angle up to 10° in all directions
  - Individual resistance adjustment (from grade 1-4) on each side
- Universal mount for accessories
- Can be used in anterior and posterior positions
- 4 castors (75 mm) with brakes
- Base plate made of birch with profiles to accommodate the heel edges
- Optional mounting point for assisted communication devices.



**HMV-No.: 28.29.01.1061**

**Dimensions**



|          |                                   | Size 0 (mini) | Size 1         | Size 2          |
|----------|-----------------------------------|---------------|----------------|-----------------|
|          | Body size                         | 65 cm - 95 cm | 85 cm - 115 cm | 105 cm - 140 cm |
| <b>A</b> | Sternum height                    | 42 cm - 62 cm | 58 cm - 86 cm  | 76 cm - 100 cm  |
| <b>B</b> | Table height from the base plate  | 48 cm - 72 cm | 63 cm - 100 cm | 80 cm - 115 cm  |
| <b>C</b> | Total length                      | 80 cm         | 88 cm          | 98 cm           |
| <b>D</b> | Total width                       | 50 cm         | 60 cm          | 70 cm           |
|          | Max. load                         | 20 kg         | 40 kg          | 60 kg           |
|          | Weight with dynamic centre column | 16 kg         | 19 kg          | 22 kg           |
|          | Weight with fixed centre column   | 13 kg         | 15 kg          | 18 kg           |

| Basic model „rigid“                                              |                                                                   |                                                                    |
|------------------------------------------------------------------|-------------------------------------------------------------------|--------------------------------------------------------------------|
| <b>trixi. size 0</b><br>for users with a height of<br>65 – 95 cm | <b>trixi. size 1</b><br>for users with a height of<br>85 – 115 cm | <b>trixi. size 2</b><br>for users with a height of<br>105 – 140 cm |
| Price / piece                                                    | Price / piece                                                     | Price / piece                                                      |
| Size 0 ( <b>63 00 000</b> ) 1.350 €                              | Size 1 ( <b>63 01 000</b> ) 1.450 €                               | Size 2 ( <b>63 02 000</b> ) 1.550 €                                |

| Basic model „dynamic“                                            |                                                                   |                                                                    |
|------------------------------------------------------------------|-------------------------------------------------------------------|--------------------------------------------------------------------|
| <b>trixi. size 0</b><br>for users with a height of<br>65 – 95 cm | <b>trixi. size 1</b><br>for users with a height of<br>85 – 115 cm | <b>trixi. size 2</b><br>for users with a height of<br>105 – 140 cm |
| Price / piece                                                    | Price / piece                                                     | Price / piece                                                      |
| Size 0 ( <b>63 00 100</b> ) 1.780 €                              | Size 1 ( <b>63 01 100</b> ) 1.880 €                               | Size 2 ( <b>63 02 100</b> ) 1.980 €                                |

| Pelvic pads                                                                                                                                                                                                                                       |                                                                                   |                                                                                     |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| <b>Pelvic pad (flex)</b> – adjustable in height and width for rotation control (standard, unless otherwise specified:<br>for basic model size 0 = strap size 0<br>for basic model size 1 = strap size 1<br>for basic model size 2 = strap size 2) | <b>Pelvic pad</b> – height and width adjustable with side guides                  | <b>Retaining strap for pelvic pad with side support</b>                             |
|                                                                                                                                                                  |  |  |
| Price / piece                                                                                                                                                                                                                                     | Price / piece                                                                     |                                                                                     |
| Size 0 ( <b>63 00 004</b> ) 298 €                                                                                                                                                                                                                 | Size 0 ( <b>63 00 003</b> ) 426                                                   | for size 0 ( <b>63 00 015</b> ) 172 €                                               |
| Size 1 ( <b>63 01 004</b> ) 298 €                                                                                                                                                                                                                 | Size 1 ( <b>63 01 003</b> ) 426                                                   | for size 1 ( <b>63 01 015</b> ) 172 €                                               |
| Size 2 ( <b>63 02 004</b> ) 298 €                                                                                                                                                                                                                 | Size 2 ( <b>63 02 003</b> ) 426                                                   | for size 2 ( <b>63 02 015</b> ) 172 €                                               |

| Lateral pads                                                                        |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| <b>Chest pad (flexible)</b> – adjustable in height and width for rotation control   | <b>Chest pad</b> – height and width adjustable with side guides                     | <b>Retaining strap for chest pad with side support</b>                                |
|  |  |  |
| Price / piece                                                                       | Price / piece                                                                       |                                                                                       |
| Size 0 ( <b>63 00 002</b> ) 298 €                                                   | Size 0 ( <b>63 00 001</b> ) 426 €                                                   | for sizes 0 ( <b>63 00 014</b> ) 172 €                                                |
| Size 1 ( <b>63 01 002</b> ) 298 €                                                   | Size 1 ( <b>63 01 001</b> ) 426 €                                                   | for sizes 1 ( <b>63 01 014</b> ) 172 €                                                |
| Size 2 ( <b>63 02 002</b> ) 298 €                                                   | Size 2 ( <b>63 02 001</b> ) 426 €                                                   | for sizes 2 ( <b>63 02 014</b> ) 172 €                                                |

| Covers for chest and pelvic pads                                                                                |                                                                                                                        |          |
|-----------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|----------|
| <b>Slate grey (No. 124)</b>  | <b>Flexi Light Grey (No. 50)</b>  |          |
| Price                                                                                                           |                                                                                                                        | Upcharge |
| <b>63 00 124</b> incl.                                                                                          | <b>63 00 050</b>                                                                                                       | 272 €    |

| Accessories                                                                                                                     |                                        |                                         |
|---------------------------------------------------------------------------------------------------------------------------------|----------------------------------------|-----------------------------------------|
| <b>Vest</b> – adjustable in size and height, including safety belt. Supports sensory perception. (only available in slate grey) | <b>Depth adjustment for chest pads</b> | <b>Depth adjustment for pelvic pads</b> |
|                                              |                                        |                                         |
| Price / piece                                                                                                                   | Price / piece                          | Price / piece                           |
| Size 0 ( <b>63 00 013</b> ) 480 €                                                                                               | up to 6 cm ( <b>63 00 005</b> ) 98 €   | up to 6 cm ( <b>63 00 017</b> ) 98 €    |
| Size 1 ( <b>63 01 013</b> ) 480 €                                                                                               | up to 10 cm ( <b>63 01 005</b> ) 98 €  | up to 10 cm ( <b>63 01 017</b> ) 98 €   |
| Size 2 ( <b>63 02 013</b> ) 480 €                                                                                               |                                        |                                         |

### Accessories

#### Flexible Knee Positioner (FKP) Standard



#### Heel edges incl. strap



Price / pair

Price / pair

Size 0 **(63 00 106)**

490 €

Size 0 **(63 00 012)**

225 €

Size 1 **(63 01 106)**

490 €

Size 1 **(63 01 012)**

225 €

Size 2 **(63 02 106)**

490 €

Size 2 **(63 02 012)**

225 €

### Accessories

**Tray** – height and depth adjustable,  
with edge


Price / piece

for sizes 0 **(63 00 008)**

425 €

for sizes 1+2 **(63 01 008)**

425 €

#### Armrests for tray



Price / piece

for size 0 **(63 00 010)**

215 €

for size 1+2 **(63 01 010)**

215 €

### Accessories

**Bracket for assisted communication**  
for mounting on the tray


Price / piece

**63 00 011**

110 €

**Bracket for assisted communication**  
for mounting on the centre column


Price / piece

**63 01 011**

110 €

**Dynamic foot cushion for sensory stimulation** – incl. cover and fastening  
material


Price / piece

for size 0 **(63 00 018)**

725 €

for size 1 **(63 01 018)**

725 €

for size 2 **(63 02 018)**

725 €

### Accessories

**Angle adjustment of the base plate**  
tilt adjustable in increments to  
2.5°, 5° or 7.5°


Price / piece

**63 00 019**

380 €

**stixi.** – The innovative app with a standing  
diary for documentation and child-friendly  
game.

The standing diary enables structured  
recording of standing times.

The integrated game offers playful motivation  
for children.

EN-Version coming soon in the Apple App  
Store and Google Play Store for smartphones  
and tablets.
