

GMFCS
Levels II–IV



Backrest
can be adjusted
from 90° to 180°.

Pelvic harness
provides reliable retention.

Seat surface
moulds to fit body
proportions.

I am your portable bath seat. I really want to go bathing with you. But that's not all – I can't wait to join you when you're building sandcastles too. Hair full of sand, hands covered in dirt... we definitely deserve a shower afterwards! My backrest can be adjusted so you can sit upright or lean back and relax. And by the way: I'm so light and compact that I can be folded up easily. That makes me perfect for trips away or visits to friends.



Back angle adjustment
Bring your back into the
desired position.





data	bibi.	bibi. with a long back
Seat depth	30 cm	
Seat width	36 cm	
Height of back	45 cm	55 cm
Total length	60 cm	
Back angle	90° - 180°	
Max. load	30 kg	
Weight	2.3 kg	
HMV-No.	04.40.04.1xxx	

Aluminium frame

I can be folded up to save space and with minimal effort.



You can find more information on accessories & harnesses in the price list.

