

GMFCS
Level III-V



I am the therapy chair for when you are young – I will give you security from the start.

I will place you on a level with your friends. My comfortable seat can be easily adjusted in height so that you can sit at any table and then be lowered back down to a height suitable for play. I will offer you just the right amount of support to help you grow and learn through your youngest years.





madita-fun mini.

Push handle
facilitates transfer.

Armrests
can be adapted in
height and width.

Pelvic harness
provides the necessary support.

Headrest
is available in a three-dimensionally
adjustable or height adjustable version.

Lateral pads
facilitate a stable and
upright sitting position.

Foot support
padded.

Seat height
can be adapted with
support from a gas pressure spring.

Castors
facilitate transfer.

**As individual as you are: in addition to a
standard seat pad, you can select a pad with
anatomic leg guidance or a dysplasia seat.**



Ideal for initial treatment.



data	mädita-fun mini.
Seat depth	12 - 19 cm
Seat width	14 - 22 cm
Height, short back	28 - 33 cm
Height, long back	30 - 35 cm
Seat height from floor	25 - 62 cm
Lower leg length	12 - 22 cm
Total width	52 cm
Length of base frame	58 cm
Seat angle	(-)5° - 37° / 0° - 42°
Back angle	(-)5° - 25°
Max. load	20 kg
Weight	16 kg
HMV-No.	26.11.05.0022



You can find more information on fabrics & embroidery
on pages 156 - 159.

You can find more information on accessories & harnesses
in the price list.

