

GMFCS
Level III-V

max. 50



(-) 5° – 25°



(-) 8° – 37°



15–42 cm



25–62 cm



16–38 cm

**I am your integrative therapy chair.**

Whether you're at the breakfast table at home or in the play corner at school, you can sit at the same level as your friends. I can be adjusted to suit your body and your taste – and I'm available in lots of great colours!

Did you know I have a unique adjustable pelvic support? It can be positioned exactly where you need it, giving you stable, comfortable support so you can sit upright and feel confident while you play, learn and join in.





Headrest

can be adjusted in height, depth and angle.

Lateral pads

allow a stable and upright sitting position.

Armrests

are individually adjustable.

Seat surface

is individually adaptable.

Push handle

for easy manoeuvrability.

Back unit

optionally removable, can be adjusted in angle and height.

Pelvic support**Foot support**

is available in three versions:

- without lock
- with lock
- with progressive angle adjustment.

Castors

available in different sizes, facilitate movement.

Base frame

compact and sturdily built.



In just one single step, my back unit can simply be removed. This gives you more freedom to move.

**Frame colours**

- Anthracite metallic
- Telegrey 4



My headrest with angle-adjustable side support is flexible, adaptable, practical and offers great therapeutic benefits.

	Size 0	Size 1	Size 1b	Size 2
Seat depth*	15 – 25 cm	21 – 30 cm	26 – 36 cm	30 – 42 cm
Seat width	16 – 26 cm	20 – 30 cm		28 – 38 cm
Back height with extra short back	32–35 cm		35–43 cm	
Back height with short back	35 – 43 cm		42 – 52 cm	
Back height with long back	42 – 52 cm		53 – 63 cm	
Back height with active backrest	35 – 43 cm			42 – 52 cm
Seat height	25 – 62 cm			
Seat height with electric height adjustment	32 – 62 cm			
Seat height with hydraulic height adjustment	27 – 63 cm			
Lower leg length	12 – 22 cm** / 15 – 25 cm	16 – 29 cm** / 19 – 32 cm	19 – 33 cm** / 22 – 36 cm	25 – 43 cm
Total width	54 cm		64 cm	
Total width with comb subframe	66 cm		76 cm	
Length of base frame	61 cm		75 cm	
Total length with combi subframe	76 cm		86 cm	
Tilt-in-space	(-)8° – 37°			
Back angle	(-)5° – 25°			
Max load.	30 kg		45 kg	50 kg
Weight	19 kg	20 kg	22 kg	24 kg
HMV-No.	26.11.05.0022			
* reduced by 1 cm with adjustable pelvic guidance ** this adjustment range can be achieved by repositioning the bracket				

* reduced by 1 cm with adjustable pelvic guidance | ** this adjustment range can be achieved by repositioning the bracket



You can find more information on fabrics & embroidery
on pages 156 - 159.

You can find more information on accessories & harnesses
in the price list.

