

GMFCS
Level III-V

kg
max. 16.5



61-105 cm

43 cm

61 cm

360°

34,5°-74,5°



ECE R129/03
ZERTIFIZIERT

i-size
universal
ISOFIX

I am your safe companion on the go.

And I am your best friend for relaxed journeys.

With my ergonomic seat, you can sit comfortably even on long trips. My unique EPP shell structure provides enhanced impact protection, and together with the side impact buddy, which absorbs up to 40% of impact energy, you are optimally protected. To make sure you can enjoy car journeys together for as long as possible, you can use me in both forward-facing and rear-facing positions. Thanks to my Isofix function, installing me on the back seat is child's play.

And your parents will love my 360-degree rotation – they can simply turn me to lift you in and out with ease.





The EPP shell structure
of the seat offers improved
impact protection.

Lateral pads
are adjustable in width.

Isofix mount
for use in
in both directions.

360° Rotation
for comfortable
installation and removal.

Headrest
can be adjusted in height.

With the 5-point harness
your child will remain safely in place.

Abduction wedge
enables optimum positioning.

Adjustable tilt angle
for 4 different sitting and reclining positions.

Foot rest
to relieve the legs,
the back and neck.

Height adjustment of the footrest
the car seat grows with the child.

**Rear facing
travel position.**

**Forward facing
travel position.**

**Comfortable and
easy to get into
and out of.**

**0 to
16.5 kg**

**9 to
16.5 kg**





data	marie integral.
Recommended body size	61–105 cm
Length	76 cm
Width	43 cm
Height	61 cm
Seat height	12.5 cm (lowest point to the horizontal plane in the seat)
Seat depth	28 cm
Seat width	28 cm
Arm rest height	34 cm
Back rest height	Max. 41 cm
Tilt-in-space	Rear-facing: 34,5°–49,5°; Forward-facing: 59,5° – 74,5° measured to horizontal plane.
Headrest height	19 cm
Footrest height	–2 cm to + 16 cm from the horizontal surface
Weight	15,6 kg
Max. load	16,5 kg
Fastening the safety device	5-point harness
HMV-No.	26.11.06.0035



Optional abduction wedge for better leg positioning.



You can find more information on accessories in the price list.

