

GMFCS
Level I-III

max. 50



29 cm



90° – 135°



26–35 cm

**I am your comfortable, foldable rehab buggy.**

I adapt to your needs with my height-adjustable and detachable foot support. My wheels cannot only be locked with a centrally operated brake for the rear wheels, but are also removable for transport. The front wheels can be steered, and are also provided with a directional lock. And the angle-adjustable backrest, the 5-point harness padded in the shoulders and crotch, the protective cover and the practical storage basket complete my convenient features.



Push handle

for easy handling.

Backrest

can be adjusted by 45° in total in five steps.

Aluminium frame

for easy transfer.

Storage basket

for the important things of everyday life.

Rear wheels

detachable and with centrally operated parking brake.

5-point strap

with padded shoulder and crotch harness facilitate optimum hold.

Calf support

can be adjusted in angle.

Foot support

can be adjusted in height.

Wheels

detachable, can be steered and with direction locks.



Super easy: With a kick on the brake lever both rear wheels are braked simultaneously.

Simply super

is also the variable backrest, which not only allows a comfortable seating position, but also protects against crosswind.





data	pipe.
Backrest height	68 cm
Seat depth	29 cm
Seat width	26 / 29 / 32 / 35 cm
Lower leg length	25 - 37,5 cm
Backrest angle	90° - 135°
Seat angle	26°
Total length	108 cm
Total width	59 cm
Total height	100 cm
Weight	50 kg
max. load	14,5 kg
Length (folded)	100 cm
Width (folded)	31 cm
Height (folded)	30 cm
HMV No.	18.99.01.1052

For easy transport I can easily be folded down.



You can find more information on accessories in the price list.

