

GMFCS
Levels I-V
kg
max. 140
38–91 cm**I'm the trike that's made just for you.**

I'm here to accompany you during your training. I help you become more mobile while also improving your balance, your spatial and visual perception, and the strength of your legs. As you'd expect from me, all my adjustments can be made without tools. My wide range of accessories means I can be customised to your individual needs. I offer four different headrests, as well as various steering handles and brakes for your companion. My backrest and pedal resistance can also be adjusted to suit you perfectly. Let's go!





rifton tricycle.

Headrests

is available in three different versions.

Handlebar

available as a straight version or as a curved handlebar for comfortable support.

Backrest

can be adjusted to 0°, 7 1/2° and 15° and can also be removed and folded for transport.

Wheels with quick-release fastener

can be removed in a flash thanks to the quick release and are also absolutely puncture-proof.

Pedal straps

can be released and moved to any position around the pedal.



As a special innovation, I'd like to introduce you to my stable stand. It transforms your tricycle into an exercise bike, so you can learn to cycle from the comfort of your own home.



Frame colour small



Turquoise (TEA)

Frame colour medium



Red (RED)

Frame colour large



Blue (BLU)



data		rifton tricycle. small	rifton tricycle. medium	Rifton tricycle. Large
Total height		54 cm		76 cm
Total length		112 cm	135 cm	165 cm
Total width		67 cm	70 cm	80 cm
Overall length without rear handle (wheels and front section removed)		91 cm	104 cm	124 cm
Overall length with rear handle (wheels and front section removed)		122 cm	130 cm	152 cm
Seat up to the top edge of the torso support (max. height)		51	61 cm	
Pedals to top edge of saddle		38 - 51 cm	48 - 71 cm	64 - 91 cm
Adjustment range of the side panels	Small torso support	15 - 30 cm		-
	Large torso support	-	24 - 39 cm	
Max. load		57 kg	113 kg	140 kg
Weight		22 kg	25 kg	29 kg