

GMFCS
Level IV-V

kg
max. 84



71–165 cm

0° – 90°



I am your versatile companion!

Because I combine many functions in one aid: with me, you get support from the front, from behind or just the way you need it. This makes me your universal companion in all situations. Meet your friends at eye level while you are standing and enjoy playing together – I'll hold you!





rifton stander. supine

Padding
with fire-retardant cover.

Armrest
for positioning the arms.

Main pillar
angle adjustable via gas pressure spring.

Knee braces
with patella recess.

8" tires
lockable.

Headrest
adjustable in depth.

Thoracic pads
can be folded down
for easy transfer.


Lap belt
provides the necessary safety.

Footrests
can be swiveled and tilted
and have a magnetic
quick-release fastener.



Variable adjustment options

Legs can be individually adjusted
in length, hip and knee angle.



Fold-down chest pads

offer upright, stable support and
enable easy transfer.



Frame colour
Size 1



Red

Frame colour
Size 2+3



Blue



data	Size 1	Size 2	Size 3 ^{new}
Total length	81,5 cm	99 cm	109 cm
Total width	53,5 cm	68,5 cm	73,5 cm
Height of the horizontal body support above the floor	56 - 71 cm	56 - 91,5 cm	60 - 96,5 cm
Height foot to knee	15 - 28 cm	25,5 - 38 cm	31 - 48 cm
Height knee to hip	15 - 25,5 cm	20,5 - 33 cm	30,5 - 44 cm
Height foot to upper edge of headrest	78,5 - 109 cm	104 - 139,5 cm	124,5 - 169 cm
Width between the fuselage side panels	12,5 - 25,5 cm	18 - 28 cm	23 - 35,5 cm
Width between medial/lateral thigh support	5 - 10 cm	7,5 - 12,5 cm	10 - 14 cm
Weight of the basic model	23,5 kg	35,5 kg	52 kg
Body height	71 - 106 cm	96 - 134 cm	122 - 165 cm
Max. load	34 kg	45 kg	84 kg

The abduction, which provides up to 30° for each leg, promotes healthy hip development.



More information about accessories & harnesses
can be found in the price list.

