

Positions you well.

GMFCS
Level IV-V



I am your adaptable reclining standing system.

With me, your transfer from lying to standing is safe and comfortable – even if you're not able to take an active part. Thanks to my wide range of adjustment options, I can be adapted precisely to your body shape from head to toe. Regular standing has a positive effect on your body's structures and functions, and it promotes your overall activity. I support you in this by giving you both safety and weight-bearing assistance at the same time. With me, you gain a new perspective – you're at eye level with your friends and can be moved around while standing. This means you can take part in different activities and situations throughout your standing training.


Height adjustment

is easily achieved via motor or with gas spring support.



Growth adjustments

can be easily adjusted by changing the length of the backrest and legs of the.











My height and angle can be easily adjusted. So nothing stands in the way of playing with your friends.



Precise tilting of the lying surface
is possible from 0 - 90° via
gas pressure spring.





Headrest with angle-adjustable side guides


is flexible, practical, individually adjustable and offers great therapeutic benefits.



Flexible knee positioners (FKP)


can be adjusted individually in three dimensions. Length-adjustable support straps adapt perfectly to the shape of the knee.

An optimal foot position
is individually 3D-adjustable. Footrests reached. These can be lowered to the floor.



Perfect leg guidance
is made possible by legs that are adjustable in length and angle.





Precise abduction of the legs
Can be adjusted from 0-30°
using the adjustment screw



data		Size 1	Size 2	Size 3	Size 4
Height		75 - 115 cm	100 - 135 cm	120 cm - 160 cm	145 - 185 cm
Back	Hip/thorax width	14 - 27 cm	17 - 31 cm	20 - 37 cm	27 - 46 cm
	Length of lying surface	30,5 - 43 cm	38 - 48 cm	49 - 61 cm	58 - 69 cm
Legs	Sole to knee	20,5 - 30,5 cm	27,5 - 37 cm	32,5 - 45 cm	42 - 54 cm
	Knee to crotch	12 - 18,5 cm	16,5 - 24 cm	20 - 29,5 cm	24,5 - 34 cm
	Knee angle	0° - 65°			
	Hip angle	0° - 45°			
	Abduction	0° - 30°			
	Level compensation	0 - 5 cm			
Feet	Equinus	0° - 20°			
	Flat foot	0° - 15°			
	Rotation	0° - 15° each direction			
Height of lying surface (for transfer)		53,5 - 73,5 cm	60 - 85 cm	71,5 - 93,5 cm	82,5 cm - 104,5 cm
cantling		0° - 90°			
Transport dimensions (L x W x H)		95 x 51 x 75 cm	105 x 61 x 85 cm	125 x 76 x 97 cm	152 x 77 x 110 cm
Max. load		30 kg	60 kg	80 kg	100 kg
Weight		30 kg	35 kg	80 kg	90 kg
HMV-No.		28.2902.1057			



Further information on accessories
can be found in the price list.

