

Always recline correctly.

GMFCS
Level IV-V



I am your modular cradle system.

I'll make sure you're always as comfortable as possible. For the most complex of cases my optional vacuum mattress relieves pressure points and pain, while my optional cushioning provides comfort for less demanding needs. I'm designed for both indoor and outdoor use, giving you the freedom to explore. With my electronic height adjustment and customizable seat, back, and knee angles, you'll always find the perfect position—whether sitting or lying down. And thanks to my smooth-rolling castor wheels, I'm easy to push and maneuver, so you can take me anywhere with ease.







You can recline on me, but also sit comfortably too so that you can always see the world around you.



Fall-out protection
for maximum safety.



Electric pump
for the vacuum mattress.



Data tina.	three-piece cover system
Height of cradle support from floor	53 - 79 cm
Adjustable back angle	130° - 180°
Adjustable seat angle	0° - 40°
Adjustable knee angle	110° - 180°
Wide base frame	64 cm
Max. load	100 kg
Weight	70 kg



Cradle system	Size 1	Size 2	Size 3
Width	59 cm (with fall protection + 12 cm)		
Length of head section	60 cm	75 cm	90 cm
Length seat section	35 - 50 cm		
Length leg section	45 cm	55 cm	65 cm
Total length	140 - 155 cm	165 - 180 cm	190 - 205 cm
Height	6 cm (with fall protection + 9 cm)		



You can find more information on accessories & harnesses
in the price list.

