

GMFCS
Level III-IV



I am your mobile standing trainer.



Using me, you can practice standing on your own. Regular standing trains your muscles and improves bone density. My body supports, which can be adapted in height, width and depth, can adapt to your body measurements. My two large wheels round me off perfectly. Using these, you can go anywhere – completely independently.



Flexible knee positioners (FKP)

are individually three-dimensionally adjustable. Length-adjustable contact straps adapt optimally to the shape of the knee.







data		Size 1		Size 2			Size 3		Size 4		
Body size		75 - 110cm		100-130 cm			120 - 155 cm		150 - 175 cm		
Wheel size		24"	28"	30"	32"	36"	32"	36"	36"	40"	42"
Grip height: footplate - hand rims at top		40 - 48 cm	50 - 58 cm	57 - 65 cm	64 - 69 cm	72 - 80 cm	60 - 68 cm	71 - 79 cm	82 - 90 cm	87 - 95 cm	
Recommended from elbow height		50 cm	60 cm	70 cm	80 - cm	92 cm	82 cm	92 cm		97 cm	99 cm
Width of hand rims at top	0°	-	-	-	-	-	65 cm		67 cm		
	3°	-	-	54 cm	53 cm	52 cm	61 cm	60 cm	63 cm	62 cm	61 cm
	6°	-	-	51 cm	49 cm	47 cm	57 cm	55 cm	59 cm	57 cm	56 cm
	9°	45 cm	44 cm	48 cm	45 cm	42 cm	54 cm	51 cm	55 cm	53 cm	51 cm
	12°	42 cm	41 cm	45 cm	41 cm	38 cm	-	-	-	-	-
	15°	40 cm	-	-	-		-	-	-	-	-
Width of hand rims at bottom (min. - max.)		68 - 74 cm	70 - 73 cm	71 - 83 cm		71 - 84 cm	73 - 83 cm	73 - 85 cm	76 - 86 cm		76 - 90 cm
Footplate height – centre of knee support pads		19 - 30 cm		28 - 36 cm			34 - 44 cm		42 - 51 cm		
Footplate height – buttocks support pad support		26 - 42 cm		43 - 58 cm			53 - 69 cm		65 - 80 cm		
Recommended for inside leg length		30 - 47 cm		48 - 63 cm			58 - 74 cm		70 - 84 cm		
Abduction split footplate		0° - 6°					0° - 30°				
Abduction split footplate (with abducted basic model)		0° - 40°					-		-		
Pelvis depth*		11 - 24 cm					15 - 28,5 cm				
Pelvic depth** (adjustment range of the pelvic position by hand crank max. 6 cm)		11 - 20 cm (+6 cm)					15 - 27 cm (+6 cm)				
Sternum height (Footplate height - centre of chest support pad)		54 - 80 cm		68 - 101 cm			92 - 123 cm		111 - 140 cm		
Tray height from footplate		43 - 66 cm		61 - 79 cm			72 - 96 cm		92 - 110 cm		
Max. load		30 kg		45 kg			65 kg		80 kg		
Weight		28 kg		31 kg			38 kg		42 kg		
HMF-No.		28.2901.3011									

* with bracket for rear support pad (standard) ** with bracket for rear support pad (flip away)



You can find more information on accessories in the price list.

