

Gives you flexibility.

GMFCS
Level II-III

kg
max. 60



65-140 cm

up to 10°



I am your dynamic standing trainer.

With me, you can comfortably explore your full range of movement. In my dynamic version, I tilt and move with you. And to make sure our joint standing training is really effective and lots of fun at the same time, I'm available via the stixi. app. There, you can not only track your progress with your therapy, but also play a cool movement game. This way, you can naturally extend your standing times while enjoying every moment.


Prefer something more classic? No problem. I'm also available in a non-dynamic version. The choice is yours!

new



Individual adjustability of resistance

(from grade 1-4) using four easy-to-reach dials.



stixi. – the brilliant new app
with interactive
game and standing diary.

EN
version
coming
soon



Therapy table
offers plenty of space
for activities.

Armrests

Mount for assisted communication

Hip support

Centre column
incl. protective cover.

Flexible knee positioners (FKP)
straps adapt perfectly to the
shape of the knee.

Integrated rails
for precise positioning of the heel edges.

Castors (75 mm)
with brakes.

Dynamic adjustment of the centre column
up to 10° in all directions.

**Slate grey vest
with manni embroidery**
is adjustable in size and height
and supports sensory perception.



**Protective pad for covering
the rotary handle**





Different foot positioning options allow for different therapeutic goals.

data	Size 0 (mini)	Size 1	Size 2
Body size	65 cm - 95 cm	85 cm - 115 cm	105 cm - 140 cm
Sternum height	42 cm - 62 cm	58 cm - 86 cm	76 cm - 100 cm
Table height from the base plate	48 cm - 72 cm	63 cm - 100 cm	80 cm - 115 cm
Total length	80 cm	88 cm	98 cm
Total width	50 cm	60 cm	70 cm
Max. load	20 kg	40 kg	60 kg
Weight with dynamic centre column	16 kg	19 kg	22 kg
Weight with fixed centre column	13 kg	15 kg	18 kg