

On the same wavelength.

GMFCS
Level III-V

kg
max. 114



I am your universal bathing and shower cradle.

You can relax beautifully on my comfortable seating surface. Here you can adjust my back so that it is best for you – even directly whilst bathing.

So that you also feel safe, there are different straps available for your chest, pelvis or legs. Of course, there is a headrest too. You can have me with or without a calf pad. My transfer sub-frame offers you ideal support when showering over the bathtub.

Bath sub-frame

puts you at the right height and relieves the strain on your parents' backs.



Shower sub-frame

with wide castors which can be braked; can also be used for transfer.





Headrest

can be adjusted in height and width.

Removable cover

for fast and uncomplicated cleaning in the washing machine.

Plastic sub-frame

can be adjusted in height.

Support feet

can be adjusted in height.

**Chest strap**

provides security and support. This is available with and without the adjustable side guidance and can also be used as a pelvic strap.

Optional leg strap

is adapted appropriately using hook and eye tape.

Optional calf pad

can be adjusted during use.

Back angle

can be adjusted with one hand.

**The angle of the calf support and seat**

can be adjusted independently of each other, even during care, so that you can sit or lie down in a stable and comfortable position.



Frame colour

☐ White



Transfer sub-frame
can be folded to save space.

In combination with the transfer frame, I can easily position you over the bathtub and even swivel by 90°. That makes showers such fun!



data	Size 1	Size 2	Size 3
Seat depth	27 cm	33 cm	43 cm
Seat width	38 cm		44 cm
Back height	58 cm	67 cm	81 cm
Seat height	5 cm / 18 cm		
Seat height with bath sub-frame	32 - 47 cm		
Seat height with shower sub-frame	70 cm		
Seat height with transfer sub-frame	53 - 75 cm		
Length of bathing chair, completely flat, without calf support	89 cm	104 cm	128 cm
Length of sun lounger, completely flat, with calf support	108 cm (43 inches)	124 cm	155 cm
Width	46 cm		53 cm
Length, completely flat	75 cm		
Height of bath sub-frame	20 cm		
Length of bath sub-frame	99 cm		
Width of bath sub-frame	44 cm		
Height of shower sub-frame	56 cm		
Width of shower sub-frame	69 cm		
Length of shower sub-frame	86 cm		
Max. load	34 kg	68 kg	114 kg
Back and seat angle	0°, 15°, 30°, 45°, 60°, 75°		
HMV-No.	04.40.04.0007		



You can find more information on accessories & harnesses
in the price list.

