



Pricelist 2023

Sitting + Positioning.



smilla.
04-07



madita.
08-14



madita-fun.
16-23



madita-fun mini.
24-27



madita hd.
28-31



madita maxi.
32-34



sofie.
36-42



activity chair.
44-47



tina.
48-50

Transport + Travelling.



sascha.
52-53



siggi combi isofix.
54-55



siggi 3 fix.
56-57



freddy.
58-61



fridolin maxi.
62-63



franka.
64-66



corzo xcountry go.
68-69



piper.
70-71

72–108

Walking + Riding.



malte.
72–74



malte outdoor.
75–77



marcy.
78–79



dynamic pacer.
80–85



momo dreirad.
86–91



follow-me momo.
92–93



momo motion.
94–100



momo therapierad.
102–108

110–129

Standing + Transfer.



till.
110–115



todd.
116–121



tim.
122–126



tram.
128–129

130–140

Bathroom + Toilet.



boris.
130–131



hts.
132–135



wave.
136–137



basti.
138–139



bibi.
140

141–147

General information.



Name embroidery.
141



Fabric collection.
142–143



Cleaning + disinfection.
144–145

smilla.
Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The convenient seating aid for everyday use.

- Sturdy wooden frame with integrated aluminium profiles
- Incl. height adjustable backrest and adjustable back angle
- Incl. seat and back pad
- Incl. height adjustable armrests
- Incl. height, angle and depth adjustable seat and footplate
- Incl. positioning castors
- Material: varnished beech



ATTENTION!

Available from **Q2** as **version 2023**
Including quick height adjustment
and even greater colour selection!


HMV-No.: 26.11.05.0034
Dimensions

		Size 0 Nursery	Size 0	Size 1 Nursery	Size 1	Size 2
A	Seat width	17 - 28 cm	17 - 28 cm	20 - 35 cm	20 - 35 cm	30 - 40 cm
B	Seat depth	16 - 26 cm	16 - 26 cm	25 - 35 cm	25 - 35 cm	30 - 40 cm
C	Max. seat height	45 cm	65 cm	45 cm	65 cm	65 cm
D	Back height	25 cm	25 cm	29 cm	29 cm	34 cm
E	Back angle	(-) 3° - 17,5°				
F	Tilt-in-space	(-) 5°, 0°, 7,5°, 15°				
G	Foot angle	(-) 5°, 0°, 7,5°, 15°				
H*	Min. footplate height	Standard	13 cm			
		with 50 mm dual castors	18 cm			
		with 75 mm dual castors	21 cm			
H*	Min. footplate height (with quick angle adjustment)	Standard	8 cm			
		with 50 mm dual castors	13 cm			
		with 75 mm dual castors	16 cm			
	Max. load	30 kg	30 kg	45 kg	45 kg	60 kg
	Weight	9,2 kg	9,7 kg	9,6 kg	10,2 kg	12 kg

* Max. adjustable lower leg length = selected / required seat height - min. footplate height (H)





Basic model		
smilla. size 0	smilla. size 1	smilla. size 2
Price / piece	Price / piece	Price / piece
☐ Size 0 Nursery (13 00 100) 1.205 €	☐ Size 1 Nursery (13 01 100) 1.205 €	☐ Size 2 (13 02 000) 1.445 €
☐ Size 0 (13 00 000) 1.205 €	☐ Size 1 (13 01 000) 1.205 €	

Frame colours	
Beech natural	Beech white
	
Price	Price
☐ 92 00 300 incl.	☐ 92 00 301 incl.

Covers					
Pink 2 (No. 115)	Azur (No. 110)	Green (No. 89)	Dark Grey (No. 92)	Flexi Blue (No. 44)	Name embroidery on the back pad please fill out the form on page 155! An embroidery of Flexi fabrics is not possible.
					
Price	Price	Price	Price	Upcharge	
☐ 91 00 115 incl.	☐ 91 00 110 incl.	☐ 91 00 089 incl.	☐ 91 00 092 incl.	☐ 13 00 044 250 €	

Accessories						
Headrest – padded and height adjustable		Lateral pads – padded and adjustable		Thigh guides – padded and adjustable		
						
	Adjustable height (in 3 cm steps)		Width 1	Width 2		Length
for size 0+1	12 – 21cm	for size 0	13 – 19cm	17 – 23cm	Size 0	17.5cm
for size 2	20.5 – 32.5cm	for size 1	20 – 26cm	24 – 30cm	Size 1	20cm
		for size 2	24 – 30cm	28 – 34cm	Size 2	22.5cm
	Price / piece		Price / pair			Price / pair
☐ for size 0 (13 00 001)	270 €	☐ Width 1 (13 01 002)	175 €	☐ Size 0 (13 00 004)		226 €
☐ for size 1 (13 01 001)	270 €	☐ Width 2 (13 01 003)	175 €	☐ Size 1 (13 01 004)		226 €
☐ for size 2 (13 02 001)	270 €			☐ Size 2 (13 02 004)		226 €

Configuration options - seat	
Standard seat pad – incl. cover (standard configuration, if nothing else stated)	Seat pad – incl. cover with anatomic leg guidance
	
Price / piece	Upcharge / piece
☐ for size 0 (13 00 021) incl.	☐ for size 0 (13 00 022) 222 €
☐ for size 1 (13 01 021) incl.	☐ for size 1 (13 01 022) 222 €
☐ for size 2 (13 02 021) incl.	☐ for size 2 (13 02 022) 222 €

Accessories

Quick angle adjustment of the footplate and seat – allows quick and easy angle adjustment of the footplate and seat to (-) 5°, 0°, 7,5° and 15° by straps (not possible in connection with quick height adjustment)



Upcharge / set

- ☐ for size 0 **(13 00 090)**
- ☐ for size 1 **(13 01 090)**
- ☐ for size 2 **(13 02 090)**

275 €
275 €
275 €

Quick height adjustment for footplate, seat and back – allows quick and easy height adjustment of the footplate, seat and back (not possible in connection with quick angle adjustment)



Upcharge / set

- ☐ **13 00 020**

68 €

Accessories

Wooden therapy tray – plug-on




	Width	Depth
for size 0	48 cm	37 cm
for size 1	52 cm	39 cm
for size 2	57 cm	39 cm

Price / piece

- ☐ for size 0 **(19 00 005)** 263 €
- ☐ for size 1 **(19 01 005)** 263 €
- ☐ for size 2 **(19 02 005)** 263 €

Large Wooden therapy tray – plug-on, with edge



	Width	Depth
for size 0	48 cm	45 cm
for size 1	52 cm	49 cm
for size 2	57 cm	49 cm

Price / piece

- ☐ for size 0 **(19 00 009)** 291 €
- ☐ for size 1 **(19 01 009)** 291 €
- ☐ for size 2 **(19 02 009)** 291 €

Abduction wedge – padded and pluggable



Price / piece

- ☐ for size 0 **(19 00 008)** 221 €
- ☐ for size 1+2 **(19 01 008)** 221 €

Accessories

Dual castors – 4 pieces



Price / set

- ☐ 50 mm – 2 with locking devices **(13 01 006)** 291 €
- ☐ 50 mm – all with locking devices **(13 01 016)** 328 €
- ☐ 75 mm – 2 with locking devices **(13 00 006)** 291 €
- ☐ 75 mm – all with locking devices **(13 00 016)** 328 €

Harness

The **smilla.** harness is assembled with the following harness assembly sets.

Further positioning aids and configuration options can be found from page 156.

Harness assembly kit for the seat
consisting of 2 brackets with 2 clip buckle



Price / set

- ☐ 20 mm **(19 00 010)** 44 €
- ☐ 25 mm **(19 01 010)** 44 €

Harness assembly kit for the back
consisting of belt holder with T-groove incl. 2 clip buckles and 2 clip buckles (drilled)




Price / set

- ☐ 20 mm **(19 00 011)** 149 €
- ☐ 25 mm **(19 01 011)** 149 €

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder, colour: black

	Size 0	Size 1	Size 2	Size 3
A	13cm	13cm	18cm	18cm
B	17.5cm	19cm	20cm	21cm
C	15cm	15cm	19cm	28cm
D	3.5cm	5cm	5cm	5cm
E	25.5cm	28cm	34cm	44cm
Strap width	2cm	2cm	2.5cm	2.5cm




Price / piece

☐ Size 0 (1300025)	177 €
☐ Size 1 (1301025)	177 €
☐ Size 2 (1302025)	177 €
☐ Size 3 (1303025)	177 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips, colour: black

	Size 0	Size 1	Size 2
A	11cm	11cm	18cm
B	21cm	21cm	22cm
C	12cm	12cm	15cm
D	2.5cm	3.5cm	3.5cm
E	25.5cm	28cm	34cm
Strap width	2cm	2cm	2.5cm




Price / piece

☐ Size 0 (1300026)	206 €
☐ Size 1 (1301026)	206 €
☐ Size 2 (1302026)	206 €

Harness

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back, colour: black

	Size 0	Size 1	Size 2
A	23cm	25cm	28cm
B	15cm	17cm	20cm
C	33cm	35cm	38cm
D	24cm	26cm	28cm
Strap width	2cm	2cm	2.5cm




Price / piece

☐ Size 0 (1300027)	148 €
☐ Size 1 (1301027)	148 €
☐ Size 2 (1302027)	148 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back, colour: black

	Size 1	Size 2	Size 3
A	24cm	29cm	33cm
B	14cm	17cm	22cm
C	35cm	43cm	49cm
Strap width	2cm	2.5cm	2.5cm



Price / piece

☐ Size 1 (1301028)	148 €
☐ Size 2 (1302028)	148 €
☐ Size 3 (1303028)	148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides, colour: black

	Size 0	Size 1	Size 2
Pad width	4cm	4cm	5cm
Pad length per side	16cm	18cm	28cm
Strap length per side	30cm	30cm	40cm
Strap width	2cm	2.5cm	2.5cm

Price / piece

☐ Size 0 (1300029)	89 €
☐ Size 1 (1301029)	89 €
☐ Size 2 (1302029)	89 €




Harness

4-point pelvic harness – with pad and adjustable clip on both sides, colour: black

	Size 0	Size 1	Size 2
Pad width	4cm	4cm	5cm
Pad length per side	16cm	18cm	28cm
Strap length per side	30cm	30cm	40cm
Strap width	2cm	2.5cm	2.5cm

Price / piece

☐ Size 0 (1300030)	155 €
☐ Size 1 (1301030)	155 €
☐ Size 2 (1302030)	155 €



Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The mobile therapy chair which grows with the child.

- Incl. armrests, footplate with slots for Harness, push handle, standard seat pad with cover
- Adjustable seat height with gas pressure spring support
- 4 castors (75 mm), 2 with locking devices
- Adjustment options:
 - Tilt-in-space (via gas pressure spring 500N)
 - Back angle
 - Seat depth and seat width
 - Angle, height and width adjustable armrests
 - Height adjustable backrest
 - Folding footplate, knee angle and height adjustable

HMV-No.: 26.11.05.0021


Basic model

madita. size 0	madita. size 1	madita. size 1b	madita. size 2
Price / piece	Price / piece	Price / piece	Price / piece
☐ 78 00 000 1.735 €	☐ 78 01 000 1.735 €	☐ 78 01 200 1.860 €	☐ 78 02 000 1.985 €

Frame colours

Anthracite grey		Telegrey 4 (white grey)	
Price		Price	
☐ 92 00 017 incl.		☐ 92 00 011 incl.	

Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:


A Axilla height	_____	D Pelvic width	_____
B Chest width	_____	E Seat depth	_____
C Back height	_____	F Lower leg length	_____

☐ **73 00 080**

Price
incl.

Dimensions

		Size 0	Size 1	Size 1b	Size 2
A	Seat depth*	15–25 cm	21–30 cm	26–36 cm	30–42 cm
B	Seat width	16–26 cm	20–30 cm	20–30 cm	28–38 cm
C	Back height with short back	35–43 cm	35–43 cm	42–52 cm	42–52 cm
C	Back height with long back	42–52 cm	42–52 cm	53–63 cm	53–63 cm
D	Seat height	47,5–58 cm	47,5–58 cm	47,5–58 cm	47,5–58 cm
E	Lower leg length	15–25 cm	19–32 cm	27–49 cm	27–49 cm
F	Total width	45 cm	45 cm	53 cm	53 cm
G	Length of base frame	66 cm	66 cm	77 cm	77 cm
	Tilt-in-space	(-)5° – 37°	(-)5° – 37°	(-)5° – 37°	(-)5° – 37°
	Back angle	(-)5° – 25°	(-)5° – 25°	(-)5° – 25°	(-)5° – 25°
	Max. load	30 kg	30 kg	45 kg	50 kg
	Weight	16 kg	17 kg	18 kg	19 kg



* reduced by 1 cm in combination with adjustable pelvic guide

Covers

Pink 2 (No. 115) 	Deep Sea Blue (No. 111) 	Leaf-Green (No. 123) 	Slate grey (No. 124) 
Price	Price	Price	Price
☐ 91 00 115	☐ 91 00 111	☐ 91 00 123	☐ 91 00 124
incl.	incl.	incl.	incl.

Covers

Name embroidery

Minor position

Flexi Blue (No. 44) 	Cover according to customer request	Name embroidery on the back pad please fill out the form on page 141! An embroidery of Flexi fabrics is not possible. 	Without seat and back upholstery and cover for the installation of custom cushions
Upcharge	Upcharge		Price reduction
☐ 91 00 044	☐ 91 00 999		☐ 73 00 059
250 €	180 €		Net -104 €

Back – incl. pad and cover		for size 0	for size 1	for size 1b	for size 2
	extra short back:	32 – 35 cm	–	35 – 43 cm	–
	short back:	35 – 43 cm	35 – 43 cm	42 – 52 cm	42 – 52 cm
	long back:	42 – 52 cm	42 – 52 cm	53 – 63 cm	53 – 63 cm
		Price / piece	Price / piece	Price / piece	Price / piece
	extra short back:	☐ 78 00 001 185 €		☐ 78 01 201 185 €	
	short back:	☐ 78 00 002 185 €	☐ 78 01 002 185 €	☐ 78 01 202 185 €	☐ 78 02 002 208 €
	long back:	☐ 78 00 003 185 €	☐ 78 01 003 185 €	☐ 78 01 203 185 €	☐ 78 02 003 208 €

		for size 0	for size 1	for size 1b	for size 2
		Price / piece	Price / piece	Price / piece	Price / piece
	for extra short back:	☐ 78 00 004 120 €		☐ 78 01 204 120 €	
	for short back:	☐ 78 00 005 120 €	☐ 78 01 005 120 €	☐ 78 01 205 120 €	☐ 78 02 005 120 €
	for long back:	☐ 78 00 006 120 €	☐ 78 01 006 120 €	☐ 78 01 206 120 €	☐ 78 02 006 120 €

Active back incl. pad and cover		for size 0	for size 1	for size 1b	for size 2
		Price / piece		Price / piece	
		☐ 78 00 007 235 €	☐ 78 01 007 235 €	☐ 78 02 007 235 €	☐ 78 03 007 235 €

Back unit – fixed (standard configuration = without push handles, if nothing else stated)			Back unit – removable			Adjustable limitation of the back angle – to determine the individual max. back angle		
Price / piece			Price / piece			Price / Set		
☐ for size 0+1 with push handles	(73 01 004)	incl.	☐ for size 0+1 with push handles	(73 01 006)	370 €	☐ 73 01 111		62 €
☐ for size 2 with push handles	(73 02 004)	incl.	☐ for size 2 with push handles	(73 02 006)	370 €			
☐ for size 0+1 with push loop	(73 01 053)	incl.						

	Armrest pad size 1 – (Standard configuration for basic model size 0, if nothing else stated) Length: 20 cm Width: 6 cm	Armrest pad size 2 – (Standard configuration for basic model size 1+2, if nothing else stated) Length: 27 cm Width: 7 cm	Armrest pad size 3 Length: 34 cm Width: 8 cm
	Price / pair	Price / pair	Price / pair
	☐ 78 01 008 incl.	☐ 78 02 008 incl.	☐ 78 03 008 incl.


Headrest (standard) padded, incl. cover, height and angle adjustable 	Headrest (Shell-shaped) padded, incl. cover, height and angle adjustable 	Headrest (Lateral guide) padded, incl. cover, height and angle adjustable 	Headrest (Lateral guide, angle adjustable) padded, incl. cover, height, depth and angle adjustable side guides are each angle adjustable A = Width contact surface B = Height headrest C = Width lateral guide																																
Price / piece	Price / piece	Price / piece	<table border="1"> <thead> <tr> <th></th> <th>A</th> <th>B</th> <th>C</th> </tr> </thead> <tbody> <tr> <td>Size 1</td> <td>13,5cm</td> <td>10,5cm</td> <td>10cm</td> </tr> <tr> <td>Size 2</td> <td>15cm</td> <td>10,5cm</td> <td>13,5cm</td> </tr> <tr> <td>Size 3</td> <td>17cm</td> <td>12,5cm</td> <td>16,5cm</td> </tr> <tr> <td></td> <td colspan="3">Price / piece</td> </tr> <tr> <td>☐ Size 1 (7801012)</td> <td colspan="2"></td> <td>588 €</td> </tr> <tr> <td>☐ Size 2 (7802012)</td> <td colspan="2"></td> <td>588 €</td> </tr> <tr> <td>☐ Size 3 (7803012)</td> <td colspan="2"></td> <td>588 €</td> </tr> </tbody> </table>		A	B	C	Size 1	13,5cm	10,5cm	10cm	Size 2	15cm	10,5cm	13,5cm	Size 3	17cm	12,5cm	16,5cm		Price / piece			☐ Size 1 (7801012)			588 €	☐ Size 2 (7802012)			588 €	☐ Size 3 (7803012)			588 €
	A	B	C																																
Size 1	13,5cm	10,5cm	10cm																																
Size 2	15cm	10,5cm	13,5cm																																
Size 3	17cm	12,5cm	16,5cm																																
	Price / piece																																		
☐ Size 1 (7801012)			588 €																																
☐ Size 2 (7802012)			588 €																																
☐ Size 3 (7803012)			588 €																																
☐ for size 0 (7800009) 254 € ☐ for size 1 (7801009) 254 € ☐ for size 2 (7802009) 254 €	☐ for size 0+1 (7801010) 360 € ☐ for size 2 (7802010) 360 €	☐ for size 0+1 (7801011) 340 € ☐ for size 2 (7802011) 340 €																																	

Headrests		Seat angle adjustment		
Universal adapter for tailored headrests in the version 15 x 15 mm	Trigger lever – mounted directly on the gas spring (standard configuration, if nothing else stated)	Trigger lever – mounted on the push bar/handle (not possible with removable back unit)	Trigger lever with detachable Bowden cable (for the removable back unit)	Safety bolt – for release lever, prevents unintentional release of height and angle adjustments
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
☐ 78 00 013 85 €	☐ 73 00 085 incl.	☐ 73 01 085 incl.	☐ 73 02 085 35 €	☐ 73 01 086 75 €

Configuration possibilities - lateral pads (Step 1: Select the bracket)

Chest width – mounted on madita. Size 0	Chest width – mounted on madita. Size 1 + 1b	Chest width – mounted on madita. Size 2	Bracket	Upcharge / pair
13 – 27 cm	16 – 30 cm	26 – 38 cm	☐ Fixed (78 00 014)	incl.
			☐ Flip away (78 01 014)	259 €

(Step 2: Select the lateral pads)

	Lateral pads size 0 padded, incl. cover (Standard configuration for madita. size 0, if nothing else stated) Height: 7 cm Depth: 9 cm	Lateral pads size 1 padded, incl. cover (Standard configuration for madita. size 1, if nothing else stated) Height: 8 cm Depth: 13 cm	Lateral pads size 2 padded, incl. cover (Standard configuration for madita. size 2, if nothing else stated) Height: 11 cm Depth: 16 cm
	Price / pair	Price / pair	Price / pair
	☐ 78 00 015 106 €	☐ 78 01 015 106 €	☐ 78 02 015 106 €


(Step 2: Select the lateral pads)


	Lateral pads size 0, long (Trunk guidance) padded, incl. cover Height: 13 cm Depth: 7 cm Min. adjustable anderarm height in combination with adjustable pelvic guide: 22 cm	Lateral pads size 1, long (Trunk guidance) padded, incl. cover Height: 15 cm Depth: 8 cm Min. adjustable anderarm height in combination with adjustable pelvic guide: 24 cm	Lateral pads size 2, long (Trunk guidance) padded, incl. cover Height: 17 cm Depth: 9 cm Min. adjustable anderarm height in combination with adjustable pelvic guide: 26 cm
	i Lateral pads (long) with flip away bracket can be folded 20° to the side.		
	Price / pair	Price / pair	Price / pair
	☐ 78 00 016 186 €	☐ 78 01 016 186 €	☐ 78 02 016 186 €

Tailored lateral pads
Tailored lateral pads, long


Height: ____ cm Depth: ____ cm	Height: ____ cm Depth: ____ cm
Price / pair	Price / pair
☐ from the Schuchmann company , padded, incl. cover (78 08 015) 250 € ☐ from the medical supply store , blank without upholstery + cover (78 09 015) 180 €	☐ from the Schuchmann company , padded, incl. cover (78 08 016) 375 € ☐ from the medical supply store , blank without upholstery + cover (78 09 016) 250 €
The blank consisting of 3 mm aluminium serves as a template for self-cutting!	



Seat option

Standard seat pad – incl. cover (standard configuration, if nothing else stated) 	Seat pad – incl. cover with anatomic leg guidance 	Adjustable pelvic guidance – padded, incl. cover  
Price / piece	Upcharge / piece	Price / piece
☐ for size 0 (78 00 017) incl. ☐ for size 1 (78 01 017) incl. ☐ for size 1b (78 01 217) incl. ☐ for size 2 (78 02 017) incl.	☐ for size 0 (78 00 018) 218 € ☐ for size 1 (78 01 018) 218 € ☐ for size 1b (78 01 218) 218 € ☐ for size 2 (78 02 018) 218 €	☐ for size 0 – seat width A = 15 – 24 cm (78 00 020) 384 € ☐ for size 1 + 1b – seat width A = 18 – 28 cm (78 01 020) 384 € ☐ for size 2 – seat width A = 26 – 36 cm (78 02 020) 408 €

Seat option				Accessories	
Abductable leg guidance with leg length compensation, padded, incl. cover (already includes thigh guides)				Knee abduction wedge pluggable using quick-release axle, padded, incl. cover (not usable with abductable leg guidance)	
					
	for size 0	for size 1+1b	for size 2		
A	0° – 20°	0° – 20°	0° – 20°	Size 0	Width
B	7 – 10 cm	9 – 12 cm	10 – 14 cm	Size 1	Depth
C	3 cm	3 cm	3 cm	Size 2	
Price / piece				Price / piece	
☐ for size 0 (7800019) 520 € ☐ for size 1 (7801019) 520 € ☐ for size 1 b (7801219) 520 € ☐ for size 2 (7802019) 540 €				☐ Size 0 (7300023) 225 € ☐ Size 1 (7301023) 225 € ☐ Size 2 (7302023) 225 €	

Thigh guides			
Thigh guides for size 0 padded, incl. cover, adjustable angle, depth and width		Thigh guides for size 1 + 1b padded, incl. cover, adjustable angle, depth and width	
12 cm = standard configuration for size 0		14 cm = standard configuration for size 1 + 1b	
Price / pair		Price / pair	
☐ 12 cm (7300114)	192 €	☐ 14 cm (7301114)	192 €
☐ 14 cm (7300115)	192 €	☐ 17 cm (7301115)	192 €
		Thigh guides for size 2 padded, incl. cover, adjustable angle, depth and width	
		17 cm = standard configuration for size 2	
		Price / pair	
		☐ 17 cm (7302114)	197 €
		☐ 22 cm (7302115)	197 €

Footplate				Calf plant							
				Calf plant for Footplate covered with black flexi-fabric							
											
Footplate standard – without fixed position, foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps		Footplate with lock – foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps.		Footplate with angle adjustment foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps.							
Price / piece		Upcharge / piece		Upcharge / piece		Price / piece					
☐ for size 0	(7800021)	incl.	☐ for size 0	(7800022)	52 €	☐ for size 0	(7800023)	96 €	☐ for size 0	(7800024)	206 €
☐ for size 1	(7801021)	incl.	☐ for size 1	(7801022)	52 €	☐ for size 1	(7801023)	96 €	☐ for size 1 + 1b	(7801024)	206 €
☐ for size 1b	(7801221)	incl.	☐ for size 1b	(7801222)	52 €	☐ for size 1b	(7801223)	96 €	☐ for size 2	(7802024)	206 €
☐ for size 2	(7802021)	incl.	☐ for size 2	(7802022)	52 €	☐ for size 2	(7802023)	96 €			

Configuration footplates												
Knee angle to accommodate the footplates individually adjustable in height, width and angle			Split footplates – foldable, fixed mounted on the footplate holder, with non-slip surface and each 4 slots for the for the fixation of straps			Foot cups Foldable, fixed mounted on the footplate holder, with magnetic lock						
												
							Min. front width	Max. front width	Min. rear width	Max. rear width	Length	
			Size 1	11,4 cm	10 cm	18,4 cm	Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm
			Size 2	13,8 cm	13,2 cm	21,4 cm	Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm
							Size 3	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm
Price / piece			Price / pair			Upcharge / pair						
☐ for size 0	(7800 025)	285 €	☐ Size 1	(7801 026)		incl.	☐ Size 1	(7801 027)				168 €
☐ for size 1	(7801 025)	285 €	☐ Size 2	(7802 026)		incl.	☐ Size 2	(7802 027)				168 €
☐ for size 1b	(7801 225)	285 €					☐ Size 3	(7803 027)				168 €
☐ for size 2	(7802 025)	285 €										

Accessories

Therapy tray – made of wood with edge, flip away 			Therapy tray – transparent, with edge, flip away 			Therapy tray – made of wood with edge, prefixable, fixed on both sides 			Therapy tray – transparent, with edge, prefixable, fixed on both sides 		
	Width	Depth		Width	Depth		Width	Depth		Width	Depth
for size 0+1	50 cm	40 cm	for size 0+1	50 cm	40 cm	for size 0+1	50 cm	40 cm	for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm	for size 2	60 cm	50 cm	for size 2	60 cm	50 cm	for size 2	60 cm	50 cm
Price / piece			Price / piece			Price / piece			Price / piece		
☐ for size 0 (78 00 028)		249 €	☐ for size 0 (78 00 029)		340 €	☐ for size 0 (78 00 038)		265 €	☐ for size 0 (78 00 039)		380 €
☐ for size 1 (78 01 028)		249 €	☐ for size 1 (78 01 029)		340 €	☐ for size 1 (78 01 038)		265 €	☐ for size 1 (78 01 039)		380 €
☐ for size 2 (78 02 028)		249 €	☐ for size 2 (78 02 029)		340 €	☐ for size 2 (78 02 038)		265 €	☐ for size 2 (78 02 039)		380 €
Add-on option:											
☐ Protection against unintentional folding up				(78 01 036)	98 €						

Accessories

Central parking brake for rear castors incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system 		4 Castors – 2 with locking devices 		4 Castors – all with locking devices 	
		Upcharge / set		Upcharge / set	
☐ Parking lever mounted left		(78 00 030)	425 €	☐ 75 mm (73 00 058)	
☐ Parking lever mounted right		(78 01 030)	425 €	☐ 100 mm (73 01 058)	
☐ Parking lever mounted left; castors front with brakes		(78 00 031)	459 €	☐ 125 mm (73 02 058)	
☐ Parking lever mounted right; castors front with brakes		(78 01 031)	459 €	☐ 75 mm (73 00 060)	
				☐ 100 mm (73 01 060)	
				☐ 125 mm (73 02 060)	

Accessories

Wheel base extension inserts at rear, length = 10 cm (not possible together with central parking brake)	Padded cover for footplate covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)	Padded cover for split footplate covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)	Table upholstery – covered with black Flexi cover, height: 1 cm 
Price / pair	Price / piece	Price / pair	Price / piece
☐ 73 01 028 168 €	☐ for size 0 (78 00 032) 118 € ☐ for size 1 & 1b (78 01 032) 118 € ☐ for size 2 (78 02 032) 118 €	☐ for size 1 (78 01 033) 169 € ☐ for size 2 (78 02 033) 169 €	☐ for size 0 (73 00 054) 124 € ☐ for size 1 (73 01 054) 124 € ☐ for size 2 (73 02 054) 124 €

Accessories

Incontinence coating in the seat area	Harness assembly kit for the seat consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)	Harness assembly kit for the back consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)
		
Price / piece	Upcharge / set	Upcharge / set
☐ 73 00 020	☐ 20 mm (73 00 025)	☐ 20 mm (73 00 026)
196 €	44 €	44 €
	☐ 25 mm (73 01 025)	☐ 25 mm (73 01 026)
	44 €	44 €


Harness

The **madita.**-harness is assembled with the harness assembly sets (see page 13).

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 0	Size 1	Size 2	Size 3
A	13 cm	13 cm	18 cm	18 cm
B	17,5 cm	19 cm	20 cm	21 cm
C	15 cm	15 cm	19 cm	28 cm
D	3,5 cm	5 cm	5 cm	5 cm
E	25,5 cm	28 cm	34 cm	44 cm
Strap width	2 cm	2 cm	2,5 cm	2,5 cm




Price / piece

☐ Size 0 (73 00 032)	177 €
☐ Size 1 (73 01 032)	177 €
☐ Size 2 (73 02 032)	177 €
☐ Size 3 (73 03 032)	177 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 0	Size 1	Size 2
A	11 cm	11 cm	18 cm
B	21 cm	21 cm	22 cm
C	12 cm	12 cm	15 cm
D	2,5 cm	3,5 cm	3,5 cm
E	25,5 cm	28 cm	34 cm
Strap width	2 cm	2 cm	2,5 cm




Price / piece

☐ Size 0 (73 00 056)	206 €
☐ Size 1 (73 01 056)	206 €
☐ Size 2 (73 02 056)	206 €

Harness

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 0	Size 1	Size 2
A	23 cm	25 cm	28 cm
B	15 cm	17 cm	20 cm
C	33 cm	35 cm	38 cm
D	24 cm	26 cm	28 cm
Strap width	2 cm	2 cm	2,5 cm




Price / piece

☐ Size 0 (73 00 033)	148 €
☐ Size 1 (73 01 033)	148 €
☐ Size 2 (73 02 033)	148 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2,5 cm	2,5 cm



Price / piece

☐ Size 1 (73 01 034)	148 €
☐ Size 2 (73 02 034)	148 €
☐ Size 3 (73 03 034)	148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2,5 cm	2,5 cm

Price / piece

☐ Size 0 (73 00 035)	89 €
☐ Size 1 (73 01 035)	89 €
☐ Size 2 (73 02 035)	89 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2,5 cm	2,5 cm

Price / piece

☐ Size 0 (73 00 036)	155 €
☐ Size 1 (73 01 036)	155 €
☐ Size 2 (73 02 036)	155 €

Harness

Ankle cuffs – with velcro fastener and idler pulley colour: black (not possible in combination with padded cover for footplate/split footplate)




	Ankle circumference	Strap width
Size 1	12 - 19 cm	2,5 cm
Size 2	15 - 23 cm	2,5 cm
Size 3	19 - 27 cm	4 cm

Price / pair

☐ Size 1 (73 01 061)	212 €
☐ Size 2 (73 02 061)	212 €
☐ Size 3 (73 03 061)	212 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



Velcro strip to open

	A	Strap width
Size 1	6,5 cm	2,5 cm
Size 2	8 cm	2,5 cm
Size 3	9,5 cm	2,5 cm

Price / pair

☐ Size 1 (78 01 034)	196 €
☐ Size 2 (78 02 034)	196 €
☐ Size 3 (78 03 034)	196 €



flexible knee positioning

The new Flexible knee pad.

Can be selected directly from our aids **till.** (page 110),
todd. (page 116) and **tim.** (page 122).

madita-fun.

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The modern and integrative therapy chair.

- Incl. armrests, footplate with slots for harness, push handle, standard seat pad with cover
- Height adjustment via gas pressure spring 1200 N
- 4 castors (75 mm), 2 with locking devices
- Adjustment options:
 - Tilt-in-space (via gas pressure spring 500 N)
 - Back angle
 - Seat depth and seat width
 - Angle, height and width adjustable armrests
 - Height adjustable backrest
 - Folding footplate, knee angle and height adjustable

HMV-No.: 26.11.05.0022


Basic model

madita-fun. size 0	madita-fun. size 1	madita-fun. size 1b	madita-fun. size 2
Price / piece	Price / piece	Price / piece	Price / piece
☐ 79 00 000 2.445 €	☐ 79 01 000 2.445 €	☐ 79 01 200 2.635 €	☐ 79 02 000 2.750 €

Frame colours

Anthracite grey	Telegrey 4 (white grey)
☐	☐
Price	Price
☐ 92 00 017 incl.	☐ 92 00 011 incl.

Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:


A Axilla height	_____	D Pelvic width	_____
B Chest width	_____	E Seat depth	_____
C Back height	_____	F Lower leg length	_____

☐ **74 00 080**

Price
incl.

Dimensions

		Size 0	Size 1	Size 1b	Size 2
A	Seat depth*	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
B	Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
C	Back height with short back	35 - 43 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
C	Back height with long back	42 - 52 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
D	Seat height	25 - 62 cm	25 - 62 cm	25 - 62 cm	25 - 62 cm
	Seat height with electric height adjustment	32 - 62 cm	32 - 62 cm	32 - 62 cm	32 - 62 cm
	Seat height with hydraulic height adjustment	27 - 63 cm	27 - 63 cm	27 - 63 cm	27 - 63 cm
E	Lower leg length	15 - 25 cm	19 - 32 cm	27 - 49 cm	27 - 49 cm
F	Total width	54 cm	54 cm	64 cm	64 cm
G	Length of base frame	61 cm	61 cm	75 cm	75 cm
	Total length with combi subframe	76 cm	76 cm	86 cm	86 cm
	Tilt-in-space	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°
	Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
	Max. load	30 kg	30 kg	45 kg	50 kg
	Weight	19 kg	20 kg	22 kg	24 kg



* reduced by 1 cm in combination with adjustable pelvic guide

Covers

Pink 2 (No. 115)		Deep Sea Blue (No. 111)		Leaf-Green (No. 123)		Slate grey (No. 124)	
	Price		Price		Price		Price
☐ 9100 115	incl.	☐ 9100 111	incl.	☐ 9100 123	incl.	☐ 9100 124	incl.

Covers

Name embroidery

Minor position

Flexi Blue (No. 44)		Cover according to customer request	Name embroidery on the back pad please fill out the form on page 141! An embroidery of Flexi fabrics is not possible.		Without seat and back upholstery and cover – for the installation of custom cushions
	Upcharge		Upcharge		Price reduction
☐ 9100 044	250 €.	☐ 9100 999	180 €		☐ 7300 059 Net -104 €

Back – incl. pad and cover		for size 0	for size 1	for size 1b	for size 2
	extra short back:	32 – 35 cm	–	35 – 43 cm	–
	short back:	35 – 43 cm	35 – 43 cm	42 – 52 cm	42 – 52 cm
	long back:	42 – 52 cm	42 – 52 cm	53 – 63 cm	53 – 63 cm
		Price / piece	Price / piece	Price / piece	Price / piece
	extra short back:	☐ 78 00 001 185 €		☐ 78 01 201 185 €	
	short back:	☐ 78 00 002 185 €	☐ 78 01 002 185 €	☐ 78 01 202 185 €	☐ 78 02 002 208 €
	long back:	☐ 78 00 003 185 €	☐ 78 01 003 185 €	☐ 78 01 203 185 €	☐ 78 02 003 208 €

		for size 0	for size 1	for size 1b	for size 2
		Price / piece	Price / piece	Price / piece	Price / piece
	for extra short back:	☐ 7800004 120 €		☐ 7801204 120 €	
	for short back:	☐ 7800005 120 €	☐ 7801005 120 €	☐ 7801205 120 €	☐ 7802005 120 €
	for long back:	☐ 7800006 120 €	☐ 7801006 120 €	☐ 7801206 120 €	☐ 7802006 120 €

Active back incl. pad and cover		for size 0	for size 1	for size 1b	for size 2
		Price / piece		Price / piece	
		☐ 78 00 007 235 €	☐ 78 01 007 235 €	☐ 78 02 007 235 €	

	Armrest pad size 1 – (standard configuration for Basic model Size 0, if nothing else stated) Length: 20 cm Width: 6 cm		Armrest pad size 2 – (standard configuration for Basic model Size 1+2, if nothing else stated) Length: 27 cm Width: 7 cm		Armrest pad size 3	
	Price / pair		Price / pair		Price / pair	
	☐ 78 01 008 incl.		☐ 78 02 008 incl.		☐ 78 03 008 incl.	

Back unit – fixed (standard configuration = without push handles, if nothing else stated)			Back unit – removable			Adjustable limitation of the back angle – to determine the individual max. back angle		
Price / piece			Price / piece			Price / Set		
☐ for size 0+1 with push handles	(7301004)	incl.	☐ for size 0+1 with push handles	(7301006)	370 €	☐ 7301111		62 €
☐ for size 2 with push handles	(7302004)	incl.	☐ for size 2 with push handles	(7302006)	370 €			
☐ for size 0+1 with push loop	(7301053)	incl.						

Headrest (Standard) padded, incl. cover, height and angle adjustable	Headrest (Shell-shaped) padded, incl. cover, height and angle adjustable	Headrest (Lateral guide) padded, incl. cover, height and angle adjustable	Headrest (lateral guide, angle adjustable) padded, incl. cover, height, depth and angle adjustable side guides are each angle adjustable																																
Price / piece	Price / piece	Price / piece	<table border="1"> <thead> <tr> <th></th> <th>A</th> <th>B</th> <th>C</th> </tr> </thead> <tbody> <tr> <td>Size 1</td> <td>13,5cm</td> <td>10,5cm</td> <td>10cm</td> </tr> <tr> <td>Size 2</td> <td>15cm</td> <td>10,5cm</td> <td>13,5cm</td> </tr> <tr> <td>Size 3</td> <td>17cm</td> <td>12,5cm</td> <td>16,5cm</td> </tr> <tr> <td></td> <td colspan="3">Price / piece</td> </tr> <tr> <td>☐ Size 1</td> <td>(78 01 012)</td> <td>588 €</td> <td></td> </tr> <tr> <td>☐ Size 2</td> <td>(78 02 012)</td> <td>588 €</td> <td></td> </tr> <tr> <td>☐ Size 3</td> <td>(78 03 012)</td> <td>588 €</td> <td></td> </tr> </tbody> </table>		A	B	C	Size 1	13,5cm	10,5cm	10cm	Size 2	15cm	10,5cm	13,5cm	Size 3	17cm	12,5cm	16,5cm		Price / piece			☐ Size 1	(78 01 012)	588 €		☐ Size 2	(78 02 012)	588 €		☐ Size 3	(78 03 012)	588 €	
	A	B	C																																
Size 1	13,5cm	10,5cm	10cm																																
Size 2	15cm	10,5cm	13,5cm																																
Size 3	17cm	12,5cm	16,5cm																																
	Price / piece																																		
☐ Size 1	(78 01 012)	588 €																																	
☐ Size 2	(78 02 012)	588 €																																	
☐ Size 3	(78 03 012)	588 €																																	
☐ for size 0 (78 00 009) 254 € ☐ for size 1 (78 01 009) 254 € ☐ for size 2 (78 02 009) 254 €	☐ for size 0+1 (78 01 010) 360 € ☐ for size 2 (78 02 010) 360 €	☐ for size 0+1 (78 01 011) 340 € ☐ for size 2 (78 02 011) 340 €	A = Width contact surface B = Height headrest C = Width lateral guide																																


Headrests		Seat angle adjustment		
Universal adapter for tailored headrests in the version 15 x 15 mm	Trigger lever – mounted directly on the gas spring (standard configuration, if nothing else stated)	Trigger lever – mounted on the push bar/handle (not possible with removable back unit)	Trigger lever with detachable Bowden cable (for the removable back unit)	Safety bolt – for release lever, prevents unintentional release of height and angle adjustments
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
☐ 78 00 013 85 €	☐ 74 00 085 incl.	☐ 74 01 085 incl.	☐ 74 02 085 35 €	☐ 73 01 086 75 €

Configuration possibilities - lateral pads (Step 1: Select the bracket)				
Chest width – mounted on madita. Size 0	Chest width – mounted on madita. Size 1 + 1b	Chest width – mounted on madita. Size 2	Bracket	Upcharge / pair
13 – 27 cm	16 – 30 cm	26 – 38 cm	☐ Fixed (78 00 014)	incl.
			☐ Flip away (78 01 014)	259 €

(Step 2: Select the lateral pads)			
	Lateral pads size 0 padded, incl. cover (Standard configuration for madita-fun. Size 0 , if nothing else stated) Height: 7 cm Depth: 9 cm	Lateral pads size 1 padded, incl. cover (Standard configuration for madita-fun. Size 1 , if nothing else stated) Height: 8 cm Depth: 13 cm	Lateral pads size 2 padded, incl. cover (Standard configuration for madita-fun. Size 2 , if nothing else stated) Height: 11 cm Depth: 16 cm
	Price / pair	Price / pair	Price / pair
	☐ 78 00 015 106 €	☐ 78 01 015 106 €	☐ 78 02 015 106 €

(Step 2: Select the lateral pads)			
	Lateral pads size 0, long (Trunk guidance) padded, incl. cover Height: 13 cm Depth: 7 cm Min. adjustable underarm height in combination with adjustable pelvic guide: 22 cm	Lateral pads size 1, long (Trunk guidance) padded, incl. cover Height: 15 cm Depth: 8 cm Min. adjustable underarm height in combination with adjustable pelvic guide: 24 cm	Lateral pads size 2, long (Trunk guidance) padded, incl. cover Height: 17 cm Depth: 9 cm Min. adjustable underarm height in combination with adjustable pelvic guide: 26 cm
	i Lateral pads (long) with flip away bracket can be folded 20° to the side.		
	Price / pair	Price / pair	Price / pair
	☐ 78 00 016 186 €	☐ 78 01 016 186 €	☐ 78 02 016 186 €

Tailored lateral pads	Tailored lateral pads, long
Height: ____ cm Depth: ____ cm	Height: ____ cm Depth: ____ cm
Price / pair	Price / pair
☐ from the Schuchmann company , padded, incl. cover (78 08 015) 250 €	☐ from the Schuchmann company , padded, incl. cover (78 08 016) 375 €
☐ from the medical supply store , blank without upholstery + cover (78 09 015) 180 €	☐ from the medical supply store , blank without upholstery + cover (78 09 016) 250 €
The blank consisting of 3 mm aluminium serves as a template for self-cutting!	


Seat option		
Standard seat pad – incl. cover (Standard configuration, if nothing else stated)	Seat pad – incl. cover with anatomic leg guidance	Adjustable pelvic guidance – padded, incl. cover
		 
Price / piece	Upcharge / piece	Price / piece
☐ for size 0 (78 00 017) incl.	☐ for size 0 (78 00 018) 218 €	☐ for size 0 – seat width A = 15 – 24 cm (78 00 020) 384 €
☐ for size 1 (78 01 017) incl.	☐ for size 1 (78 01 018) 218 €	☐ for size 1 + 1 b – seat width A = 18 – 28 cm (78 01 020) 384 €
☐ for size 1 b (78 01 217) incl.	☐ for size 1 b (78 01 218) 218 €	☐ for size 2 – seat width A = 26 – 36 cm (78 02 020) 408 €
☐ for size 2 (78 02 017) incl.	☐ for size 2 (78 02 018) 218 €	

Accessories

with leg length compensation, padded,
incl. cover (already includes thigh guides)

Price / piece

☐ for size 0	(7800019)	520 €
☐ for size 1	(7801019)	520 €
☐ for size 1 b	(7801219)	520 €
☐ for size 2	(7802019)	540 €



A = Abduction angle
B = Thigh width
C = Leg length compensation

pluggable using quick-release axle, padded, incl. cover (not usable with abductible leg guidance)



Case	15.000	10.000	5.000	Price / piece
------	--------	--------	-------	---------------

☐ Size 0	(73 00 023)	225 €
☐ Size 1	(73 01 023)	225 €
☐ Size 2	(73 02 023)	225 €

padding, incl. cover, adjustable angle,
depth and width

12 cm = standard configuration for Size 0

		Price / pair
☐ 12 cm	(73 00 114)	192 €
☐ 14 cm	(73 00 115)	192 €

padding, incl. cover, adjustable angle,
depth and width

14 cm = standard configuration for Size 1 + 1 b

	Price / pair
☐ 14 cm (7301 114)	192 €
☐ 17 cm (7301 115)	192 €

padding, incl. cover, adjustable angle,
depth and width

17 cm = Standard configuration for Size 2

	Price / pair
☐ 17 cm (73 02 114)	197 €
☐ 22 cm (73 02 115)	197 €

	Width	Lower leg length
Size 0	28,4 cm	15 - 25 cm
Size 1	33,4 cm	19 - 32 cm
Size 1b	33,4 cm	27 - 49 cm
Size 2	40,3 cm	27 - 49 cm


Footplate standard – without fixed position, foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps

Footplate with lock – foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps

Footplate with angle adjustment
foldable, fixed mounted on the
seat unit, with non-slip surface
and 4 slots each for fixing straps.

Calf plant for Footplate

covered with black flexi-fabric




	for size 0	for size 1+1b	for size 2
Width	27 cm	32 cm	38 cm
Height	6,5 cm	6,5 cm	6,5 cm

Price / piece			Upcharge / piece			Upcharge / piece			Price / piece		
☐ for size 0	(7800021)	incl.	☐ for size 0	(7800022)	52 €	☐ for size 0	(7800023)	96 €	☐ for size 0	(7800024)	206 €
☐ for size 1	(7801021)	incl.	☐ for size 1	(7801022)	52 €	☐ for size 1	(7801023)	96 €	☐ for size 1 + 1b	(7801024)	206 €
☐ for size 1b	(7801221)	incl.	☐ for size 1b	(7801222)	52 €	☐ for size 1b	(7801223)	96 €	☐ for size 2	(7802024)	206 €
☐ for size 2	(7802021)	incl.	☐ for size 2	(7802022)	52 €	☐ for size 2	(7802023)	96 €			


individually adjustable in height, width and angle



Split footplates – foldable, fixed mounted on the footplate holder, with non-slip surface and each 4 slots for the for the fixation of straps



Foldable, fixed mounted
on the footplate holder,
with magnetic lock

[illegible]

Price / piece		Price / pair		Upcharge / pair	
☐ for size 0 (7800 025)	285 €	☐ Size 1 (7801 026)	incl.	☐ Size 1 (7801 027)	168 €
☐ for size 1 (7801 025)	285 €	☐ Size 2 (7802 026)	incl.	☐ Size 2 (7802 027)	168 €
☐ for size 1 b (7801 225)	285 €			☐ Size 3 (7803 027)	168 €
☐ for size 2 (7802 025)	285 €				

Accessories

Therapy tray – made of wood with edge, flip away



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Price / piece
☐ for size 0 (78 00 028)	249 €
☐ for size 1 (78 01 028)	249 €
☐ for size 2 (78 02 028)	249 €

Add-on option:

☐ Folding up protection

Therapy tray – transparent, with edge, flip away



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Price / piece
☐ for size 0 (78 00 029)	340 €
☐ for size 1 (78 01 029)	340 €
☐ for size 2 (78 02 029)	340 €

(78 01 036) 98 €

Therapy tray – made of wood with edge, prefixable, fixed on both sides



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Price / piece
☐ for size 0 (78 00 038)	265 €
☐ for size 1 (78 01 038)	265 €
☐ for size 2 (78 02 038)	265 €

Therapy tray – transparent, with edge, prefixable, fixed on both sides



	Width	Depth
for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm

	Price / piece
☐ for size 0 (78 00 039)	380 €
☐ for size 1 (78 01 039)	380 €
☐ for size 2 (78 02 039)	380 €

Accessories

4 Castors – 2 with locking devices



	Upcharge / set
☐ 75 mm (73 00 058)	incl.
☐ 100 mm (73 01 058)	58 €
☐ 125 mm (73 02 058)	92 €

4 Castors – all with locking devices



	Upcharge / set
☐ 75 mm (73 00 060)	35 €
☐ 100 mm (73 01 060)	92 €
☐ 125 mm (73 02 060)	127 €

Accessories

Central parking brake for rear castors – incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system



Conversion kit for combi subframe

complete with 12 1/2" wheels, screw-in axles, drum brake, steering wheels 6" and push bar



	Price / set
☐ for size 0+1 Parking lever mounted left (79 00 030)	425 €
☐ for size 0+1 Parking lever mounted right (79 01 030)	425 €
☐ for size 1b+2 Parking lever mounted left (79 02 030)	425 €
☐ for size 1b+2 Parking lever mounted right (79 03 030)	425 €
☐ for size 0+1 Parking lever mounted left – castors front with brakes (79 00 031)	459 €
☐ for size 0+1 Parking lever mounted right – castors front with brakes (79 01 031)	459 €
☐ for size 1b+2 Parking lever mounted left – castors front with brakes (79 02 031)	459 €
☐ for size 1b+2 Parking lever mounted right – castors front with brakes (79 03 031)	459 €

	Price / set
☐ for size 0+1 (74 01 001)	853 €
☐ for size 1b+2 (74 02 001)	853 €

Accessories

Hydraulic height adjustment

Retrofitting not possible! Not possible together with the conversion kit for combi subframe!



Electric height adjustment

by a 24 Volt electric motor; incl. battery and hand control



Additional gas dampers (2 pcs)

to support the lift function

(not possible in combination with the hydraulic or electric height adjustment!)

2x 300 N	recommended from 30 kg load
2x 500 N	recommended from 40 kg load

Anti tip – adjustable in height



	Upcharge / set
☐ for size 0+1 (74 01 050)	537 €
☐ for size 1b+2 (74 02 050)	537 €

	Price / set
☐ 74 00 102	939 €

	Price / set
☐ 2x 300 N (74 00 003)	119 €
☐ 2x 500 N (74 01 003)	119 €

	Price / piece
☐ for standard & electric height adjustment (74 00 004)	256 €
☐ for hydr. height adjustment (74 01 004)	256 €

Accessories

Wheel base extension insertable, inserts at front, length = 10cm 	Padded cover for footplate – cover- ed with black Flexi cover, height: 3 cm (not possible in combination with foot straps)	Padded cover for split footplate – covered with black Flexi cover, height: 3 cm (not possible in com- bination with foot straps)	Table upholstery – covered with black Flexi cover, height: 1 cm 
Price / pair	Price / piece	Price / pair	Price / piece
☐ 74 00 005 176	☐ for size 0 (78 00 032) 118 € ☐ for size 1 & 1b (78 01 032) 118 € ☐ for size 2 (78 02 032) 118 €	☐ for size 1 (78 01 033) 169 € ☐ for size 2 (78 02 033) 169 €	☐ for size 0 (73 00 054) 124 € ☐ for size 1 (73 01 054) 124 € ☐ for size 2 (73 02 054) 124 €

Accessories

Incontinence coating in the seat area 	Harness assembly kit for the seat consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)	Harness assembly kit for the back consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)
Price / piece	Upcharge / set	Upcharge / set
☐ 73 00 020 196 €	☐ 20 mm (73 00 025) 44 € ☐ 25 mm (73 01 025) 44 €	☐ 20 mm (73 00 026) 44 € ☐ 25 mm (73 01 026) 44 €


Harness

The **madita-fun.**-harness is assembled with the harness assembly sets (see page 22).

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 0	Size 1	Size 2	Size 3
A	13 cm	13 cm	18 cm	18 cm
B	17,5 cm	19 cm	20 cm	21 cm
C	15 cm	15 cm	19 cm	28 cm
D	3,5 cm	5 cm	5 cm	5 cm
E	25,5 cm	28 cm	34 cm	44 cm
Strap width	2 cm	2 cm	2,5 cm	2,5 cm




Price / piece

☐ Size 0 (73 00 032)	177 €
☐ Size 1 (73 01 032)	177 €
☐ Size 2 (73 02 032)	177 €
☐ Size 3 (73 03 032)	177 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 0	Size 1	Size 2
A	11 cm	11 cm	18 cm
B	21 cm	21 cm	22 cm
C	12 cm	12 cm	15 cm
D	2,5 cm	3,5 cm	3,5 cm
E	25,5 cm	28 cm	34 cm
Strap width	2 cm	2 cm	2,5 cm




Price / piece

☐ Size 0 (73 00 056)	206 €
☐ Size 1 (73 01 056)	206 €
☐ Size 2 (73 02 056)	206 €

Harness

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 0	Size 1	Size 2
A	23 cm	25 cm	28 cm
B	15 cm	17 cm	20 cm
C	33 cm	35 cm	38 cm
D	24 cm	26 cm	28 cm
Strap width	2 cm	2 cm	2,5 cm




Price / piece

☐ Size 0 (73 00 033)	148 €
☐ Size 1 (73 01 033)	148 €
☐ Size 2 (73 02 033)	148 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2,5 cm	2,5 cm



Price / piece

☐ Size 1 (73 01 034)	148 €
☐ Size 2 (73 02 034)	148 €
☐ Size 3 (73 03 034)	148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2,5 cm	2,5 cm

Price / piece

☐ Size 0 (73 00 035)	89 €
☐ Size 1 (73 01 035)	89 €
☐ Size 2 (73 02 035)	89 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2,5 cm	2,5 cm

Price / piece

☐ Size 0 (73 00 036)	155 €
☐ Size 1 (73 01 036)	155 €
☐ Size 2 (73 02 036)	155 €

Harness

Ankle cuffs – with velcro fastener and idler pulley, colour: black (not possible in combination with padded cover for footplate/split footplate)




	Ankle circumference	Strap width
Size 1	12 – 19 cm	25 mm
Size 2	15 – 23 cm	25 mm
Size 3	19 – 27 cm	40 mm

Price / pair

☐ Size 1 (73 01 061)	212 €
☐ Size 2 (73 02 061)	212 €
☐ Size 3 (73 03 061)	212 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



Klettstreifen zum Öffnen

	A	Strap width
Size 1	6,5 cm	2,5 cm
Size 2	8 cm	2,5 cm
Size 3	9,5 cm	2,5 cm

Price / pair

☐ Size 1 (78 01 034)	196 €
☐ Size 2 (78 02 034)	196 €
☐ Size 3 (78 03 034)	196 €

madita-fun mini.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The therapy chair for early intervention.



- Incl. armrests, footplate with slots for harness, upholstery for footplate, ergonomically shaped push handle, standard seat pad incl. cover,
- 4 Castors (75 mm) – 2 with locking devices
- Height adjustable per gas pressure spring (800 N)
- Frame colour: Telegrey 4 (white grey)
- Adjustment options:
 - Tilt-in-space (by gas pressure spring 250 N, (-)5° - 37° / 0° - 42°)
the trigger lever is mounted on the push handle
 - Back angle
 - Seat depth and seat width
 - Armrest, height, width and angle adjustable
 - Height adjustable backrest
 - Footplate, knee angle and height adjustable

HMV-No.: 26.11.05.0022



Dimensions

A	Seat depth	12 - 19 cm
B	Seat width	14 - 22 cm
C	Height short back	28 - 33 cm
C	Height long back	30 - 35 cm
D	Seat height	25 - 62 cm
E	Lower leg length	12 - 22 cm
F	Total width	52 cm
G	Length of base frame	58 cm
	Tilt-in-space	(-)5° - 37° / 0° - 42°
	Back angle	(-)5° - 25°
	Max. load	20 kg
	Weight	16 kg



Basic model
madita-fun mini.

Price / piece

☐ **3100 000** 2.410 €

Tilt-in-spacevoreinstellung

Please select a default setting ((-5° - 37° = standard configuration, if nothing else stated)

Price / piece

☐ (-5° - 37° **(3100 001)** incl.

Price / piece

☐ 0° - 42° **(3100 101)** incl.

Covers
Flamingo-Rosé
(No. 125)


Price

☐ **9100 125** incl.

Pazifik-Mint
(No. 126)



Price

☐ **9100 126** incl.

Cover according to customer request

Upcharge

☐ **9100 999** 180 €


Configuration options - back
Back incl. pad and cover
Standard


Height Width

Back (short)	28 - 33 cm	24 cm
Back (long)	30 - 35 cm	24 cm

Price / piece

☐ Back (short)	(3100 002)	185 €
☐ Back (long)	(3100 003)	185 €

Back with integrated upper arm guidance incl. pad and cover – with adjustable hol-low form in the lower trunk area through upholstery wedges (height approx. 3cm)


Height

Back (short)	28 - 33 cm
Back (long)	30 - 35 cm

Price / piece


☐ Back (short)	(3100 005)	208 €
☐ Back (long)	(3100 006)	208 €

Lateral pads – padded and covered, incl. bracket


Height	7 cm
Depth	9 cm
Chest width	14 - 20 cm

Price / pair

☐ **3100 009** 108 €

Lateral pads – long (trunk guidance) padded and covered, incl. bracket


Height	14 cm
Depth	6 cm
Trunk width	14 - 20 cm


Price / pair

☐ **3100 007** 170 €





Back without pad and cover – for the installation of custom cushions

Price reduction


☐ **3100 008** Net -52 €

Name embroidery on the back pad
please fill out the form on page 14!! An embroidery of Flexi fabrics is not possible.


Configuration options - seat

Seat pad standard, incl. cover 	Seat pad – with anatomic leg guidance , incl. cover 		Dysplasia seat – with a seat plate (up to 10 cm abductable) in the front sitting area, internal and external leg guidance (approx. 30° abduction per side), incl. pad and cover and a 10 cm wider adjustable footplate
Price / piece	Upcharge / piece		Upcharge / piece
☐ 3100011 incl.	☐ 3100012 218 €		☐ 3100013 324 €
Thigh guides – padded and covered, angle, width, depth adjustable (included in dysplasia seat) 	Internal guide – abductable (included in dysplasia seat) 	Incontinence coating in the seat area 	Seat without pad and cover for the installation of custom cushions
Price / pair	Price / set	Price / piece	Price reduction
☐ 13 cm (3100016) 195 € ☐ 16 cm (3100017) 195 €	☐ 3100015 208 €	☐ 3100014 196 €	☐ 3100018 Net -52 €

Accessories

Headrest (standard) padded, height, depth and angle adjustable 	Headrest (Shell-shaped) padded, height, depth and angle adjustable 	Headrest (Lateral guide, angle adjustable) padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable  A = 13,5 cm B = 10,5 cm C = 10,0 cm
Price / piece	Price / piece	Price / piece
☐ 3100021 259 €	☐ 3100124 366 €	☐ 3100026 588 €

Accessories

Therapy tray – made of wood with edge, pluggable 	Table upholstery covered with black Flexi cover, height: 1cm 	4 Castors – 2 with locking devices 	4 Castors – all with locking devices 
Price / piece	Price / piece	Upcharge / set	Upcharge / set
☐ 3100022 253 €	☐ 3100023 127 €	☐ 75 mm (7300058) incl. ☐ 100 mm (7301058) 58 €	☐ 75 mm (7300060) 35 € ☐ 100 mm (7301060) 92 €

madita hd.

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The strong therapy chair.

- Incl. padding (with increased safety) and cover
- Incl. armrests, footplate and push handle
- Frame colour: Anthracit metallic
- Reinforced base incl. 1x weight-bar (ca. 8 kg)
- 4 castors (75 mm) with locking devices
- Reinforced back angle adjustment
- Double-supported tilt
- Adjustment options:
 - Height adjustment via gas pressure spring (350 N)
 - Tilt-in-space
 - Back angle
 - Seat depth and width
 - Angle and width adjustable armrests
 - Knee angle and height adjustable footplate

HMV-No.: 26.11.05.0021

Basic model
madita hd.

Price

☐ **3002 100**

2.399 €

Covers

Black / Dark Grey
(No. 105 + 92)

Price

☐ **9100 205**

incl.


Dark Grey / Deep orange
(No. 92 + 106)

Price

☐ **9100 214**

incl.

Coated cover

Flexi Black / Light grey
(No. 43 + 50)


Upcharge

☐ **9100 226**

250 €

Dimensions

		madita hd.
A	Seat depth	30 – 42 cm
B	Seat width	28 – 38 cm
C	Back height with short back	42 – 52 cm
C	Back height with long back	53 – 63 cm
D	Seat height	47,5 – 58 cm
E	Lower leg length	30 – 48 cm
F	Total width	53 cm
G	Length of base frame	77 cm
	Tilt-in-space	(-)5° – 20°
	Back angle	(-)5° – 25°
	Max. load	65 kg
	Weight	29 kg




Headrest

Headrest

with reinforced bracket,
padded incl. cover,
height and
angle adjustable



Headrest (lateral guide, angle adjustable) – padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable



A = Width contact surface
B = Height headrest
C = Width lateral guide

	A	B	C
Size 1	13,5cm	10,5 cm	10 cm
Size 2	15 cm	10,5 cm	13,5 cm
Size 3	17 cm	12,5 cm	16,5 cm

Price / piece

☐ **3002 104**

Price / piece

383 €

☐ **Size 1 (3001030)**

594 €

☐ **Size 2 (3002030)**


594 €

☐ **Size 3 (3003030)**

594 €

Back

Back – reinforced, incl. pad and cover



short back	42 – 52 cm
long back	53 – 63 cm

Price / piece

☐ **short (3002 102)**

259 €

☐ **long (3002 103)**

259 €

Lateral pads

Lateral pads – with reinforced holder, padded, incl. cover



	Height	Depth	Chest width
Size 1	8 cm	13 cm	21 – 38 cm
Size 2	11 cm	16 cm	21 – 38 cm

Price / pair

☐ **Size 1 (3001025)**


174 €

☐ **Size 2 (3002025)**

174 €

Configuration options - seat

Standard seat pad – incl. cover
(standard configuration, if nothing
else stated)



Price / piece

☐ **3002 018**

incl.

Seat pad – incl. cover with anatomi-
c leg guidance



Upcharge / piece

☐ **3002 019**

219 €

Incontinence coating
in the seat area



Price / piece

☐ **3000 001**

198 €

Minor position


**Without seat and back upholstery
and cover** – for the installation of
custom cushions

Price reduction

☐ **7300 059**

Net -104 €

Seat angle adjustment		Thigh guides	
Trigger lever mounted directly on the gas spring (standard configuration, if nothing else stated)	Trigger lever mounted on the push bar	Thigh guides padded, incl. cover, adjustable angle, depth and width	
Price / piece	Price / piece	Price / pair	
☐ 73 00 085 incl.	☐ 73 01 085 incl.	☐ 19 cm (73 02 114) ☐ 24 cm (73 02 115)	197 € 197 €

Accessories			
Adjustable pelvic guidance – padded, incl. cover, seat width A = 26 - 36 cm  	Knee abduction wedge – pluggable, padded, incl. cover 	Therapy tray – made of wood with edge, prefixable, fixed on both sides 	
	Upcharge / piece	Price / piece	Price / piece
☐ 73 02 021	407 €	☐ 73 02 023	225 €
		☐ 30 02 023	315 €


Accessories			
Weight ca. 8 kg (max. 3 extra weights possible)		4 Castors – 100 mm, all with locking devices 	Wheel base extension Length = 10 cm, inserts at back
Price / piece	Quantity	Upcharge / set	Price / set
☐ 30 02 021	152 €	☐ 73 01 060	92 €
		☐ 73 01 028	168 €

Polster		Accessories	
Padded cover for footplate – cover- ed with black Flexi cover, height: 3 cm (not possible in combination with foot straps)		Padded tray cover – covered with black Flexi cover, height: 1 cm 	Harness assembly kit for the seat consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)
Price / piece	Price / piece	Upcharge / set	Upcharge / set
☐ 73 02 038	110 €	☐ 20 mm (73 00 025) 44 € ☐ 25 mm (73 01 025) 44 €	☐ 20 mm (73 00 026) 44 € ☐ 25 mm (73 01 026) 44 €

Harness

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2
A	25 cm	28 cm
B	17 cm	20 cm
C	35 cm	38 cm
D	26 cm	28 cm
Strap width	2 cm	2,5 cm




Price / piece

- ☐ Size 1 (**73 01 033**)
☐ Size 2 (**73 02 033**)

148 €
148 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2,5 cm	2,5 cm



Price / piece

- ☐ Size 1 (**73 01 034**)
☐ Size 2 (**73 02 034**)
☐ Size 3 (**73 03 034**)

148 €
148 €
148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece

- ☐ Size 1 (**73 01 035**)
☐ Size 2 (**73 02 035**)

89 €
89 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece


- ☐ Size 1 (**73 01 036**)
☐ Size 2 (**73 02 036**)

155 €
155 €

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 1	Size 2	Size 3
A	13 cm	18 cm	18 cm
B	19 cm	20 cm	21 cm
C	15 cm	19 cm	28 cm
D	5 cm	5 cm	5 cm
E	28 cm	34 cm	44 cm
Strap width	2 cm	2,5 cm	2,5 cm




Price / piece

- ☐ Size 1 (**73 01 032**)
☐ Size 2 (**73 02 032**)
☐ Size 3 (**73 03 032**)

177 €
177 €
177 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 1	Size 2
A	11 cm	18 cm
B	21 cm	22 cm
C	12 cm	15 cm
D	3,5 cm	3,5 cm
E	28 cm	34 cm
Strap width	2 cm	2,5 cm




Price / piece

- ☐ Size 1 (**73 01 056**)
☐ Size 2 (**73 02 056**)

206 €
206 €

Harness

Ankle cuffs – with velcro fastener and idler pulley colour: black (not possible in combination with padded cover for footplate)




	Ankle circumference	Strap width
Size 1	12 – 19 cm	25 mm
Size 2	15 – 23 cm	25 mm
Size 3	19 – 27 cm	40 mm

Price / pair

- ☐ Size 1 (**73 01 061**)
☐ Size 2 (**73 02 061**)
☐ Size 3 (**73 03 061**)

212 €
212 €
212 €

Foot positioner – with straps, incl. velcro fastener (not possible in combination with padded cover for footplate)



Velcro strip to open

	A	Strap width
Size 1	6,5 cm	2,5 cm
Size 2	8 cm	2,5 cm
Size 3	9,5 cm	2,5 cm

Price / pair

- ☐ Size 1 (**78 01 034**)
☐ Size 2 (**78 02 034**)
☐ Size 3 (**78 03 034**)

196 €
196 €
196 €

madita maxi.

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____


Commission:
The solid therapy chair.

- Seat pad with anatomic shaping and cover
- The seat height is adjustable using a crank, and the angular adjustment can be operated through two gas pressure springs à 200 N via a foot pedal
- Incl. push handle and armrests with adjustable height, width and angle with extra side support
- Incl. 4 dual castors (100 mm), 2 with locking devices


HMV-No.: 26.11.05.0021
Dimensions

A	Seat depth	36 - 47 cm
B	Seat width	34 - 44 cm
C*	Height short back	40 cm
C*	Height long back	48 cm
C**	Height short back	48 cm
C**	Height long back	56 cm
E	Seat height	48 - 60 cm
F	Lower leg length	37 - 49 cm
G	Total width	64 cm
H	Length of base frame	80 cm
	Tilt-in-space	(-) 5° - 30°
	Back angle	(-) 5° - 25°
	Max. load	120 kg
	Weight	26 kg

* made of 3 mm aluminium sheet metal | ** made of 4 mm aluminium sheet metal


Basic model
madita maxi.

Price / piece

☐ **30 03 000** 2.555 €

Frame colours
Anthrazit metallic


Price


☐ **92 00 014** incl.

Telegrey 4 (white grey)


Price

☐ **92 00 011** incl.

Covers		Coated cover		Minor position
	Petrol/Dark Grey (No. 112 + 92)		Black/Dark Grey (No. 105 + 92)	
	Price		Price	Upcharge
☐ 9100235	incl.	☐ 9100236	incl.	☐ 9100226 250 €
				☐ 3003001 Net -104 €

Configuration options - back				
	Standard-Back – incl. individually adjustable back pad, anatomic shape, back angle adjustment with mechanical locking option		Maxi back made of 4 mm aluminium sheet metal incl. individually adjustable back pad, anatomic shape, back angle adjustment with mechanical locking option	
	Short back	40 cm	Short back	48 cm
	Long back	48 cm	Long back	56 cm
	Price / piece		Price / piece	
	☐ short (30 02 002)	232 €	☐ short (30 03 002)	270 €
	☐ long (30 02 003)	232 €	☐ long (30 03 003)	270 €

Headrests			
Headrest (Standard) padded, height adjustable		Headrest (Lateral guide) padded, height and angle adjustment	
	Price / piece		Price / piece
☐ 30 03 004	283 €	☐ 30 02 011	377 €
		☐ 30 03 012	398 €

Accessories			
Reinforced lateral pads – padded		Thigh guides – padded and angle adjustable	
			
Price / pair		Price / pair	
☐ Chest width 24 – 32 cm (3003206)	142 €	☐ Standard 37 cm (3003021)	302 €
☐ Chest width 32 – 38 cm (3003006)	142 €	☐ Long 46 cm (3003022)	302 €
☐ Chest width 38 – 44 cm (3003106)	142 €		
		Knee abduction wedge pluggable and padded	
		Price / piece	
		☐ 3002008	186 €

Accessories			
Central parking brake for rear castors incl. 4 dual castors (100 mm), operated by using the parking lever which is mounted to the seat system		Therapy tray – made of wood with edge, pluggable, folding and removable	Incontinence coating in the seat area 
	Price / set	Price / piece	Price / piece
☐ Parking lever mounted left	(3003051) 440 €	☐ 3003023	☐ 3000001
☐ Parking lever mounted right	(3003551) 440 €	276 €	198 €
☐ Parking lever mounted left; Castors front with brakes	(3003151) 475 €		
☐ Parking lever mounted right; Castors front with brakes	(3003651) 475 €		

Accessories		Footplate	
4 Castors – 100 mm, all with locking devices		Wheel base extension – insertable, inserts at rear, length = 10 cm (not possible together with central parking brake)	
	Upcharge / set		Price / set
☐ 30 03 050	35 €	☐ 30 03 024	170 €
		Wooden footplate – height and foot angle adjustable, also folding	
			Price / piece
		☐ 30 03 015	440 €
		Plug-on footplate height approx. 8 cm above the floor	
			Price / piece
		☐ 30 03 016	357 €

Harnessmontagesätze

Harness assembly kit for the seat 25 mm – consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)

Harness assembly kit for the back 25 mm – consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

Upcharge / set

Upcharge / set

☐ 30 03 070

44 €


☐ 30 03 071

44 €

Harness


Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 3	Size 4
A	18 cm	20 cm
B	21 cm	22 cm
C	28 cm	30 cm
D	5 cm	6 cm
E	44 cm	52 cm
Strap width	2,5 cm	2,5 cm



Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 4
A	36 cm
B	26 cm
C	55 cm
Strap width	2,5 cm



Price / piece

Price / piece

☐ Size 3 (30 03 055)

177 €

☐ Size 4 (30 04 055)


177 €

☐ Size 4 (30 04 058)

148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 2	Size 3
Pad width	5 cm	6 cm
Pad length per side	28 cm	30 cm
Strap length per side	40 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece


☐ Size 2 (30 02 059)

89 €

☐ Size 3 (30 03 059)

89 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 2	Size 3
Pad width	5 cm	6 cm
Pad length per side	28 cm	30 cm
Strap length per side	40 cm	40 cm
Strap width	2,5 cm	2,5 cm

Price / piece

☐ Size 2 (30 02 060)

155 €

☐ Size 3 (30 03 060)

155 €



ATTENTION!

Available from **Q2** as version 2023!
Incl. quick height adjustment and
even more colour choice!
QR code to the current fitting sheet



smilla.

The convenient seating aid for everyday use.

You can find more information starting on page
8 or in the video by following the QR code:



sofie.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The individual seat unit that grows with the child.


- Incl. armrests with adjustable height, width and angle
- Incl. standard seat pad with cover
- Adjustment options: back angle, back height, seat depth, seat width
- Adaptable to all standard exterior and interior chassis

HMV-No.: 26.11.04.0xxx



Presettings

Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:



A Axilla height _____

B Chest width _____

C Back height _____

D Pelvic width _____


E Seat depth _____

F Lower leg length _____

☐ **17 00080**

Price

incl.

Dimensions						
		Size 0	Size 1	Size 1 b	Size 2	
A	Seat depth*	15 – 25 cm	21 – 30 cm	26 – 36 cm	30 – 42 cm	
B	Seat width	16 – 26 cm	20 – 30 cm	20 – 30 cm	28 – 38 cm	
C	Back height (short back)	35 – 43 cm	35 – 43 cm	42 – 52 cm	42 – 52 cm	
C	Back height (long back)	42 – 52 cm	42 – 52 cm	53 – 63 cm	53 – 63 cm	
	Back angle	(-)5° – 25°	(-)5° – 25°	(-)5° – 25°	(-)5° – 25°	
	Max. load	30 kg	30 kg	45 kg	50 kg	
	Weight	9 kg	10 kg	10 kg	12 kg	
	Recommended for subframes	freddy. Size 0	freddy. Size 0	freddy. size 1 + anti tip bar, franka. size 2, & fridolin maxi.	freddy. size 2 + anti tip bar, franka. size 2 & fridolin maxi.	

* reduced by 1 cm in combination with adjustable pelvic guide

Basic model			
sofie. size 0	sofie. size 1	sofie. size 1b	sofie. size 2
Price / piece	Price / piece	Price / piece	Price / piece
☐ 17 00 100 1.455 €	☐ 17 01 100 1.455 €	☐ 17 01 200 1.535 €	☐ 17 02 100 1.625 €

Frame colours					
Lasurit blue metallic 	White aluminium (silver) 	Dormant Red 	Pure white 	Anthracit metallic 	Coated with an individual colour
Price	Price	Price	Price	Price	Upcharge
☐ 92 00 010 incl.	☐ 92 00 002 incl.	☐ 92 00 003 incl.	☐ 92 00 100 incl.	☐ 92 00 014 incl.	☐ 17 00 090 180 €

Covers		
Dark Grey / Pink (No. 92 + 114) 	Dark Grey / Deep orange (No. 92 + 107) 	Dark Grey / Light blue (No. 92 + 106) 
Price	Price	Price
☐ 91 00 222 incl.	☐ 91 00 213 incl.	☐ 91 00 214 incl.


Covers	Name embroidery	Minor position
Black / Dark Grey (No. 105 + 92) 	Name embroidery on the back pad please fill out the form on page 141! An embroidery of Flexi fabrics is not possible. 	Without seat and back upholstery and cover for the installation of custom cushions.
Price		Price reduction
☐ 91 00 205 incl.		☐ 60 00 062 Net -104 €

Coated covers			
Flexi Black / Light grey (No. 43 + 50) 	Flexi Black / Red (No. 43 + 45) 	Flexi Blue (No. 44) 	Cover according to customer request
Upcharge	Upcharge	Upcharge	Upcharge
☐ 17 00 226 250 €	☐ 17 00 227 250 €	☐ 17 00 044 250 €	☐ 17 00 999 180 €

Headrests

Headrest (Standard) padded, incl. cover, height and angle adjustable	Headrest (Shell-shaped) padded, incl. cover, height and angle adjustable	Headrest (Lateral guide) padded, incl. cover, height and angle adjustable
		
Price / piece ☐ for size 0 (1700013) 254 € ☐ for size 1 (1701013) 254 € ☐ for size 2 (1702013) 254 €	Price / piece ☐ for size 0+1 (1701015) 360 € ☐ for size 2 (1702015) 360 €	Price / piece ☐ for size 0+1 (1701014) 340 € ☐ for size 2 (1702014) 340 €

Headrests

Headrest (Lateral guide angle adjustable) – padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable			Universal adapter for tailored headrests in the version 15 x 15 mm	
 A = Width contact surface B = Height headrest C = Width lateral guide		A	B	C
	Size 1	13,5 cm	10,5 cm	10 cm
	Size 2	15 cm	10,5 cm	13,5 cm
	Size 3	17 cm	12,5 cm	16,5 cm
	Price / piece			
☐ Size 1 (1701050)		588 €		☐ 1701016
☐ Size 2 (1702050)		588 €		85 €
☐ Size 3 (1703050)		588 €		

Configuration options - back (Standard-Ausführung oder mit integrierter Oberarmführung)

Back – incl. pad and cover	for size 0	for size 1	for size 1b	for size 2
short back: long back:	35 – 43 cm 42 – 52 cm	35 – 43 cm 42 – 52 cm	42 – 52 cm 53 – 63 cm	42 – 52 cm 53 – 63 cm
	Price / piece	Price / piece	Price / piece	Price / piece
 Short standard design: Long standard design:	☐ 1700003 185 € ☐ 1700004 185 €	☐ 1701003 185 € ☐ 1701004 185 €	☐ 1701103 185 € ☐ 1701104 185 €	☐ 1702003 208 € ☐ 1702004 208 €
 Short with upper arm guide: Long with upper arm guide:	☐ 1700005 295 € ☐ 1700006 295 €	☐ 1701005 295 € ☐ 1701006 295 €	☐ 1701105 295 € ☐ 1701106 295 €	☐ 1702005 314 € ☐ 1702006 314 €

Configuration options - backeinheit

Back unit – fixed (standard configuration = without push handles, if nothing else stated)			Back unit – removable			Adjustable limitation of the back angle – to determine the individual max. back angle	
							
Price / piece			Price / piece			Price / Set	
☐ for size 0+1 with push handles	(1701009)	incl.	☐ for size 0+1 with push handles	(1701011)	370 €	☐ 1700046	62 €
☐ for size 2 with push handles	(1702009)	incl.	☐ for size 2 with push handles	(1702011)	370 €		
☐ for size 0+1 with push loop	(1701007)	incl.					

Configuration possibilities - lateral pads (Step 1: Select the bracket)

Chest width – mounted on sofie. Size 0	Chest width – mounted on sofie. Size 1 + 1b	Chest width – mounted on sofie. Size 2	Bracket	Upcharge / pair
13 – 20 cm	16 – 23 cm	26 – 32 cm	☐ Fixed (1700017)	incl.
			☐ Flip away (1700018)	259 €
20 – 27 cm	23 – 30 cm	32 – 38 cm	☐ Fixed (1701017)	incl.
			☐ Flip away (1701018)	259 €
___ – ___ cm	___ – ___ cm	___ – ___ cm	☐ Fixed and tailored (1709017)	70 €
			☐ Flip away and tailored (1709018)	330 €

(Step 2: Select the lateral pads)

	Lateral pads Size 0 padded, incl. cover (Standard configuration for sofie. size 0, if nothing else stated) Height: 7 cm Depth: 9 cm	Lateral pads Size 1 padded, incl. cover (Standard configuration for sofie. size 1, if nothing else stated) Height: 8 cm Depth: 13 cm	Lateral pads Size 2 padded, incl. cover (Standard configuration for sofie. size 2, if nothing else stated) Height: 11 cm Depth: 16 cm
	Price / pair	Price / pair	Price / pair
	☐ 1700019 106 €	☐ 1701019 106 €	☐ 1702019 106 €

Thigh guides

	Thigh guides for size 0 padded, incl. cover, adjustable angle, depth and width 12 cm = Standard configuration for size 0	Thigh guides for size 1 + 1b padded, incl. cover, adjustable angle, depth and width 14 cm = Standard configuration for size 1 + 1b	Thigh guides for size 2 padded, incl. cover, adjustable angle, depth and width 19 cm = Standard configuration for size 2
	Price / pair	Price / pair	Price / pair
	☐ 12 cm (1700021) 192 € ☐ 14 cm (1700022) 192 €	☐ 14 cm (1701021) 192 € ☐ 19 cm (1701022) 192 €	☐ 19 cm (1702021) 197 € ☐ 24 cm (1702022) 197 €

Seat option

Standard seat pad – incl. cover (Standard configuration, if nothing else stated) 	Seat pad – incl. cover with anatomic leg guidance 	Adjustable pelvic guidance – padded, incl. cover  
Price / piece	Upcharge / piece	Price / piece
☐ for size 0 (1700023) incl. ☐ for size 1 (1701023) incl. ☐ for size 1 b (1701123) incl. ☐ for size 2 (1702023) incl.	☐ for size 0 (1700024) 217 € ☐ for size 1 (1701024) 217 € ☐ for size 1 b (1701124) 217 € ☐ for size 2 (1702024) 217 €	☐ for size 0 – Seat width A = 15 – 24 cm (1700026) 384 € ☐ for size 1 + 1 b – Seat width A = 18 – 28 cm (1701026) 384 € ☐ for size 2 – Seat width A = 26 – 36 cm (1702026) 407 €

Seat option

Abductable leg guidance – with leg length compensation, padded, incl. cover (already includes thigh guides)

	for size 0	for size 1+1b	for size 2
A	0° – 20°	0° – 20°	0° – 20°
B	7 – 10 cm	9 – 12 cm	10 – 14 cm
C	3 cm	3 cm	3 cm


Price / piece

☐ for size 0	(1700025)	520 €
☐ for size 1	(1701025)	520 €
☐ for size 1 b	(1701125)	520 €
☐ for size 2	(1702025)	540 €

A = Abduction angle

B = Thigh width

C = Leg length compensation



Footplate

Footplate – height and knee angle adjustable, foldable, fixed to the seat unit, with non-slip surface and 8 slots for attaching straps



Price / piece

☐ for size 0	(1700102)	339 €
☐ for size 1	(1701102)	339 €
☐ for size 1b	(1701112)	339 €
☐ for size 2	(1702102)	339 €

Knee abduction wedge – pluggable, padded, incl. cover



Price / piece

☐ Size 0	(1700108)	225 €
☐ Size 1	(1701108)	225 €
☐ Size 2	(1702108)	225 €

Accessories

Therapy tray – made of wood with edge



	Width	Depth
for size 0 + 1	50 cm	40 cm
for size 2	60 cm	50 cm

Price / piece

☐ for size 0	(1700027)	249 €
☐ for size 1	(1701027)	249 €
☐ for size 2	(1702027)	249 €

Therapy tray – transparent, with edge



	Width	Depth
for size 0 + 1	50 cm	40 cm
for size 2	60 cm	50 cm

Price / piece

☐ for size 0	(1700028)	340 €
☐ for size 1	(1701028)	340 €
☐ for size 2	(1702028)	340 €

Padding

Padded cover for footplate – covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)



Price / piece

☐ for size 0	(1700035)	110 €
☐ for size 1	(1701035)	110 €
☐ for size 2	(1702035)	110 €

Table upholstery – covered with black Flexi cover, height: 1 cm



Price / piece

☐ for size 0	(1700037)	124 €
☐ for size 1	(1701037)	124 €
☐ for size 2	(1702037)	124 €

Accessories

Incontinence coating – in the seat area



Price / piece

☐ 1700029	196 €
----------------------------------	-------

Harness assembly kit for the seat

consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)

Upcharge / set


☐ 20 mm	(1700033)	44 €
☐ 25 mm	(1701033)	44 €


Harness assembly kit for the back

consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled)

Upcharge / set

☐ 20 mm	(1700034)	44 €
☐ 25 mm	(1701034)	44 €

Adapter									
Trapezoidal adapter lower part			Trapezoidal adapter top part			Assembly of one of the adjacent adapters for seat depth			
	to suit			to suit		Price			
Size 0	sofie. size 0		Size 0	sofie. size 0		☐ Size 0:	15 – 21cm	(1700030)	incl.
Size 1	sofie. size 0-2		Size 1	sofie. size 0-2		☐ Size 0:	19 – 25 cm	(1700031)	incl.
Price / piece			Price / piece			☐ Size 1:	21 – 26 cm	(1701030)	incl.
☐ Size 0 (1700085)		199 €	☐ Size 0 (1700086)		130 €	☐ Size 1:	25 – 30 cm	(1701031)	incl.
☐ Size 1 (1701085)		199 €	☐ Size 1 (1701086)		130 €	☐ Size 2:	30 – 36 cm	(1702030)	incl.
						☐ Size 2:	37 – 42cm	(1702031)	incl.

Footplate for subframes					Calf plant				
					Calf plant for Footplate covered with black flexi-fabric 				
				Width	Lower leg length				
			Size 1	28,4 cm	15 – 25 cm				
			Size 2	33,4 cm	19 – 32 cm				
			Size 3	40,3 cm	27 – 49 cm				
Footplate standard – without fixed position, foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps			Footplate with lock – foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps		Footplate with angle adjustment foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps				
Price / piece			Price / piece		Upcharge / piece		Price / piece		
☐ Size 1 (1701090)		346 €	☐ Size 1 (1701091)	398 €	☐ Size 1 (1701092)	442 €	☐ Size 1 (1701093)	250 €	
☐ Size 2 (1702090)		346 €	☐ Size 2 (1702091)	398 €	☐ Size 2 (1702092)	442 €	☐ Size 2 (1702093)	250 €	
☐ Size 3 (1703090)		346 €	☐ Size 3 (1703091)	398 €	☐ Size 3 (1703092)	442 €	☐ Size 3 (1703093)	250 €	

Padding	
Padded cover for footplate (Footplate for subframes) covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)	
Price / piece	
☐ for size 0 (1700094)	118 €
☐ for size 1 (1701094)	118 €
☐ for size 2 (1702094)	118 €


Harness

The **sofie**-harness is assembled with the harness assembly sets (see page 40).

Harness

Chest-shoulder strap – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder; colour: black

	Size 0	Size 1	Size 2	Size 3
A	13 cm	13 cm	18 cm	18 cm
B	17,5 cm	19 cm	20 cm	21 cm
C	15 cm	15 cm	19 cm	28 cm
D	3,5 cm	5 cm	5 cm	5 cm
E	25,5 cm	28 cm	34 cm	44 cm
Strap width	2 cm	2 cm	2,5 cm	2,5 cm




Price / piece

☐ Size 0 (1700038)	177 €
☐ Size 1 (1701038)	177 €
☐ Size 2 (1702038)	177 €
☐ Size 3 (1703038)	177 €

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips; colour: black

	Size 0	Size 1	Size 2
A	11 cm	11 cm	18 cm
B	21 cm	21 cm	22 cm
C	12 cm	12 cm	15 cm
D	2,5 cm	3,5 cm	3,5 cm
E	25,5 cm	28 cm	34 cm
Strap width	2 cm	2 cm	2,5 cm




Price / piece

☐ Size 0 (1700039)	206 €
☐ Size 1 (1701039)	206 €
☐ Size 2 (1702039)	206 €

Harness

T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 0	Size 1	Size 2
A	23 cm	25 cm	28 cm
B	15 cm	17 cm	20 cm
C	33 cm	35 cm	38 cm
D	24 cm	26 cm	28 cm
Strap width	2 cm	2 cm	2,5 cm




Price / piece

☐ Size 0 (1700040)	148 €
☐ Size 1 (1701040)	148 €
☐ Size 2 (1702040)	148 €

Groin harness (standard) – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back; colour: black

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2,5 cm	2,5 cm



Price / piece

☐ Size 1 (1701041)	148 €
☐ Size 2 (1702041)	148 €
☐ Size 3 (1703041)	148 €

Harness

Pelvic harness – with pad and adjustable clip on both sides; colour: black



	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2,5 cm	2,5 cm

Price / piece

☐ Size 0 (1700042)	89 €
☐ Size 1 (1701042)	89 €
☐ Size 2 (1702042)	89 €

4-point pelvic harness – with pad and adjustable clip on both sides; colour: black




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2,5 cm	2,5 cm

Price / piece

☐ Size 0 (1700043)	155 €
☐ Size 1 (1701043)	155 €
☐ Size 2 (1702043)	155 €

Harness

Ankle cuffs – with velcro fastener and idler pulley colour: black (not possible in combination with padded cover for footplate/split footplate)




	Ankle circumference	Strap width
Size 1	12 – 19 cm	25 mm
Size 2	15 – 23 cm	25 mm
Size 3	19 – 27 cm	40 mm

Price / pair

☐ Size 1 (1701044)	212 €
☐ Size 2 (1702044)	212 €
☐ Size 3 (1703044)	212 €

Foot positioner – with straps, incl. velcro fastener, colour: black (not possible in combination with padded cover for footplate/split footplate)



Velcro strip to open

	A	Strap width
Size 1	6,5 cm	2,5 cm
Size 2	8 cm	2,5 cm
Size 3	9,5 cm	2,5 cm

Price / pair

☐ Size 1 (1701046)	196 €
☐ Size 2 (1702046)	196 €
☐ Size 3 (1703046)	196 €



hts.

The universal toilet seat system.

Further information can be found from page 132.



activity chair.

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The dynamic therapy chair.

- Seat and back unit incl. pad
- Back angle adjustment via gas pressure spring from (-)5° to 20°
- Back height and seat depth adjustment
- Incl. 2-point pelvic harness
- Tilt-in-space adjustment via gas pressure spring
 - for the standard base (-) 15° to 15°
 - for the Hi/Lo base (-) 15° to 25°
- Frame colour: Effect black

HMV-No.: 26.11.05.0xxx

Dimensions

		Size 0	Size 1	Size 2
Seat depth		20 – 30 cm	28 – 41 cm	38 – 51 cm
Seat width		30 cm	36 cm	43 cm
Seat width (with hip support pads)		18 – 23 cm	22 – 29 cm	28 – 36 cm
Chest width		17 – 29 cm	17 – 30 cm	24 – 36 cm
Back height		32 – 39 cm	39 – 48 cm	48 – 61 cm
Height headrest – seat		37 – 53 cm	44 – 62 cm	50 – 75 cm
Height Armlehne – seat		13 – 19 cm	18 – 27 cm	18 – 27 cm
Lower leg length		18 – 30 cm	20 – 38 cm	29 – 54 cm
Seat height above floor	with short / long feet	24 – 32 cm / 47 – 55 cm	32 – 42 cm / 48 – 58 cm	42 – 52 cm / 51 – 61 cm
	with short / long feet and dual castors	34 – 42 cm / 47 – 55 cm	41 – 51 cm / 48 – 58 cm	42 – 52 cm / 51 – 61 cm
	with Hi/Lo base	25,5 – 63,5 cm	30,5 – 58,5 cm	40,5 – 61 cm
Frame width	with short / long feet	53 / 58 cm	59 / 62 cm	68 / 68 cm
	with short / long feet and dual castors	57 / 60 cm	62 / 65 cm	71 / 71 cm
	with Hi/Lo base	66 cm	68,5 cm	73,5 cm
Tilt-in-space	with standard base / with Hi/Lo base	(-)15° – 15° / (-)15° – 25°		
Back angle		(-)5° – 20°		
Max. load		34 kg	68 kg	113 kg

Basic model

activity chair.

Seat unit (standard)

Tilt-in-space adjustment via gas pressure spring (-) 15° to 20°



	Price / piece
☐ Size 0 (72 00 200) (R821)	960 €
☐ Size 1 (72 01 200) (R841)	980 €
☐ Size 2 (72 02 200) (R861)	1.120 €

activity chair.

Seat unit (dynamic)

Back angle adjustment via gas pressure spring from (-)5° to 20°; a gas pressure spring enables self-generated movement to be accommodated - the spring option offers the user 10° of movement.



	Price / piece
☐ Size 0 (72 00 300) (R822)	1.095 €
☐ Size 1 (72 01 300) (R842)	1.120 €
☐ Size 2 (72 02 300) (R862)	1.295 €

Untergestelle

Standard base activity chair.

incl. height-adjustable, non-slip feet; seat angle adjustment by gas pressure spring from (-) 15° to 15°



	Price / piece
☐ Size 0 (72 00 400) (R823)	682 €
☐ Size 1 (72 01 400) (R843)	693 €
☐ Size 2 (72 02 400) (R873)	813 €

Dynamic base activity chair.

incl. height-adjustable, non-slip feet; seat angle adjustment by gas pressure spring from (-) 15° to 15°; a gas pressure spring enables self-generated movement to be accommodated - the spring option offers the user room to move.



	Price / piece
☐ Size 0 (72 00 500) (R824)	786 €
☐ Size 1 (72 01 500) (R844)	813 €
☐ Size 2 (72 02 500) (R874)	948 €

Hi/Lo base activity chair.

incl. 4 double castors (100 mm, braked, 1x with directional lock); height adjustment (by gas pressure spring for size 0 and hydraulically for sizes 1+2); seat angle adjustment by gas pressure spring from (-) 25° to 15°



	Price / piece
☐ Size 0 (72 00 600) (R833)	2.119 €
☐ Size 1 (72 01 600) (R853)	2.975 €
☐ Size 2 (72 02 600) (R883)	3.157 €

Pad colour

Blue



	Price
☐ for size 0 (72 05 201) (R816)	incl.
☐ for size 1 (72 00 201) (R836)	incl.
☐ for size 2 (72 02 201) (R876)	incl.

Red



	Price
☐ for size 0 (72 05 301) (R817)	incl.
☐ for size 1 (72 00 301) (R837)	incl.
☐ for size 2 (72 02 301) (R817)	incl.

Feet



Short feet – standard configuration, if nothing else stated

	Price / set
☐ for size 0 (72 04 004) (R803)	incl.
☐ for size 1 (72 00 004) (R893)	incl.
☐ for size 2 (72 05 004) (R808)	incl.

Long feet – for all sizes

	Price / set
☐ 72 01 004 (R894)	incl.

Feet



4 Dual castors (75 mm)

short version, with brakes, 1 x with direction lock

	Upcharge / set
☐ for size 0 (72 06 004) (R805)	367 €
☐ for size 1+2 (72 02 004) (R895)	367 €



4 Dual castors (75 mm)

long version, with brakes, 1 x with direction lock, for all sizes

	Upcharge / set
☐ 72 03 004 (R896)	367 €

Accessories

Headrest (Standard) – depth, angle and height adjustable



Price / piece

☐ **72 00 010 (R898)**

329 €

Headrest (Shell-shaped) – depth, angle and height adjustable



Price / piece

☐ **72 00 011 (R897)**

353 €

Headrest (Lateral guide) – depth, angle and height adjustable, wing pads can be adjusted from almost flat to 90°



Price / piece

☐ **72 00 018 (R881)**

453 €

Accessories

Lateral pads – height and width adjustable


Depth
Height
for size 0-2

13 cm

10 cm

for size 1-2

16,5 cm

13 cm

Price / pair

☐ for size 0-2 **(72 01 005) (R828)**

289 €

☐ for size 1-2 **(72 02 005) (R838)**

312 €

Lateral pads – height and width adjustable, with chest strap


Depth
Height
for size 0+1

13 cm

10 cm

for size 2

16,5 cm

13 cm

Price / set

☐ for size 0 **(72 00 006) (R819)**

367 €

☐ for size 1 **(72 01 006) (R829)**

367 €

☐ for size 2 **(72 02 006) (R839)**

393 €

Push handles – recommended for all **activity chair**. with dual castors and models with Hi/Lo base



Price / pair

☐ **72 00 012 (R878)**

234 €

Accessories

Arm rests – height adjustable


Width
Length
for size 0

6 cm

22 cm

for size 1+2

7,5 cm

30,5 cm

Price / pair

☐ for size 0 **(72 05 007) (R827)**

367 €

☐ for size 1+2 **(72 00 007) (R868)**

367 €

Hip support pads – padded and angle adjustable (only for use together with the armrests)


Height
Length
for size 0

8 cm

14 cm

for size 1+2

10 cm

23,5 cm

Price / pair

☐ for size 0 **(72 05 008) (R814)**

367 €

☐ for size 1+2 **(72 00 008) (R867)**

367 €

Therapy tray – only for use together with the armrests



Price / piece

☐ for size 0 **(72 00 009) (R667)**

432 €

☐ for size 1 **(72 01 009) (R668)**

432 €

☐ for size 2 **(72 02 009) (R669)**

432 €

Accessories

Backrest extension



Price / piece

☐ for size 0+1 **(72 01 016) (R849)**

186 €

☐ for size 2 **(72 02 016) (R856)**

186 €

Spinal support and tube edge



Price / set

☐ for size 0 **(72 00 017) (R886)**

161 €

☐ for size 1 **(72 01 017) (R887)**

161 €

☐ for size 2 **(72 02 017) (R888)**

161 €

Abduction wedge
depth adjustable


Width
Depth
Height
for size 0+1

7,5 cm

15 cm

9 cm

for size 2

10 cm

19 cm

10 cm

Price / piece

☐ for size 0+1 **(72 01 013) (R884)**

198 €

☐ for size 2 **(72 02 013) (R891)**

198 €

Accessories

Knee lateral guide



Price / pair

☐ 72 00 014 (R892)

289 €

Knee support pad

with key lock,
depth adjustable


Price / piece

☐ for size 0 (72 00 015) (R854)

325 €

☐ for size 1 (72 01 015) (R858)

325 €

☐ for size 2 (72 02 015) (R859)

353 €

Harness

Chest strap



Price / piece

☐ for size 0 (72 00 020) (R802)

143 €

☐ for size 1 (72 01 020) (R863)

143 €

☐ for size 2 (72 02 020) (R807)

143 €

Chest-shoulder strap



normal

thin

Price / piece

☐ normal for size 0 (72 00 021) (R834)

212 €

☐ normal for size 1 (72 01 021) (R864)

212 €

☐ normal for size 2 (72 02 021) (R889)

212 €

☐ dünn for size 1 (72 04 021) (R882)

212 €

☐ dünn for size 2 (72 03 021) (R826)

212 €

Thigh straps



Price / pair

☐ for size 0 (72 00 022) (R811)

143 €

☐ for size 1 (72 01 022) (R865)

143 €

☐ for size 2 (72 02 022) (R813)

143 €

Groin harness



Price / piece

☐ for size 0 (72 00 023) (R813)

250 €

☐ for size 1 (72 01 023) (R866)

250 €

☐ for size 2 (72 02 023) (R846)

250 €

Accessories

Footplate – height adjustable, can be folded up and swung away for transfer, recommended for models with Hi/Lo base.


for size 0

for size 1

for size 2

Lower leg length

23 – 30 cm

30 – 42 cm

38 – 51 cm

Price / piece

☐ for size 0 (72 00 029) (R847)

656 €

☐ for size 1 (72 01 029) (R848)

656 €

☐ for size 2 (72 02 029) (R857)

656 €

Ankle positioner – with socket lock, only for use together with the footplate


Price / pair

☐ for size 0+1 (72 01 026) (R899)

143 €

☐ for size 2 (72 02 026) (R855)

143 €

Accessories

Sandals with foot straps

only for use together
with the footplate


Price / pair

☐ for size 0 (72 00 027) (R944)

143 €

☐ for size 1 (72 01 027) (R954)

143 €

☐ for size 2 (72 02 027) (R964)

143 €

Wedges – only for use together with the footplate and sandals


Price / pair

☐ for size 0 (72 00 028) (R697)

113 €

☐ for size 1 (72 01 028) (R698)

113 €

☐ for size 2 (72 02 028) (R699)

113 €

tina.
Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____


Signature: _____

Commission:
The cradle subframe with cradle system.

- Electrical height adjustment via hand switch
- 12" inflated rear wheels with drum brake, 8" PU front wheels
- Universal mounting points for all types of cradles
- Frame colour: white aluminium (silver)


Dimensions - cradle subframe


		for three-piece cradle system	for one-piece cradle system
A	Height of cradle support from floor	53 - 79 cm	53 - 79 cm
B	Adjustable back angle	0° - 50°	-
C	Seat angle (adjustable in steps)	0° - 18°	-
D	Adjustable foot angle	0° - 50°	-
E	Inclination angle lying surface	-	0° - 20°
F	Receptive length for cradles	-	95 cm
	Base frame width	64 cm	64 cm
	Max. load	100 kg	100 kg
	Weight	44 kg	44 kg


Basic model cradle subframe

tina. (for three-piece cradle system) with electrical back adjustment, adjustable seat and foot angle		tina. (for one-piece cradle system) with electrical inclination angle adjustment of the lying surface	
Price / piece		Price / piece	
☐ 4000 101	5.440 €	☐ 4000 201	5.440 €

Dimensions - cradle system

		Size 2	Size 3	Size 4
A	Width	70 cm	70 cm	70 cm
B	Length of head section	60 cm	75 cm	90 cm
C	Centre section length	35 cm	35 cm	35 cm
D	Length of foot section	45 cm	55 cm	65 cm
E	Total length	140 cm	165 cm	190 cm
F	Height	12 cm	12 cm	12 cm



Cradle system

Head section



Note: Size 4 includes an extended push handle and anti-tip device!

	Length	Width
Size 2	60 cm	70 cm
Size 3	75 cm	70 cm
Size 4	90 cm	70 cm

i Custom modification is possible.

	Price / piece
☐ Size 2 (40 02 102)	395 €
☐ Size 3 (40 03 102)	395 €
☐ Size 4 (40 04 102)	838 €

Seat section



	Length	Width
Size 2-4	35 cm	70 cm

i Custom modification is possible.

	Price / piece
☐ Size 2-4 (40 02 103)	253 €

Leg section



	Length	Width
Size 2	45 cm	70 cm
Size 3	55 cm	70 cm
Size 4	65 cm	70 cm

i Custom modification is possible.

	Price / piece
☐ Size 2 (40 02 104)	395 €
☐ Size 3 (40 03 104)	395 €
☐ Size 4 (40 04 104)	395 €

Matratzen

Vacuum mattress – for **tina**. incl. valve and electric pump with two-way switch



	Length	Width
Size 2	ca. 150 cm	ca. 75 cm
Size 3	ca. 175 cm	ca. 75 cm
Size 4	ca. 200 cm	ca. 80 cm

Cover made of breathable fabric with insert made of spacer fabric (for adjacent vacuum mattress)

i Custom modification is possible!

	Price / piece		Price / piece
☐ Size 2 (41 02 202)	3.302 €	☐ for size 2 (41 02 201)	397 €
☐ Size 3 (41 03 202)	3.302 €	☐ for size 3 (41 03 201)	397 €
☐ Size 4 (41 04 202)	3.302 €	☐ for size 4 (41 04 201)	397 €

Support mattress (contourable with lateral guide)

Filling made of polystyrene beads and soft foam

Cover in bi-elastic fabric with polyurethane coating:

- waterproof
- breathable
- washable and wipe disinfectable
- anti-slip on the underside

	Length	Width
Size 2	ca. 150 cm	ca. 75 cm
Size 3	ca. 175 cm	ca. 75 cm
Size 4	ca. 200 cm	ca. 80 cm


	Price / piece
☐ Size 2 (41 02 207)	900 €
☐ Size 3 (41 03 207)	900 €
☐ Size 4 (41 04 207)	1.015 €

Accessories

Extended push handle – for head sections larger than 75 cm (already included with head section size 4)



Anti tip – pivoting, adjustable in height (already included with head section size 4)



Upcharge / piece

Price / piece

□ 40 00 120

164 €

□ 40 00 140

281 €

Accessories

Tray for medical devices – height adjustable, with slotted holes for strap (25mm) to secure the medical devices (The tray for medical devices is only available in conjunction with the extended push handle)



Width 43cm

Depth 40cm

Price / piece

□ 40 09 011

408 €

Infusion holder – telescopic



Price / piece

□ 40 00 010

261 €

Accessories

Preparation of the frame for mounting an E-Fixdrive 12"

Installation of the delivered E-Fix-drive

Large support cushion

Wide seat belt – strap with sliding unit to fix on the tube frame

Net price (not discountable)

Net price (not discountable)

Price / piece

Price / piece

□ 40 00 130

234 €

□ 40 00 131

177 €

□ 40 00 005

180 €

□ 40 00 009

253 €

new



todd.

Swivel arm pelvic frame.

Further information can be found from page 116.



Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The pivoting car seat.


- Seat and back plate made of aluminium with back angle adjustment
- Removable pads and fixed lateral pads
- Incl. height adjustable and padded headrest and 5-point brace straps
- Three-point strap guidance, frame colour: black – silky matt, RAL 9005
- Tested and approved based on ECE-44/03

HMV-No.: 26.11.06.0009

Dimensions

		Size 1		Size 2		Size 3	
		a	b	a	b	a	b
A	Seat depth	25 cm	30 cm	30 cm	35 cm	35 cm	40 cm
B	Seat width	28 cm		32 cm		36 cm	
C	Back height	40 cm		48 cm		56 cm	
D	Headrests height	59 - 69 cm		65 - 77 cm		76 - 86 cm	
E	Lateral pad height	25 - 33 cm		32 - 40 cm		36 - 44 cm	
	Tilt-in-space*	0° / 15°		0° / 10°		0° / 10°	
F	Chest width	14 - 25 cm		16 - 29 cm		16 - 32 cm	
G	Back angle	90° - 108°		90° - 108°		90° - 108°	
	Max. load	30 kg		45 kg		60 kg	

* When using the footplate and tilt adjustment



Basic model

sascha. size 1		sascha. size 2		sascha. size 3	
Price / piece		Price / piece		Price / piece	
☐ Size 1a (14 01 001)	1.995 €	☐ Size 2a (14 02 001)	2.070 €	☐ Size 3a (14 03 001)	2.115 €
☐ Size 1b (14 01 100)	1.995 €	☐ Size 2b (14 02 100)	2.070 €	☐ Size 3b (14 03 100)	2.115 €

Covers		
Fabric design Black / Red (No. 105 + 91) 	Fabric design Black / Light blue (No. 105 + 106) 	Fabric design Black / Spacer fabric (No. 105 + 47) 
Price	Price	Price
☐ 9100206 incl.	☐ 9100230 incl.	☐ 9100231 incl.

Lateral pads		
Standard setting range for fixed lateral pads The max. chest width reduces by 2 cm on each side when using the pluggable lateral pads For size 1 a / b Height from seat to armpit 25 – 33 cm Width of chest area 14 – 25 cm For size 2 a / b Height from seat to armpit 32 – 40 cm Width of chest area 16 – 29 cm For size 3 a / b Height from seat to armpit 36 – 44 cm Width of chest area 16 – 32 cm	Setting range for tailored lateral pads outside of the adjacent dimensionse Width of chest area _____ cm Width of chest area outside of the adjacent dimensions on request.	Lateral pads – pluggable, easier to get in and out, incl. cover  Upcharge / piece ☐ for size 1, mounted on right (1401007) 226 € ☐ for size 1, mounted on left (1401008) 226 € ☐ for size 2, mounted on right (1402007) 226 € ☐ for size 2, mounted on left (1402008) 226 € ☐ for size 3, mounted on right (1403007) 226 € ☐ for size 3, mounted on left (1403008) 226 €
	Upcharge	
	☐ 1400009 231 €	

Strap guide		Accessories	
Shoulder strap guide 1400002: Shoulder strap guide mounted on right (standard configuration, if nothing else stated) 1400003: Shoulder strap guide mounted on left 1400004: Shoulder strap guide mounted on left	Lateral guide and torso support pad as safety impact protection, incl. cover and fixing strap 	Abduction wedge – pluggable, with bracket and cover 	ISOFIX bracket 
Price / piece	Price / piece	Price / piece	Price / set
☐ rechts (1400002) incl. ☐ links (1400003) incl. ☐ beidseitig (1400004) 90 €	☐ for size 1 (1401005) 736 € ☐ for size 2 (1402005) 736 € ☐ for size 3 (1403005) 736 €	☐ for size 1 (1401006) 334 € ☐ for size 2 (1402006) 334 € ☐ for size 3 (1403006) 334 €	☐ for size 1 (1401010) 237 € ☐ for size 2 (1402010) 237 € ☐ for size 3 (1403010) 237 €

Pivoting device		
Pivoting device with tilt adjustment incl. footplate support and retainer bar   (Measurements in mm)	Right pivoting	Left pivoting
	Price / set	Price / set
	☐ X = 394 mm for size 1 (1401002) 458 € ☐ X = 454 mm for size 2 (1402002) 458 € ☐ X = 494 mm for size 3 (1403002) 458 €	☐ X = 394 mm for size 1 (1401102) 458 € ☐ X = 454 mm for size 2 (1402102) 458 € ☐ X = 494 mm for size 3 (1403102) 458 €

Accessories								
Tilt adjustment on base plate incl. footplate support 	Footplate – removable (only for use in combination with the pivoting device or tilt adjustment) 	Footplate – for installation in the footwell, adjustable in height and angle  <table border="1"> <tr> <th>Width</th><th>Height</th><th>Depth</th></tr> <tr> <td>36 cm</td><td>19–31 cm</td><td>26 cm</td></tr> </table>	Width	Height	Depth	36 cm	19–31 cm	26 cm
Width	Height	Depth						
36 cm	19–31 cm	26 cm						
Price / piece	Price / piece	Price / piece						
☐ for size 1 (1401003) 171 € ☐ for size 2 (1402003) 171 € ☐ for size 3 (1403003) 171 €	☐ for size 1 (1401004) 161 € ☐ for size 2 (1402004) 161 € ☐ for size 3 (1403004) 161 €	☐ 1604007 216 €						

54 **siggi combi isofix.**

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The special car seat.

- Based on the BeSafe iZi Combi X4 ISOfix and was verified and approved with special accessories according to ECE-44/04
- From 0 - 18 kg to be installed facing against the flow of traffic by way of ISOFIX
- From 9 - 18 kg to be installed with the flow of traffic
- With three different sitting / resting positions
- The patented belt guide is included in the basic model
- With serial 5-point belt incl. shoulder pad
- With integrated support leg, additional side impact protection SIP+ and chock
- Height of head support and belt can be adjusted easily with one move
- High quality expanded polystyrene (EPS) in head, torso and upper thigh areas ensures particularly high security in case of side impact
- Fabric design: Fresh Black Cap

new

now also available as a set!

incl. lateral pads,
abduction wedge & positioning jacket

HMV-No.: 26.11.06.0021





**0 to
18 kg**

Dimensions

A	Seat width front	32 cm
	Seat width rear	27 cm
B	Seat surface depth	25 cm
C	Height of headrest	48 - 60 cm
	Seat angle	facing the direction of travel 30° / 37,5° / 45° facing away from the direction of travel* 25° / 32,5° / 40°
D	Total width	46 cm
E	Total height	68 cm
	Group	0 - 18 kg
	Weight	16 kg

* at a vehicle seat angle of 5°



Basic model

	Price / piece
☐ 16 02 300	1.621 €
☐ 16 02 500 Set (incl. lateral pads, abduction wedge & positioning jacket)	1.899 €

Accessories

Lateral pads – adjustable in height and width, for chest width 17 – 21 cm



Price / pair

☐ **16 00 001** 212 €

Abductor wedge – hook and loop fastener



Price / piece

☐ **16 00 002** 100 €

Arm rest cushion



Price / piece

☐ **16 01 002** 234 €

Accessories

Positioning jacket – hook and loop fastener



Price / piece

☐ **16 01 004** 89 €

Baby cushion



Price / piece

☐ **16 01 006** 89 €

Sun protection



Price / piece

☐ **16 01 005** 95 €

siggi 3 fix.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The car seat that grows with you.

- Based on the BeSafe iZi Up X3 FIX and was verified and approved with special accessories according to ECE-44/04
- With ISOFIX attachment points for installation in the direction of travel (use in case of non-existent bracket for ISOFIX is also permitted)
- With innovative side impact rotation (SIR) when installed with ISOFIX
- Height of head support and belt can be adjusted easily with one move
- Automatic adjustment of the shoulder belt guide to the back height
- Headrest padding with the possibility of customization
- Adjustable rest position, back angle adjustment by 20°
- Fabric design: Fresh Black Cap

new

now also available as a set!
incl. lateral pads & abduction wedge

HMV-No.: 26.11.06.0023


ECE R44/04
Special Needs
15-36 kg
E4
0444 3808S

**15 to
36 kg**



Dimensions

A	Seat width at front	32 cm
B	Seat width at rear	27 cm
C	Seat surface depth	25 cm
D	Height of headrest	55 - 75 cm
E	Total width	53 cm
F	Total height	68 - 87 cm
	Back angle	90° / 110°
	Group	15 - 36 kg
	Weight	8 kg



Basic model

	Price / piece
☐ 16 04 100	1.320 €
☐ 16 04 500 Set (incl. lateral pads & abduction wedge)	1.499 €

Accessories

Lateral pads – adjustable in height and width, for chest width 23 – 28 cm



Price / pair

☐ **16 04 005** 212 €

Abductor wedge – hook and loop fastener



Price / piece

☐ **16 00 002** 100 €

Lateral guide and torso support pad – as safety impact protection, incl. cover and fixing strap



Price / piece

☐ **16 04 006** 396 €

Accessories

Positioning jacket – with zipper



Price / piece

☐ **16 04 003** 256 €

Pelvic positioning harness



Price / set

☐ **16 04 004** 256 €

Footplate – for installation in the footwell, adjustable in height and angle



Width	Height	Depth
36 cm	19 - 31 cm	26 cm

Price / piece

☐ **16 04 007** 216 €

freddy.
Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The compact seat pan subframe.

- Tilt adjustment of the seat pan via gas pressure spring 300 N
- Adjustable seat depth using depth adjustable back unit
- Back and push grips, angle and height adjustable
- Incl. support head as a base plate for trapezoidal and external adapter, height and seat depth adjustment
- Incl. socket for footplate
- Incl. rear wheels with quick-release axle and drum brake

HMV-No.: 26.99.01.1033
(with 16" rear wheels)

HMV-No.: 26.99.01.1079
(with 24" rear wheels and grip ring)

Basic model


freddy. size 0	freddy. size 1	freddy. size 2	freddy. size 3
Price / piece	Price / piece	Price / piece	Price / piece
☐ 25 00 000 2.399 €	☐ 25 01 000 2.399 €	☐ 25 02 000 2.399 €	☐ 25 03 000 2.399 €

Frame colours

Turquoise blue	Lasurit blue metallic	White aluminium (silver)	Dormant Red	Pure white	Anthracite metallic	Black satin finish
Price	Price	Price	Price	Price	Price	Price
☐ 92 00 016 incl.	☐ 92 00 010 incl.	☐ 92 00 002 incl.	☐ 92 00 003 incl.	☐ 92 00 100 incl.	☐ 92 00 017 incl.	☐ 92 00 005 incl.

Dimensions

		Size 0	Size 1	Size 2	Size 3
A	Seat pan support height	16"	45 - 55 cm (47 - 57 cm*)		
		22"	47 - 57 cm (49 - 59 cm*)		
		24"	50 - 60 cm (52 - 62 cm*)		
B	Short seat depth	26 - 28,5 cm	29 - 36,5 cm	-	-
B	Standard Seat depth	30,5 - 38 cm	33 - 40,5 cm	42,5 - 50 cm	42,5 - 50 cm
	Tilt-in-space	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°
C	Width between rear wheels (with clothing guard) = Max. outer width of seat pan	39 cm	44 cm	49 cm	54 cm
D	Frame width	inside	28 cm	33 cm	38 cm
		outside	35 cm	40 cm	45 cm
E	Seat frame width	25 cm	30 cm	35 cm	35 cm
F**	Total length	63 - 79 cm	63 - 79 cm	66 - 82 cm	71 - 87 cm
G	Total width	58 cm	63 cm	68 cm	73 cm
	Max. load	90 kg	90 kg	120 kg	120 kg
	Weight	from 14 kg	from 14 kg	from 14 kg	from 14 kg





* incl. trapezoidal adapter ** depending on the wheel size

Variante



Push handles / bar

Basic frame design with additional gas pressure spring (300 N) – for sizes 2+3 (recommended from a patient weight of 90 kg incl. seat pan), retrofitting not possible		Standard push handles 	Push bar 
Upcharge / piece		Price / pair	Upcharge / piece
☐ for size 2 (25 02 020)	470 €	☐ 25 00 001 incl.	☐ for size 0 (25 00 002) 166 €
☐ for size 3 (25 03 020)	470 €		☐ for size 1 (25 01 002) 166 €
			☐ for size 2+3 (25 02 002) 166 €

Back

Back angle adjustment via snap joint (Standard version, if nothing else selected) 	Centralised fast mechanical adjustment mechanical 	Centralised fast mechanical adjustment via gas pressure spring (300 N) 
<i>Please note that one of these options must be selected.</i>		
Price / piece	Price / piece	Upcharge / set
☐ for size 0 (25 00 003) incl.	☐ for size 0 (25 00 004) incl.	☐ for size 0 (25 00 005) 104 €
☐ for size 1 (25 01 003) incl.	☐ for size 1 (25 01 004) incl.	☐ for size 1 (25 01 005) 104 €
☐ for size 2+3 (25 02 003) incl.	☐ for size 2+3 (25 02 004) incl.	☐ for size 2+3 (25 02 005) 104 €

Accessories

Back rest support of seat unit 	Armrests – mounted to the back tubes, adjustable width and angle 	Headrest bracket – to mount square tubes (15x15 mm) 
Price / set	Price / set	Price / piece
☐ for size 0 (25 00 007) 247 €	☐ for size 0 (25 00 025) 347 €	☐ 25 00 019 106 €
☐ for size 1 (25 01 007) 247 €	☐ for size 1 (25 01 025) 347 €	
☐ for size 2+3 (25 02 007) 247 €	☐ for size 2+3 (25 02 025) 347 €	


Front wheels

6" front wheels with rubber coat 150 x 30 mm, colour: black (standard configuration, if nothing else stated) 	Full rubber front wheels 140 x 40 mm, colour: black 	7" front wheels with rubber coat 180 x 45 mm, colour: black 
Price / pair	Upcharge / pair	Upcharge / pair
☐ 25 00 022 incl.	☐ 25 01 022 47 €	☐ Inflated tyres (25 01 021) 47 €

Rear wheels

Rear wheels inflated tyres with quick-release axle and drum brake, colour: black 	Puncture-resistant rear wheels with quick-release axle and drum brake, colour: black 	Clothing guard – fixed 	Clothing guard – removable 
Upcharge / set	Upcharge / set	Price / pair	Price / pair
☐ 16" (25 00 008) incl.	☐ 16" (25 00 009) incl.	☐ 16" (25 00 011) 253 €	☐ 16" (25 00 018) 351 €
☐ 22" (25 02 008) incl.	☐ 22" (25 02 009) incl.	☐ 22" (25 02 011) 253 €	☐ 22" (25 02 018) 351 €
☐ 24" (25 03 008) 116 €	☐ 24" (25 03 009) 116 €	☐ 24" (25 03 011) 253 €	☐ 24" (25 03 018) 351 €

Rear wheels

Spoke guard – for rear wheels, transparent 	Wheel base extension (not possible in combination with the bracket for Viamobil / Scalamobil and S-Max) 	Knee lever brake with bracket 
Price / pair	Price / pair	Price / pair
☐ 22" (25 02 010) 201 €	☐ 25 00 014 144 €	☐ 16" (25 00 015) 259 €
☐ 24" (25 03 010) 201 €		☐ 22" (25 02 015) 259 €
		☐ 24" (25 03 015) 259 €

Accessories

Anti tip bar – adjustable in height and length, also pivoting (not possible in combination with the bracket for Viamobil / Scalamobil and S-Max) 	Tilting device (not possible in combination with the bracket for Viamobil / Scalamobil and S-Max) 	Tilting device for dual anti tip bar mounted on the yoke (not possible in combination with the bracket for Viamobil / Scalamobil and S-Max)	Tilting device for dual anti tip bar and wheel base extension (not possible in combination with the bracket for Viamobil / Scalamobil and S-Max)
Price / piece	Price / piece	Price / piece	Price / piece
☐ Right attachment (25 01 013) 186 €	☐ Right attachment (25 01 012) 82 €	☐ Right attachment (25 03 012) 82 €	☐ Right attachment (25 05 012) 82 €
☐ Left attachment (25 02 013) 186 €	☐ Left attachment (25 02 012) 82 €	☐ Left attachment (25 04 012) 82 €	☐ Left attachment (25 06 012) 82 €


Adapter

Support head as a base plate for trapezoidal and external adapters			Trapezadapter	Lower part	Top part		
	short seat depth	stand. seat depth		Length	Front width	Rear width	
			Lower part size 0	20 cm	21 cm	14,5 cm	
for size 0	23 – 33 cm	28 – 38 cm	Lower part size 1	26 cm	26 cm	17,5 cm	
for size 1	27 – 37 cm	32 – 42 cm	Top part size 0	19 cm	16,5 cm	9,8 cm	
for size 2+3	–	40 – 50 cm	Top part size 1	25 cm	21,5 cm	12,7 cm	
			Price / piece		Price / piece		
☐ for size 0, short	(25 00 116)		incl.	☐ Lower part size 0	(25 00 034)		196 €
☐ for size 0, standard	(25 00 016)		incl.	☐ Lower part size 1	(25 01 034)		196 €
☐ for size 1, short	(25 01 116)		incl.	☐ Top part size 0	(25 00 035)		128 €
☐ for size 1, standard	(25 01 016)		incl.	☐ Top part size 1	(25 01 035)		128 €
☐ for size 2+3, standard	(25 02 016)		incl.				

Retainer systems

Securing eyes to transport the chassis in the car incl. mounting to the chassis, universally combinable e.g. with systems of Q-Straint, etc. (not compatible with E-Fix, tested according to ISO 7176/19)	Power junction system to transport the chassis by car with heavy-duty eyelets, as assembly kit, incl. red pelvic harness (freddy. is tested with the power junction system according to ISO 7176/19 and the system itself according to DIN 75078-2)	Bracket for Viamobil / Scalamobil as an assembly set – only possible for freddy. sizes 2+3 – only possible with use of „back angle adjustment by snap joint“ – please use accessories from Alber (anti tip bar, tilting device)	Bracket for S-Max as assembly set – only possible for freddy. sizes 2+3 – only possible with use of „back angle adjustment by snap joint“ – please use accessories from AAT (anti tip bar, tilting device)	Power junction system assembly at Schuchmann
Net price (not discountable) / set	Net price (not discountable) / set	Net price (not discountable) / set	Net price (not discountable) / set	Net price (not discountable)
☐ for size 0 (25 00 033) 193 €	☐ für 22+24" (25 02 030) 396 €	☐ 25 00 031 310 €	☐ 25 00 032 310 €	☐ 25 01 030 151 €
☐ for size 1-3 (25 01 033) 193 €	☐ für 16" (25 03 030) 396 €			

Footplate

							
						Calf plant for Footplate covered with black flexi-fabric	

Split footplate

Split footplates – adjustable in height and width, swing away, knee angle adjustable, with calf plates			Polster		
			Padded cover for footplate covered with black Flexi cover (not possible for split footplates)		
			Size 1	Size 2	
			Lower leg length (base plate to footrest)	16–28 cm	26–38 cm
			recommended for lower leg length	24–36 cm	34–36 cm
			Min. width	30 cm	36 cm
			Adjustable angle	100° – 165°	
			Price / set		Price / piece
			☐ Size 1 (93 01 002) 520 €		☐ Size 0 (25 00 044) 118 €
			☐ Size 2 (93 02 002) 520 €		☐ Size 1 (25 01 044) 118 €
					☐ Size 2 (25 02 044) 118 €

fridolin maxi.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The functional subframe for indoors.

- The seat height is adjustable using a crank
- The angular adjustment of the seat unit can be operated through two gas pressure springs à 200 N via a foot pedal
- Incl. 4 dual castors (100 mm), 2 with locking devices
- Incl. support head and plug-in footplate
- Frame colour: anthracite metallic
- Max. load: 100 kg

HMV-No.: 26.99.01.0032




Dimensions

A	Height of seat pan support	44 - 56 cm (46 - 58 cm*)
B	Tilt-in-space	(-) 5° - 30°
C	Total length	80 cm
D	Total width	66 cm
	Strength of gas pressure spring for tilt adjustment	2x200 N
	Max. load	100 kg
	Weight	11,5 kg

* incl. trapezoidal adapter

Basic model



	Price / piece
☐ 22 03 000	1.280 €






Accessories

Wheel base extension – insertable, inserts at rear, length = 10 cm (not possible together with central parking brake)	4 Castors – all with locking devices	Central parking brake for rear castors – incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system
		
Price / set	Upcharge / set	Upcharge / set
☐ 22 03 006 170 €	☐ 22 03 007 34 €	☐ standard (22 03 008) 440 € ☐ castors front with brakes (22 03 009) 475 €

Adapter

Maxi support head – for trapezoidal adapter size 1 and external adapter	Trapezoidal adapter lower part size 1 – for assembly on the subframe	Trapezoidal adapter top part size 1 – for assembly on the seat pan
		
Price / piece	Price / piece	Price / piece
☐ 22 03 003 incl.	☐ 22 01 034 201 €	☐ 22 01 035 131 €

Footplate

			Calf plant Calf plant for Footplate covered with black flexi-fabric
Footplate standard – without fixed position, foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps	Footplate with lock – foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps	Footplate with angle adjustment foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps	
Price / piece	Price / piece	Upcharge / piece	Price / piece
☐ Size 1 (22 01 040) 346 € ☐ Size 2 (22 02 040) 346 € ☐ Size 3 (22 03 040) 346 €	☐ Size 1 (22 01 041) 398 € ☐ Size 2 (22 02 041) 398 € ☐ Size 3 (22 03 041) 398 €	☐ Size 1 (22 01 042) 442 € ☐ Size 2 (22 02 042) 442 € ☐ Size 3 (22 03 042) 442 €	☐ Size 1 (22 01 043) 250 € ☐ Size 2 (22 02 043) 250 € ☐ Size 3 (22 03 043) 250 €

Pad

Push handle

Footplate

Padded cover for footplate – covered with black Flexi cover, Height: 3 cm (not for plug-on footplate and not possible in combination with foot straps)	Maxi push handle height and angle adjustable	Plug-on footplate – height approx. 8 cm above the floor
		
Price / piece	Price / piece	Price / set
☐ Size 1 (22 01 044) 118 € ☐ Size 2 (22 02 044) 118 € ☐ Size 3 (22 03 044) 118 €	☐ 22 03 002 340 €	☐ 22 03 004 349 €

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The universal subframe.


- Tilt adjustment of the seat pan via gas pressure spring 500 N
- Height adjustment via gas pressure spring 1200 N (800 N for size 0)
- 4 castors (75 mm), 2 with locking devices
- incl. mounting head for trapezoidal and third-party adapters and footrest insertion

HMV-No.: 26.99.01.0022

Dimensions

		Size 0	Size 1	Size 2
A*	Height of seat pan support	18 - 55 cm (20 - 57 cm**)		
	Height of seat pan support with electric height adjustment	26 - 55 cm (28 - 57 cm**)		
	Height of seat pan support with hydraulic height adjustment	20 - 56 cm (22 - 58 cm**)		
A	Height of seat pan support with combi subframe	-	25 - 64 cm (27 - 66 cm**)	
B	Total length	58 cm	61 cm	75 cm
	Total length with combi subframe	71 cm	76 cm	86 cm
C	Total width	52 cm	54 cm	64 cm
	Total width with combi subframe	62 cm	66 cm	76 cm
D	Frame height	27 cm	27 cm	27 cm
E	Tilt-in-space	(-) 5° - 35°	(-) 5° - 35°	(-) 5° - 35°
	Max. load	30 kg	50 kg	75 kg
	Weight	8 kg	10 kg	12 kg



* with 75 mm castors ** incl. trapezoidal adapter



Basic model		
franka. size 0	franka. size 1	franka. size 2
Price / piece	Price / piece	Price / piece
☐ 2100000 1.220 €	☐ 2101000 1.330 €	☐ 2102000 1.390 €




Frame colours					
Lasurit blue metallic	White aluminium (silver) Subframes in this colour available from stock	Dormant Red	Pure white	Telegrey 4	Anthracite grey Subframes in this colour available from stock
					
Price	Price	Price	Price	Price	Price
☐ 9200010 incl.	☐ 9200002 incl.	☐ 9200003 incl.	☐ 9200100 incl.	☐ 9200011 incl.	☐ 9200017 incl.

Push bar / handle		Accessories	
Push bar – angle adjustable, mounted on the support head 	Push handle – height and angle adjustment (mounted on the base frame) 	Anti tip – adjustable in height 	Wheel base extension – pluggable, mounted in the front, length = 10 cm 
Price / piece	Price / piece	Price / piece	Price / set
☐ for support head size 0 (2100013) 242 €	☐ for size 0+1 (2100001) 329 €	☐ for stand. & electric height adjustment (2100009) 256 €	☐ 2100209 176 €
☐ for support head size 1 (2101013) 242 €	☐ for size 2 (2100201) 329 €	☐ for hydr. height adjustment (2100109) 256 €	




Seat angle adjustment			
Trigger lever – mounted directly on the gas pressure spring (standard configuration, if nothing else stated) 	Trigger lever – mounted on the push bar/handle 	Trigger lever – with fixing plate for assembly to the back of the seat pan 	Safety bolt for Trigger lever, to prevent unintentional release of the gas pressure spring 
Price / piece	Upcharge / piece	Upcharge / piece	Price / piece
☐ 2100002 incl.	☐ 2100014 35 €	☐ 2100015 57 €	☐ 2101014 75 €

Accessories			
Conversion kit for combi subframe with 12 1/2" wheels, screw-in axles, drum brake, steering wheels 6" and Push handle 	Hydraulic height adjustment retrofitting not possible! (not possible together with the conversion kit for combi subframe or electric height adjustment) 	Electric height adjustment by a 24-volt electric motor; incl. battery and hand control unit 	Additional gas dampers (2 pcs) to support the lift function (not possible in combination with the hydraulic or electric height adjustment!)
Price / set	Upcharge / set	Price / set	Price / set
☐ for size 0 (2100101) 853 €	☐ for size 1 (2101050) 537 €	☐ 2100107 939 €	☐ 2x 300 N Recommended from 30 kg load
☐ for size 1 (2101001) 853 €	☐ for size 2 (2102050) 537 €		☐ 2x 500 N Recommended from 430 kg load
☐ for size 2 (2102001) 853 €			☐ 2x 300 N (2100005) 119 €
			☐ 2x 500 N (2100006) 119 €


Adapter

Support head size 0 – for trapezoidal adapter size 0 and external adapter		Trapezoidal adapter lower part size 0 – for assembly on the subframe			Trapezoidal adapter top part size 0 – for assembly on the seat pan					
										
			Length	Front width	Rear width		Length	Front width	Rear width	
Size 0			20 cm	21cm	14,5 cm		Size 0	19 cm	16,5 cm	9,8 cm
Price / piece		Price / piece			Price / piece					
☐ 2102003 incl.		☐ 2100016 200 €			☐ 2100017 130 €					

Adapter

Support head size 1 – for trapezoidal adapter size 1 and external adapter		Trapezoidal adapter lower part size 1 – for assembly on the subframe			Trapezoidal adapter top part size 1 – for assembly on the seat pan					
										
			Length	Front width	Rear width			Length	Front width	Rear width
Size 1			26cm	26cm	17,5cm	Size 1		25cm	21,5cm	12,7cm
Price / piece					Price / piece				Price / piece	
☐ 2100003		incl.			☐ 2101016				☐ 2101017	
					200 €				130 €	

Footplate

				Calf plant for Footplate covered with black flexi-fabric	
					
			Width	Lower leg length	
		Size 1	28,4 cm	15 – 25 cm	
		Size 2	33,4 cm	19 – 32 cm	
		Size 3	40,3 cm	27 – 49 cm	
Footplate standard – without fixed position, foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps		Footplate with lock – foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps		Footplate with angle adjustment foldable, fixed mounted on the seat unit, with non-slip surface and 4 slots each for fixing straps	
Price / piece		Price / piece		Upcharge / piece	
☐ Size 1 (2101023)	346 €	☐ Size 1 (2101024)	398 €	☐ Size 1 (2101025)	442 €
☐ Size 2 (2102023)	346 €	☐ Size 2 (2102024)	398 €	☐ Size 2 (2102025)	442 €
☐ Size 3 (2103023)	346 €	☐ Size 3 (2103024)	398 €	☐ Size 3 (2103025)	442 €
				Price / piece	
				☐ Size 1 (2101026)	250 €
				☐ Size 2 (2102026)	250 €
				☐ Size 3 (2103026)	250 €

Padded cover for footplate

Padded cover for footplate – covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)		Central parking brake for rear castors incl. 4 dual castors (100 mm), operated by using the brake lever which is mounted to the seat system		Central parking brake for rear castors incl. 4 dual castors (100 mm), front with locking devices , operated by using the brake lever which is mounted to the seat system	
Price / piece		Price / set		Price / set	
☐ Size 1 (2100027)		☐ for size 0 (2100011)		☐ for size 0 (2100111)	
118 €		424 €		459 €	
☐ Size 2 (2101027)		☐ for size 1 (2101011)		☐ for size 1 (2101111)	
118 €		424 €		459 €	
☐ Size 3 (2102027)		☐ for size 2 (2102011)		☐ for size 2 (2102111)	
118 €		424 €		459 €	

Castors

		4 Castors – 2 with locking devices (125 mm not possible for size 0)		4 Castors – all with locking devices (125 mm not possible for size 0)	
		Upcharge / set		Upcharge / set	
☐ 75 mm (2104075)		incl.		☐ 75 mm (2106075)	
☐ 100 mm (2104100)		58 €		☐ 100 mm (2106100)	
☐ 125 mm (2104125)		92 €		☐ 125 mm (2106125)	
				127 €	

Fabric combinations.

for our seating systems **madita.** & **madita-fun.**



Pink 2
(No. 115)



Deep Sea Blue
(No. 111)



Leaf-Green
(No. 123)



Slate grey
(No. 124)



Flexi Blue
(No. 44)

You can select the fabric combinations shown for our covers directly from our aids **madita.** (page 8) and **madita-fun.** (page 16).

68 **corzo xcountry go.**

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

I am your foldable rehab buggy.

- Tested according to the requirements of EN ISO 7176/19
- Steel / aluminum frame with footbrake
- Steerable and puncture-resistant front wheels
- Handles (height adjustable only for sizes 1+2)
- Adjustable seat depth and back height
- Height adjustable, foldable and removable footplate
- 5-point harness with padded crotch strap
- Basket
- Soft top
- Cover: black
- Frame: anthracite

HMV-No.: 18.99.01.1032

new **New model:**
now basket & soft top
included!

patron®




Basic modele

corzo xcountry go. size 1	corzo xcountry go. size 2	corzo xcountry go. size 3	corzo xcountry go. size 4
Price / piece	Price / piece	Price / piece	Price / piece
☐ 29 01 100 / 1 2.360 €	☐ 29 02 100 / 1 2.360 €	☐ 29 03 100 / 1 2.490 €	☐ 29 04 100 / 1 2.690 €

Dimensions

		Size 1	Size 2	Size 3	Size 4
A	Backrest height	55 cm	67 cm	65 / 85 cm	69 / 88 cm
B	Seat depth	28 - 34 cm	29 - 35 cm	32 - 39 cm	33 - 41 cm
	Seat width	30 cm	34 cm	38 cm	42 cm
C	Lower leg length	29 - 42 cm	32 - 46 cm	32 - 46 cm	33 - 50 cm
	Backrest angle	90° / 89°	94° / 91°	90° / 88°	94° / 92°
	Seat angle	20°	25°	25°	21°
D	Total length	98 cm	102 cm	102 cm	105 cm
	Total width	52 cm	57 cm	62 cm	65 cm
E	Total height (with adjustable handles over for sizes 1+2)	101 cm	105 cm	117 cm	121 cm
	Weight	12,7 kg	13 kg	14 kg	14,7 kg
	Max. load	45 kg	55 kg	75 kg	75 kg
	Length (folded)	117 cm	122 cm	122 cm	126 cm
	Width (folded)	32 cm	31 cm	32 cm	37 cm
	Height (folded)	34 cm	36 cm	36 cm	36 cm



Accessories

Footplate – incl. calf strap, height adjustable, foldable and removable (standard configuration, if nothing else stated)



Split footplates – incl. calf strap, height adjustable, foldable and removable



Headrest



Backrest extension



Please choose either footplate or split footplates.

Price / piece		Price / piece		Price / piece		Price / piece	
☐ for size 1 (29 01 001)	incl.	☐ for size 1 (29 01 002)	incl.	☐ 29 00 003	185 €	☐ for size 3 (29 03 012)	125 €
☐ for size 2 (29 02 001)	incl.	☐ for size 2-4 (29 02 002)	incl.			☐ for size 4 (29 04 012)	125 €
☐ for size 3 (29 03 001)	incl.						
☐ for size 4 (29 04 001)	incl.						

Accessories

Lateral guide



Chest-shoulder strap



Positioning jacket



Pelvic harness



Price / pair		Price / piece		Price / piece		Price / piece	
☐ for size 1 (29 01 105)	89 €	☐ for size 1 (29 01 006)	148 €	☐ 29 00 007	180 €	☐ Size 1 (29 01 008)	95 €
☐ for size 2-4 (29 02 105)	89 €	☐ for size 2-4 (29 02 006)	148 €			☐ Size 2 (29 02 008)	95 €

Accessories

Arm supports



Soft top – with window for accompanying person



Rain cape – can only be used in conjunction with the soft top



Basket



Price / pair		Price / piece		Price / piece		Price / piece	
☐ 29 00 013	189 €	☐ 29 00 009	incl.	☐ for size 1+2 (29 01 010)	80 €	☐ for size 1+2 (29 01 011)	incl.
				☐ for size 3+4 (29 02 010)	80 €	☐ for size 3+4 (29 02 011)	incl.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The comfortable, foldable rehab buggy.

- Light aluminium frame, rear wheels can be locked with a centrally operated brake
- Detachable wheels, front wheels with directional lock can be steered
- Angle-adjustable backrest (90° - 135°)
- Height-adjustable and detachable foot support
- 5-point strap with padded shoulder and crotch strap
- Incl. canopy and storage basket
- Removable covers
- Frame colour: black

HMV-No.: 18.99.01.1052

patron®



Basic modele

piper.

with black covers



Price / piece

☐ **2001000** / 0 1.889 €

piper.

with blue covers



Price / piece

☐ **2001000** / 1 1.889 €

piper.

with pink covers





Price / piece

☐ **2001000** / 2 1.889 €

Dimensions

		piper.
A	Backrest height	68 cm
B	Seat depth	29 cm
C	Seat width	24 - 35 cm
D	Lower leg length	25 - 37,5 cm
	Backrest angle	90° - 135°
	Seat angle	26°
E	Total length	108 cm
F	Total width	59 cm
G	Total height	100 cm
	Weight	14,5 kg
	Max. load	50 kg
	Length (folded)	100 cm
	Width (folded)	31 cm
	Height (folded)	30 cm



Accessories

Headrest



Price / piece

□ **20 00 001**

164 €

Flexipads set – to be fastened by Velcro



Price / Set

□ **20 00 002**

194 €

Seat padding – incl. 4 side padding (velcro) and one back cushion (seat depth is reduced to 23 cm)



Price / Set

□ Black

(20 00 003/0)

253 €

□ Blue

(20 00 003/1)

253 €

□ Pink

(20 00 003/2)

253 €

Foot positioner



Price / pair

□ **20 00 004**

80 €

Accessories

Footmuff - summer



Price / piece

□ **20 00 005**

175 €

Footmuff - winter



Price / piece

□ **20 00 006**

306 €

Rain cape



Price / piece

□ **20 00 007**

73 €

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____


Commission:
The compact and versatile posterior-walker.

- Lightweight and folding aluminium frame (size 0 cannot be folded)
- Individual height and width adjustable handles
- EVA rear wheels with integrated backwards stop function in the hub, can be switched off
- EVA front wheels with direction locks

HMV-No.: 10.46.02.3029

Dimensions

		Size 0 frame can't be folded	Size 1	Size 2	Size 3	Size 4
A1	Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
A2	Height of forearm supports	-	46 - 64 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
B	Handle spacing	22 - 36 cm	22 - 36 cm	30 - 44 cm	38 - 52 cm	43 - 57 cm
C	Seat height	-	33 cm	45 cm	60 cm	70 cm
D	Total width	52 cm	51 cm	61 cm	68 cm	75 cm
E	Total length	58 cm	58 cm	70 cm	77 cm	87 cm
	Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
	Transport dimensions (L x W x H)	58x52x37 cm	58x51x25 cm	70x61x32 cm	77x68x34 cm	87x75x37 cm
	Max. load	20 kg	30 kg	50 kg	75 kg	100 kg
	Weight	3,9 kg	7,1 kg	8,3 kg	8,5 kg	10,2 kg



Basic model				
malte mini. size 0 Frame colour: Lasurit blue metallic (frame cannot be folded)	malte. size 1	malte. size 2	malte. size 3	malte. size 4
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
☐ 35 00 100 790 €	☐ 75 01 000 825 €	☐ 75 02 000 840 €	☐ 75 03 000 940 €	☐ 75 04 000 1.140 €

Frame colours			
Colour selection for size 1	Colour selection for size 2	Colour selection for size 3+4	Coated with an individual colour (on request) Please provide the requested colour here: RAL _____
			
Price	Price	Price	Price
☐ Telemagenta (Pink) (92 00 012) incl. ☐ Sky blue (92 00 013) incl.	☐ Lasurit blue metallic (92 00 010) incl. ☐ Traffic red (92 00 007) incl.	☐ Anthracit metallic (92 00 014) incl. ☐ Pure white (92 00 100) incl.	☐ 75 00 090 210 €

Accessories					
Forearm armrests with handles – adjustable in depth, height and angle, handles adjustable in depth and angle (for size 0 on request)		Forearm armrests with grip balls – adjustable in depth, height and angle, handles adjustable in depth and angle (for size 0 on request)		Grip bar – pivoting, as a cross-bar between the handles (for sizes 0+4 and in combination with adjacent forearm armrests not possible)	
					
	Length	Width		Length	Width
Size 1	18 cm	6,5 cm	Size 1	18 cm	6,5 cm
Size 2	25 cm	8 cm	Size 2	25 cm	8 cm
Price / set			Price / set		Price / piece
☐ Size 1 (75 01 001) 332 €			☐ Size 1 (75 01 002) 332 €		☐ for size 1 (75 01 003) 272 €
☐ Size 2 (75 02 001) 332 €			☐ Size 2 (75 02 002) 332 €		☐ for size 2 (75 02 003) 272 €
					☐ for size 3 (75 03 003) 272 €

Accessories			
Stops for the hand	arking handbrake (not possible for size 0)		Handbrake – with extra small brake handles at reduced grip width with separate locking function
			
		Assembly in combination with a forearm rest (not possible in combination with the grip ball)	
Price / set		Price / set	Price / piece
☐ 75 00 004 129 €	☐ for size 1 (75 01 008) 235 € ☐ for size 2 (75 02 008) 235 € ☐ for size 3 (75 03 008) 235 € ☐ for size 4 (75 04 008) 235 €		☐ for size 1 (75 01 015) 392 € ☐ for size 2 (75 02 015) 392 € ☐ for size 3 (75 03 015) 392 € ☐ for size 4 (75 04 015) 392 €

Accessories

Combined pelvic and buttock support pad for sizes 1-4 – height, width and depth adjustable, PU padded		Combined pelvic and buttock support pad for size 0 – height, width and depth adjustable, PU padded		4-point pelvic harness for width adjustable pelvic support pad				
								
	Width		Width		Outer part width	Centre section width	Outer part length	Centre section length
Size 1	20 – 28cm			Size 1	7 cm	7 cm	12 cm	12 cm
Size 2	25 – 33cm			Size 2	9 cm	8,5 cm	13 cm	14 cm
Size 3	30 – 36cm			Size 3	11,5cm	11 cm	15 cm	16 cm
Price / set		Price / set		Price / piece				
☐ Size 1 (75 01 005) 385 €		☐ for size 0 (35 01 103) 385 €		☐ Size 1 (75 01 020) 172 €				
☐ Size 2 (75 02 005) 385 €				☐ Size 2 (75 02 020) 172 €				
☐ Size 3 (75 03 005) 385 €				☐ Size 3 (75 03 020) 172 €				

Accessories

Folding seat – padded and covered (not possible for size 0)		Friction brake		Push bar for attendant (not possible for size 0)	
					
Price / piece		Price / set		Price / piece	
☐ for size 1 (75 01 007) 307 € ☐ for size 2 (75 02 007) 307 € ☐ for size 3+4 (75 03 007) 307 €		☐ for size 0+1 (75 00 009) 76 € ☐ for size 2-4 (75 01 009) 76 €		☐ for right (75 00 018) 178 € ☐ for left (75 01 018) 178 €	

Accessories

Anti tip – incl. support roller (for size 0 without support roller)		Deflector set (not possible for size 0)		Bumper wheels – for size 0		Storage basket	
							
Price / set		Price / set		Price / set		Price / piece	
☐ for size 0 (35 00 005) 129 € ☐ for size 1-4 (75 00 019) 129 €		☐ 75 00 012 190 €		☐ 35 00 010 190 €		☐ for size 0 (35 00 011) 65 € ☐ for size 1-4 (75 00 011) 65 €	

malte outdoor.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The compact and versatile posterior-walker for outdoor use.


- Lightweight and folding aluminium frame
- Individual height and width adjustable handles
- Drum brake incl. parking handbrake
- EVA back wheels
- EVA front wheels with direction locks

HMV-No.: 10.46.02.3052



Dimensions

		Size 1	Size 2	Size 3
A1	Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm
A2	Handle spacing	46 - 64 cm	63 - 87 cm	79 - 103 cm
B	Seat height	22 - 36 cm	30 - 44 cm	38 - 52 cm
C	Total width	33 cm	45 cm	60 cm
D	Total length	56 cm	66 cm	73 cm
E	Frame height (folded)	66 cm	75 cm	82 cm
	Transport dimensions (LxWxH)	25 cm	32 cm	34 cm
	Max. load	66x56x25 cm	75x66x32 cm	82x73x34 cm
	Weight	30 kg	50 kg	70 kg
	Weight	7,1kg	8,3kg	9,8kg



Basic model

malte outdoor. size 1 6" tyres (front) and 8" (rear)	malte outdoor. size 2 8" tyres (front) and 10" (rear)	malte outdoor. size 3 8" tyres (front) and 10" (rear)
Price / piece	Price / piece	Price / piece
☐ 75 01 100 965 €	☐ 75 02 100 980 €	☐ 75 03 100 1.080 €

Frame colours

Colour selection for size 1	Colour selection for size 2	Colour selection for size 3	Coated with an individual colour (on request) Please provide the requested colour here: RAL _____
			
Price	Price	Price	Price
☐ Telemagenta (Pink) (92 00 012) incl. ☐ Sky blue (92 00 013) incl.	☐ Lasurit blue metallic (92 00 010) incl. ☐ Traffic red (92 00 007) incl.	☐ Anthracit metallic (92 00 014) incl. ☐ Pure white (92 00 100) incl.	☐ 75 00 090 210 €

Accessories

Combined pelvic and buttock support pad height, width and depth adjustable, PU padded		4-point pelvic harness for width adjustable pelvic support pad				
						
	Width		Outer part width	Centre section width	Outer part length	Centre section length
Size 1	20 – 28 cm	Size 1	7 cm	7 cm	12 cm	12 cm
Size 2	25 – 33 cm	Size 2	9 cm	8,5 cm	13 cm	14 cm
Size 3	30 – 36 cm	Size 3	11,5 cm	11 cm	15 cm	16 cm
		Price / piece				
☐ Size 1 (75 01 005)	385 €	☐ Size 1 (75 01 020)	172 €			
☐ Size 2 (75 02 005)	385 €	☐ Size 2 (75 02 020)	172 €			
☐ Size 3 (75 03 005)	385 €	☐ Size 3 (75 03 020)	172 €			

Accessories

Stops for the hand		Forearm armrests with handles – adjustable in depth, height and angle, handles adjustable in depth and angle		Handbrake – with extra small brake handles at reduced grip width with separate locking function	
					
			Length	Width	
		Size 1	18cm	6,5cm	
		Size 2	25cm	8 cm	
Price / set		Price / set		Upcharge / set	
☐ 75 00 004 129 €		☐ Size 1 (75 01 001) 332 €		☐ for size 1 (75 01 016) 157 €	
		☐ Size 2 (75 02 001) 332 €		☐ for size 2 (75 02 016) 157 €	
				☐ for size 3 (75 03 016) 157 €	

Accessories

Folding seat – padded and covered



Price / piece

☐ for size 1	(75 01 007)	307 €
☐ for size 2	(75 02 007)	307 €
☐ for size 3	(75 03 007)	307 €

Push bar for attendant



Price / piece

☐ for right	(75 00 018)	178 €
☐ for left	(75 01 018)	178 €

Storage basket



Price / piece

☐ 75 00 011	65 €
---	------

Accessories

Anti tip – incl. support roller



Price / set

☐ 75 00 019	129 €
---	-------

Backwards stop – (Retrofitting not possible)



Upcharge / pair

☐ for size 1	(75 01 013)	235 €
☐ for size 2+3	(75 02 013)	235 €

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The anterior-walker for more independence and mobility.


- Lightweight aluminium frame (sizes 2+3 can be folded)
- Individual height and width adjustable handles
- EVA tyres (150 mm for sizes 0+1, 200 mm for sizes 2+3)

HMV-No.: 10.46.02.3049



Dimensions

		Size 0	Size 1	Size 2	Size 3
A1	Handle height	30 – 47 cm	38 – 59 cm	53 – 75 cm	63 – 88 cm
A2	Height of forearm supports	41 – 58 cm	49 – 70 cm	63 – 85 cm	73 – 98 cm
B	Handle spacing	18 – 32 cm	22 – 36 cm	30 – 44 cm	36 – 50 cm
C	Seat height	25 cm	30 cm	45 cm	45 cm
D	Total width	44 cm	48 cm	62 cm	71 cm
E	Total length	36,5 cm	42,5 cm	57,5 cm	61 cm
E	Total length (telescopic wheelbase)	39 – 47 cm	45 – 55 cm	60 – 74 cm	65 – 79 cm
	Max. load	20 kg	30 kg	50 kg	65 kg
	Weight	4,2 kg	5 kg	7,4 kg	9,1 kg



Basic model

marcy. size 0 - 3

	Price / piece
☐ Size 0 (76 00 000)	705 €
☐ Size 1 (76 01 000)	705 €
☐ Size 2 (76 02 000)	705 €
☐ Size 3 (76 03 000)	740 €



marcy. size 0 - 3
(telescopic)

	Price / piece
☐ Size 0 (76 00 100)	915 €
☐ Size 1 (76 01 100)	915 €
☐ Size 2 (76 02 100)	915 €
☐ Size 3 (76 03 100)	990 €

Basic features
with additional,
adjustable
wheelbase



Frame colours			
Colour selection for size 0+1	Colour selection for size 2	Colour selection for size 3	Coated with an individual colour (on request) Please provide the requested colour here: RAL _____
			
Price	Price	Price	Price
☐ Telemagenta (Pink) (92 00 012) incl. ☐ Sky blue (92 00 013) incl.	☐ Lasurit blue metallic (92 00 010) incl. ☐ Traffic red (92 00 007) incl.	☐ Anthracit metallic (92 00 014) incl. ☐ Pure white (92 00 100) incl.	☐ 76 00 090 210 €

Accessories							
Forearm armrests with handles (only possible with telescopic Basic model) adjustable in depth, height and angle - handles adjustable in depth and angle			Forearm armrests with grip balls (only possible with telescopic basic model) adjustable in depth, height and angle - grip balls adjustable in depth and angle			Grip bar – pivoting, as a cross-bar between the handles (in combination with adjacent forearm armrests not possible)	
							
	Length	Width		Length	Width		
Size 1	18 cm	6,5 cm	Size 1	18 cm	6,5 cm		
Size 2	25 cm	8 cm	Size 2	25 cm	8 cm		
Price / set			Price / set			Price / piece	
☐ Size 1	(76 01 001)	344 €	☐ Size 1	(76 01 002)	344 €	☐ for size 0	(76 00 003) 283 €
☐ Size 2	(76 02 001)	344 €	☐ Size 2	(76 02 002)	344 €	☐ for size 1	(76 01 003) 283 €
						☐ for size 2	(76 02 003) 283 €
						☐ for size 3	(76 03 003) 283 €
						☐ for size 0	(76 00 004) 172 €
						☐ for size 1	(76 01 004) 172 €
						☐ for size 2	(76 02 004) 172 €
						☐ for size 3	(76 03 004) 172 €
							
						Seat – fixed, padded and covered incl. back strap	

Accessories			
Parking handbrake 	Handbrake – with extra small brake handles at reduced grip width with separate locking function 	Friction brake 	Backwards stop 
Assembly in combination with the forearm armrest (not possible in combination with grip balls) 			
Price / set	Upcharge / set	Price / set	Upcharge / pair
☐ for size 0 (76 00 005) 244 € ☐ for size 1 (76 01 005) 244 € ☐ for size 2 (76 02 005) 244 € ☐ for size 3 (76 03 005) 244 €	☐ for size 0 (76 00 010) 407 € ☐ for size 1 (76 01 010) 407 € ☐ for size 2 (76 02 010) 407 € ☐ for size 3 (76 03 010) 407 €	☐ for size 0+1 (76 00 006) 85 € ☐ for size 2+3 (76 02 006) 85 €	☐ for size 0+1 (76 00 009) 276 € ☐ for size 2+3 (76 02 009) 276 €

Accessories			
Direction locks 	Storage basket for size 0+1  for size 2+3 	Push bar for attendant 	
Price / pair	Price / piece	Price / piece	
☐ for size 0+1 (76 00 008) 94 € ☐ for size 2+3 (76 02 008) 94 €	☐ for size 0 (76 00 007) 72 € ☐ for size 1 (76 01 007) 72 € ☐ for size 2+3 (76 02 007) 72 €	☐ for right (76 00 012) 185 € ☐ for left (76 01 012) 185 €	

dynamic pacer.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The dynamic gait trainer.

- Modular set with various options to combine different bases and upper frames (except size 4)
- Upper frame to accommodate accessories, continuous height adjustment via gas spring support (for size 1-4)



HMV-No.: 10.46.02.3005
(Size 0)

HMV-No.: 10.46.02.3044
(with standard base)

HMV-No.: 10.46.02.3046
(with utility base)



Standard upper frames



Dynamic Upper frames

Additionally, with integrated and separately adjustable vertical and horizontal dynamic for a physiological natural gait.



	Price / piece		Price / piece
☐ Size 1 (33 01 200) (K622)	1.055 €	☐ Size 1 (33 01 202) (K621)	1.790 €
☐ Size 2 (33 02 200) (K632)	1.450 €	☐ Size 2 (33 02 202) (K631)	2.455 €
☐ Size 3 (33 03 200) (K642)	1.760 €	☐ Size 3 (33 03 202) (K641)	2.765 €
☐ Size 4 (33 04 200) (K652)	2.455 €	☐ Size 4 (33 04 202) (K651)	3.485 €

Standard base

incl. four 5½" (140mm) castors with parking brake, friction brake, backwards stop and direction lock



Utility base

incl. two 8" (200 mm) front wheels with direction lock and two 11 1/2" (290 mm) rear wheels with parking brake, friction brake and backwards stop, for indoor and outdoor use (for outdoor use the function and parking handbrake must be ordered)



Treadmill/stability base

only for size 1 + 3, including four 5½" (140mm) castors with parking brake, friction brake, backwards stop and direction lock (only available in gray)



	Price / piece		Price / piece		Price / piece
☐ Size 1 without odometer (33 01 204) (K624)	680 €	☐ Size 2without odometer (33 02 206) (K636)	1.350 €	☐ standard (33 00 208) (K676)	1.710 €
☐ Size 1 with odometer (33 01 205) (K623)	950 €	☐ Size 2with odometer (33 02 207) (K635)	1.650 €		
☐ Size 2without odometer (33 02 204) (K634)	900 €	☐ Size 3without odometer (33 03 206) (K646)	1.655 €		
☐ Size 2 with odometer (33 02 205) (K633)	1.200 €	☐ Size 3with odometer (33 03 207) (K645)	1.960 €		
☐ Size 3without odometer (33 03 204) (K644)	1.205 €				
☐ Size 3with odometer (33 03 205) (K643)	1.510 €				
☐ Size 4without odometer (33 04 204) (K654)	2.130 €				
☐ Size 4with odometer (33 04 205) (K653)	2.385 €				



pacer mini. (size 0)

one-piece base frame with gripping frame, 4 swivel castors with grinding brake, recoil stop, parking brake and directional fixed spring

		Price / piece
☐ Size 0 without odometer	(33 00 200) (K614)	1.310 €
☐ Size 0 with odometer	(33 00 201) (K613)	1.605 €


Colour selection for size 0 (available from stock)	Colour selection for size 1 (available from stock)	Colour selection for size 2 (available from stock)	Colour selection for size 3 (available from stock)	Colour selection for size 4 (available from stock)
	 	 	 	
Price	Price	Price	Price	Price
☐ Lime (92 00 153) incl.	☐ Lime (92 00 153) incl. ☐ White (92 00 155) incl.	☐ Red (92 00 150) incl. ☐ White (92 00 155) incl.	☐ Blue (92 00 152) incl. ☐ Anthracite (92 00 151) incl.	☐ Anthracite (92 00 151) incl.

These frame colors are also available for all sizes with a delivery time of approx. 10 weeks:



		Price
☐ Anthracite	(92 00 151)	incl.
☐ Red	(92 00 150)	incl.
☐ Blue	(92 00 152)	incl.
☐ Lime	(92 00 153)	incl.
☐ Pink	(92 00 154)	incl.

Dimensions



			Size 0	Size 1	Size 2	Size 3	Size 4
A	Elbow height		39 – 52 cm	47 – 70 cm	61 – 89 cm	81 – 119 cm	86 – 124 cm
B	Total width	Standard base	52 cm	58 cm	66 cm	71 cm	80 cm
		Utility base	-	-	76 cm	81 cm	-
		Schmales Treadmill/stability base	-	89 cm	89 cm	89 cm	-
C	Total length	Standard base	57 cm	66 cm	76 cm	91 cm	101 cm
		Utility base	-	-	91 cm	103 cm	-
		Wide treadmill/stability base	-	104 cm	104 cm	104 cm	-
D	Frame height	Standard-/ Utility base	32 cm	41 – 53 cm	55 – 69 cm	72 – 98 cm	79 – 104 cm
		Wide treadmill/stability base	-	58 – 89 cm	60 – 93 cm	76 – 119 cm	-
E	Base height (without upper frame)	Standard-/ Utility base	-	28 cm	37 cm	41 cm	47 cm
		Wide treadmill/stability base	-	53 cm	53 cm	53 cm	-
	Weight	Standard base	4,6 kg	5,2 kg	6,8 kg	7 kg	10,2 kg
		Utility base	-	-	8,4 kg	8,7 kg	-
		Narrow treadmill/stability ba	-	12,5 kg	12,5 kg	12,5 kg	-
		Standard upper frames	-	2,7 kg	3,4 kg	4,4 kg	4,9 kg
		Dynamic Upper frames	-	3,4 kg	5,2 kg	6,1 kg	6,8 kg
	Dynamic upper frame movement	Vertical	-	3,5 cm	8 cm	8 cm	8 cm
		Horizontal	-	4 cm	5 cm	5 cm	5 cm
	Maximum treadmill width	Narrow treadmill/stability base	-	74 cm	74 cm	74 cm	-
		Wide treadmill/stability base	-	86 cm	86 cm	86 cm	-
	Maximum treadmill height	Wide treadmill/stability base	-	28 cm	28 cm	28 cm	-
Max. load			23 kg	34 kg	68 kg	91 kg	113 kg

Accessories

Multi-position saddle (MPS)

The multi-position saddle (MPS) is the ultimate positioning accessory for gait. It quickly fine tunes pelvis positioning for anterior or posterior use. Adjustable in five directions to accommodate each client's posture, the MPS easily attaches to both the dynamic and standard upper frames.



Seat pad Dimensions	Size 1+2	Size 3
Front width (a)	7,5 cm	11,5 cm
Middle width (b)	4,5 cm	6,5 cm
Rear width (c)	9 cm	13,5 cm
Length (d)	24 cm	27,5 cm

Saddle depth – adjusts to optimally position the pelvis in relation to the upper body to encourage forward motion during ambulation.

Adjustment range: 13 cm



Saddle angle – is adjustable to position the pelvis in anterior or posterior tilt for the most effective gait pattern.

Forward: 7.5° and 15°
Back: 7.5° and 15°



Saddle height – adjusts to finetune the position of the pelvis during gait.

Adjustment range: 22 cm



Hip corral – depth and height adjustments provide cushioning and support for the pelvis at the level of ASIS.

Adjustment range depth: 8 cm
Adjustment range height: 8 cm



Price / piece

☐ Size 1 for size 1	(33 01 012) (K672)	1.523 €
☐ Size 2 for size 2-4	(33 02 012) (K673)	1.576 €
☐ Size 3 for size 2-4	(33 03 012) (K674)	1.711 €

Accessories


MPS Padding cover to reduce frictional resistance

Price / piece

☐ for MPS size 1+2	(33 01 024) (K627)	123 €
☐ for MPS size 3	(33 03 024) (K628)	123 €

Accessories

Pelvic support – incl. bracket and a pair of handholds



	Size 0	Size 1	Size 2
Key user dimension*	15 – 20 cm	15 – 25 cm	20 – 36 cm
a	17 cm	22 cm	30 cm
b	30 cm	41 cm	51 cm
c	18 cm	23 cm	30 cm

* half of thigh circumference

Price / Set

☐ Size 0	(33 00 015) (K556)	688 €
☐ Size 1	(33 01 015) (K557)	750 €
☐ Size 2	(33 02 015) (K558)	814 €

Hip positioner – adjustable, incl. bracket and a pair of handholds



	Size 1	Size 2
Crotch length (a)	6 cm	9 cm
Width (b)	3 cm	4,5 cm
Length (c)	23 cm	28 cm

Price / Set

☐ Size 1	(33 01 003) (K513)	601 €
☐ Size 2	(33 02 003) (K533)	668 €

Pad – for hip positioner



Price / piece

☐ for size 1	(33 01 004) (K541)	217 €
☐ for size 2	(33 02 004) (K542)	217 €

Accessories

Chest prompt – incl. bracket

	Range of circumference	Height
Size 1	36 – 71 cm	11 cm
Size 2	56 – 102 cm	17 cm
Size 3	71 – 127 cm	20 cm
Price / piece		
☐ Size 1 (33 01 002)		800 €
☐ Size 2 (33 02 002)		865 €
☐ Size 3 (33 03 002)		865 €



Accessories

Chest unit – with back strap and holder

	Height	Range of circumference
for size 0	10 cm	36 – 71 cm
for size 1	10 cm	36 – 71 cm
for size 2	12,5 cm	56 – 102 cm
for size 3	15 cm	71 – 127 cm
Price / piece		
☐ for size 0 (33 00 020)		500 €
☐ for size 1 (33 01 020) (K629)		500 €
☐ for size 2 (33 02 020) (K639)		545 €
☐ for size 3 (33 03 020) (K649)		571 €



Handles – with bracket

	Height	Width
Size 1	12,5 cm	9 cm
Size 2	18 cm	9 cm
Price / pair		
☐ Size 1 (33 01 019) (K625)		327 €
☐ Size 2 (33 02 019) (K675)		347 €



Accessories

Arm prompts – incl. handles and bracket

	Padded surface	Adjustment range (elbow to fist)
Size 1	11x18 cm	20 – 32 cm
Size 2	13x22 cm	27 – 41 cm
Price / pair		
☐ Size 1 (33 01 001) (K512)		685 €
☐ Size 2 (33 02 001) (K532)		685 €



Arm platforms – incl. bracket (recommended for sizes 2-4)

	without handles	Padded surface
without handles		34x14 cm
with handles		34x14 cm
Price / pair		
☐ without handles (33 00 010) (K697)		643 €
☐ with handles (33 01 010) (K696)		761 €



Accessories

Functional and parking handbrake

only for use with the utility base in conjunction with arm prompts

(Retrofitting only possible through frame drilling)



Price / pair	
☐ for size 2 (33 02 011) (K687)	503 €
☐ for size 3 (33 03 011) (K688)	552 €

Grip rings – incl. bracket




Price / pair	
☐ 33 00 008 (K579)	316 €

Accessories

Additional handholds for all sizes		Therapy tray – with hollow and transparent cover		Thigh prompts – incl. bracket	
					
		Inner dimensions	22x28cm	Size 1	Length 13 - 20cm
		Max. load	5kg	Size 2	Max. leg circumference 13 - 28cm
Price / pair		Price / piece		Price / pair	
☐ 33 00 002 (K519) 194 €		☐ for size 0 (33 00 001) (K518) 654 € ☐ for size 1-4 (33 00 013) (K698) 654 €		☐ Size 1 (33 01 005) (K515) 597 € ☐ Size 2 (33 02 005) (K535) 710 €	

Accessories

Ankle prompts for size 0 	Ankle prompts for standard base 	Ankle prompts for utility base 	Ankle prompts for treadmill/stability base 
Price / set	Price / set	Price / set	Price / set
☐ for size 0 (33 00 016) (K681) 322 €	☐ for size 1 (33 01 016) (K682) 322 € ☐ for size 2 (33 02 016) (K683) 322 € ☐ for size 3+4 (33 03 016) (K684) 364 €	☐ for size 2 (33 02 017) (K693) 322 € ☐ for size 3 (33 03 017) (K694) 364 €	☐ 33 00 018 (K664) 364 €

Accessories

Attendant guide bar 	Accessories tote 35cm x 40cm x 23cm 														
<table border="1"> <tr> <td></td><td>Price / piece</td></tr> <tr> <td>☐ for size 0 (33 00 009) (K589)</td><td>262 €</td></tr> <tr> <td>☐ for size 1 (33 01 022) (K589)</td><td>262 €</td></tr> <tr> <td>☐ for size 2-3 (33 02 022) (K589)</td><td>262 €</td></tr> <tr> <td>☐ for size 4 (33 04 022)</td><td>262 €</td></tr> </table>		Price / piece	☐ for size 0 (33 00 009) (K589)	262 €	☐ for size 1 (33 01 022) (K589)	262 €	☐ for size 2-3 (33 02 022) (K589)	262 €	☐ for size 4 (33 04 022)	262 €	<table border="1"> <tr> <td></td><td>Price / piece</td></tr> <tr> <td>☐ 33 00 014 (K522)</td><td>149 €</td></tr> </table>		Price / piece	☐ 33 00 014 (K522)	149 €
	Price / piece														
☐ for size 0 (33 00 009) (K589)	262 €														
☐ for size 1 (33 01 022) (K589)	262 €														
☐ for size 2-3 (33 02 022) (K589)	262 €														
☐ for size 4 (33 04 022)	262 €														
	Price / piece														
☐ 33 00 014 (K522)	149 €														

momo tricycle.

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____


Commission:
The tricycle as therapy aid for children, teenagers and adults.

- Aluminium frame with extra-low access point
- Rim brake with separate parking brake
- Angle adjustable handlebar stem
- Steering damper for stabilisation when moving straight on (not applicable for steering for attendant)
- Drum brake in front wheel
- incl. guard plates and basket
- Lighting system according to StVZO with hub dynamo (from 24/20")

new in sizes
16/12" & 24/20"


Dimensions

		12"	16/12"	16"	20"	24/20"	24"	26"	26" (XL)
A*	Pedal to upper edge of saddle	30 - 43 cm	35 - 48 cm	42 - 62 cm	49 - 71 cm	53 - 75 cm	55 - 77 cm	65 - 87 cm	72 - 94 cm
A**	Pedal to upper edge of saddle	33 - 40 cm	38 - 45 cm	45 - 60 cm	52 - 67 cm	56 - 71 cm	58 - 81 cm	68 - 91 cm	75 - 98 cm
B	Saddle support tube to front tube	33 cm	40 cm	39 cm	45 cm	54 cm	52 cm	58 cm	58 cm
C	Wheelbase	68 cm	75 cm	80 cm	91 cm	99 cm	104 cm	112 cm	112 cm
D	Access point height	16 cm	16 cm	18 cm	20 cm	20 cm	22 cm	23 cm	23 cm
E	Total length	102 cm	110 cm	120 cm	140 cm	153 cm	165 cm	178 cm	178 cm
	Total width	58 cm	65 cm	65 cm	69 cm	75 cm	75 cm	75 cm	75 cm
	Max. load	40 kg	50 kg	60 kg	80 kg	90 kg	100 kg	120 kg	120 kg
	Weight	12 kg	13,5 kg	14,5 kg	17 kg	18,5 kg	19,5 kg	21,5 kg	21,7 kg



* Dimensions with standard saddle post / ** Dimensions with T-saddle post

Basic model

momo dreirad. 12" (HMV-No.: 22.51.02.0053)	momo dreirad. 16/12" (HMV-No.: 22.51.02.0054)	momo dreirad. 16" (HMV-No.: 22.51.02.0054)	momo dreirad. 20" (HMV-No.: 22.51.02.0055)
Basic features, without drum brake and rear guard plate, additionally with pennant	Basic features, additionally with pennant	Basic features, additionally with pennant	Basic features, additionally with pennant
Price / piece	Price / piece	Price / piece	Price / piece
☐ 37 01 000 870 €	☐ 37 02 100 1.300 €	☐ 37 02 000 1.300 €	☐ 37 03 000 1.420 €

Basic model

momo dreirad. 24/20" (HMV-No.: 22.51.02.0051)	momo dreirad. 24" (HMV-No.: 22.51.02.0051)	momo dreirad. 26" (HMV-No.: 22.51.02.0052)	momo dreirad. 26" (XL) (HMV-No.: 22.51.02.0052)
Basic features	Basic features	Basic features	Basic features
Price / piece	Price / piece	Price / piece	Price / piece
☐ 37 04 100 1.650 €	☐ 37 04 000 1.650 €	☐ 37 05 000 1.775 €	☐ 37 06 000 1.940 €

Frame colours

Colour selection for 12"	Colour selection for 16" + 16/12" + 20"	Colour selection for 24" + 24/20" + 26" + 26" XL
  	  	 
Price	Price	Price
☐ Telemagenta (Pink) (92 00 012) incl.	☐ Melon yellow (92 00 006) incl.	☐ Anthracit metallic (92 00 014) incl.
☐ Sky blue (92 00 013) incl.	☐ Marine blue (92 00 008) incl.	☐ Pure white (92 00 100) incl.
☐ Pure white (92 00 100) incl.	☐ Pure white (92 00 100) incl.	

Configuration possibilities - drive

	Art.-No.	Drive	Upcharge / piece
	☐ 37 00 001	Rigid sprocket (standard configuration, if nothing else stated)	incl.
	☐ 37 00 004	Freewheel brake hub (with back pedal brake)	171 €
	☐ 37 00 403	Change from rigid sprocket to neutral hub	390 €
	☐ 37 00 002	Extra light transmission (only for above drive types)	incl.
	☐ 37 00 303	Switchover from rigid sprocket, neutral hub and 3-gears-freewheel brake hub with reverse option	610 €
	☐ 37 00 005	3-gear freewheel brake hub (with back pedal brake)	226 €
	☐ 37 00 007	3-gear freewheel brake hub with reverse option	379 €
	☐ 37 00 008	7-gear freewheel brake hub (with back pedal brake)	429 €
	☐ 37 00 202	7-gear-freewheel hub with roll brake (brake lever mounted on the handlebar, without back pedal brake)	540 €

Saddle

Standard saddle (Standard configuration, if nothing else stated)			Gel saddle				Saddle seat with leg hollows				
											
	Width	Length		Rear width	Front width	Length		Rear width	Front width	Rear length	Total length
Size 1	15cm	21cm	Size 1	19 cm	4 cm	24 cm	Size 1	20,5cm	4 cm	9,5cm	14,5cm
Size 2	15cm	24 cm	Size 2	24 cm	7 cm	27 cm	Size 2	24 cm	4 cm	9,5cm	14,5cm
Size 3	18cm	26cm									
Price / piece			Upcharge / piece				Upcharge / piece				
☐ Size 1 (37 01 001)		incl.	☐ Size 1 (37 01 024)			57 €	☐ Size 1 (37 01 003)				216 €
☐ Size 2 (37 02 001)		incl.	☐ Size 2 (37 02 024)			57 €	☐ Size 2 (37 02 003)				216 €
☐ Size 3 (37 03 001)		incl.									


Saddle

Saddle supports

Unicycle saddle – banana shape (caution: min. inside leg length increases by 2 cm)				Moped saddle – (caution: min. inside leg length increases by 2 cm)				Standard saddle support (standard configuration, if nothing else stated)		T-saddle support – for horizontal saddle adjustment	
											
Rear width	Middle width	Front width	Length	Rear width	Front width	Rear length	Total length	for		for	
11 cm	6 cm	8 cm	25 cm	26 cm	9,5 cm	12 cm	25 cm	Size 1	12" + 16"/12"	Size 1	12" + 16"/12"
								Size 2	16"	Size 2	16" + 20"
								Size 3	20" – 26"	Size 3	24" + 24"/20" + 26"
Upcharge / piece				Upcharge / piece				Price / piece		Upcharge / piece	
☐ Size 1 (37 00 010) 106 €				☐ Size 1 (37 00 023) 158 €				☐ Size 1 (37 01 004)	incl.	☐ Size 1 (37 01 005)	83 €
								☐ Size 2 (37 02 004)	incl.	☐ Size 2 (37 02 005)	83 €
								☐ Size 3 (37 03 004)	incl.	☐ Size 3 (37 03 005)	83 €

Holding brackets

Headrest

Holding brackets incl. holder – only for use together with the T-saddle support				Headrest – only for use together with a holding bracket			
		Max. Height*	Depth				
Size 1		20 cm	12 cm				
Size 2		30 cm	12 cm				
Size 3		37 cm	12 cm				
Size 4		53 cm	12 cm				
Tailored		____ cm (max. 85 cm)	____ cm (max. 30 cm)				
						Width	Height
						Size 1	20 cm
						Size 2	23 cm
							
						Price / set	
☐ Size 1 (37 01 007)		129 €	☐ Size 3 (37 03 007)		129 €	☐ Size 1 (37 01 029)	343 €
☐ Size 2 (37 02 007)		129 €	☐ Size 4 (37 04 007)		129 €	☐ Size 2 (37 02 029)	343 €
			☐ Tailored (37 09 007)		272 €		

* Max. height: measured from the standard saddle to the top edge of the respective support pad

Torso supports

Pelvic supports

Dynamic back support pad padded with strap and bracket (only for use together with a holding bracket)		Width adjustable back support pad padded, incl. bracket (only for use together with a holding bracket)		Dynamic pelvic guide support pad padded with strap and bracket (only for use together with a holding bracket)		Width adjustable pelvic guide support pad – padded, incl. bracket (only for use together with a holding bracket)	
							
Width		Width		Width		Width	
Size 1	20 – 30 cm	Size 1	20 – 28 cm	Size 1	20 – 30 cm	Size 1	20 – 28 cm
Size 2	25 – 35 cm	Size 2	25 – 33 cm	Size 2	25 – 35 cm	Size 2	25 – 33 cm
Size 3	30 – 40 cm	Size 3	30 – 36 cm	Size 3	30 – 40 cm	Size 3	30 – 36 cm
Price / piece		Price / piece		Price / piece		Price / piece	
☐ Size 1 (37 01 008)	243 €	☐ Size 1 (37 01 025)	324 €	☐ Size 1 (37 01 009)	243 €	☐ Size 1 (37 01 026)	324 €
☐ Size 2 (37 02 008)	243 €	☐ Size 2 (37 02 025)	324 €	☐ Size 2 (37 02 009)	243 €	☐ Size 2 (37 02 026)	324 €
☐ Size 3 (37 03 008)	243 €	☐ Size 3 (37 03 025)	324 €	☐ Size 3 (37 03 009)	243 €	☐ Size 3 (37 03 026)	324 €

Harness system

Chest strap – for width adjustable back support pad			Positioning jacket – for width adjustable back support pad, incl. universal bracket			
						
	Width	Length	Inside width	Total width	Side height	Total height
Size 1	5 cm	18,5 cm	Size 1	8 cm	25 cm	10 cm
Size 2	6,5 cm	25,5 cm	Size 2	9 cm	30 cm	12 cm
Size 3	7 cm	30,5 cm	Size 3	12 cm	34 cm	16 cm
Price / piece			Price / piece			
☐ Size 1 (37 01 050)		98 €	☐ Size 1 (37 01 051)		243 €	
☐ Size 2 (37 02 050)		98 €	☐ Size 2 (37 02 051)		243 €	
☐ Size 3 (37 03 050)		98 €	☐ Size 3 (37 03 051)		243 €	

Harness system

T-shape groin harness – for width adjustable pelvic support pad						4-point pelvic harness – for width adjustable pelvic support pad					
	Seat width	Front width	Seat length	Total length			Outer part width	Centre section width	Outer part length	Centre section length	
Size 1	24 cm	23 cm	15 cm	32 cm		Size 1	7 cm	7 cm	12 cm	12 cm	
Size 2	26 cm	26 cm	17 cm	34 cm		Size 2	9 cm	8,5 cm	13 cm	14 cm	
Size 3	28 cm	28 cm	19 cm	38 cm		Size 3	11,5 cm	11 cm	15 cm	16 cm	
Price / piece					Price / piece						
☐ Size 1 (37 01 052)					166 €	☐ Size 1 (37 01 053)					170 €
☐ Size 2 (37 02 052)					166 €	☐ Size 2 (37 02 053)					170 €
☐ Size 3 (37 03 052)					166 €	☐ Size 3 (37 03 053)					170 €

Accessories

Universal bracket for harness used without support pads (width = 27.5 cm)		Hand positioner – for all handlebars (if no quantity was specified: standard = 1 piece)		Steering stop – adjustable up to direction lock (incl. for steering for attendant)		Handlebar extension	
			Front width				
		Size 0	8,5 cm				
		Size 1	10 cm				
		Size 2	12 cm				
		Size 3	14 cm				
Price / piece		Price / piece		Menge		Price / piece	
☐ 37 02 055		110 €					
		☐ Size 0 (37 00 054)	70 €			☐ Size 1 (37 01 006)	129 €
		☐ Size 1 (37 01 054)	70 €			☐ Size 2 (37 03 006)	129 €
		☐ Size 2 (37 02 054)	70 €				
		☐ Size 3 (37 03 054)	70 €				
						Price / piece	
						☐ Size 1 (37 01 022)	29 €
						☐ Size 2 (37 02 022)	29 €

Handlebar

Universal bracket – black (standard configuration, if nothing else stated)			Touring handlebar			Round handlebar			Multifunctional handlebar suitable for 20" – 26"		
		Standard for									
		Width									
	Size 1	12"					Size 1	40 cm			
	Size 2	16/12" + 16"					Size 2	43 cm			
	Size 3	20"					Size 3	48 cm			
	Size 4	24/20" – 26"		Size 1	55 cm					Size 1	61,5 cm
				Size 2	58 cm						
Price / piece			Upcharge / piece			Upcharge / piece			Upcharge / piece		
☐ Size 1 (37 01 011)			incl.			☐ Size 1 (37 01 013)			111 €		
☐ Size 2 (37 02 011)			incl.			☐ Size 2 (37 02 013)			111 €		
☐ Size 3 (37 03 011)			incl.			☐ Size 3 (37 03 013)			111 €		
☐ Size 4 (37 04 011)			incl.						☐ 37 02 014		
									162 €		

Handlebar brake

Round handlebar with brake lever ring				Round handlebar with brake function – the brake function is simply activated by pressing down the round handlebars			
	Width	Depth			Width	Depth	
Size 1	40 cm	17 cm		Size 1	40 cm	17 cm	
Size 2	43 cm	21 cm		Size 2	43 cm	21 cm	
Size 3	48 cm	25 cm		Size 3	48 cm	25 cm	
Upcharge / piece			Upcharge / piece				
☐ Size 1 (3701028)			329 €	☐ Size 1 (3701031)			329 €
☐ Size 2 (3702028)			329 €	☐ Size 2 (3702031)			329 €
☐ Size 3 (3703028)			329 €	☐ Size 3 (3703031)			329 €

Accessories

Push bar height adjustable and removable with integrated twist lock 	Parking brake for attendant (can only be used in connection with a push bar) 	Steering for attendant – via deflection bar, incl. steering stop, in this case the steering damper is omitted (can only be used in connection with a push bar)  <table border="1" data-bbox="1058 274 1404 491"> <thead> <tr> <th colspan="3">Price / piece</th> </tr> </thead> <tbody> <tr> <td>☐ for 12"</td> <td>(37 01 010)</td> <td>432 €</td> </tr> <tr> <td>☐ for 16/12"</td> <td>(37 02 110)</td> <td>432 €</td> </tr> <tr> <td>☐ for 16"</td> <td>(37 02 010)</td> <td>432 €</td> </tr> <tr> <td>☐ for 20"</td> <td>(37 03 010)</td> <td>432 €</td> </tr> <tr> <td>☐ for 24/20"</td> <td>(37 04 110)</td> <td>432 €</td> </tr> <tr> <td>☐ for 24"</td> <td>(37 04 010)</td> <td>432 €</td> </tr> <tr> <td>☐ for 26"</td> <td>(37 05 010)</td> <td>432 €</td> </tr> </tbody> </table>	Price / piece			☐ for 12"	(37 01 010)	432 €	☐ for 16/12"	(37 02 110)	432 €	☐ for 16"	(37 02 010)	432 €	☐ for 20"	(37 03 010)	432 €	☐ for 24/20"	(37 04 110)	432 €	☐ for 24"	(37 04 010)	432 €	☐ for 26"	(37 05 010)	432 €
Price / piece																										
☐ for 12"	(37 01 010)	432 €																								
☐ for 16/12"	(37 02 110)	432 €																								
☐ for 16"	(37 02 010)	432 €																								
☐ for 20"	(37 03 010)	432 €																								
☐ for 24/20"	(37 04 110)	432 €																								
☐ for 24"	(37 04 010)	432 €																								
☐ for 26"	(37 05 010)	432 €																								
Price / piece ☐ 37 00 011 179 €	Price / piece ☐ 37 00 012 119 €																									

Crank

Progressively adjustable crank reducer for 12"–20" there is no widening of the pedal spacing and for 24"+26" there is a widening of the pedal spacing of approx. 18 mm on each side 	Crank reducer – minimal pedal distance widening – reduction in size by 2.5 or 5 cm respectively 	Special crank for knee contracture for 20" – 26" 
Price / pair ☐ for 12"–20"+16/12"+24/20" (37 00 016) 176 € ☐ for 24"+26" (37 04 016) 176 €	Price / pair ☐ for 16"+20"+24/20" (37 02 015) 91 € ☐ for 24"+26" (37 05 015) 91 €	Price / piece ☐ mounted on right (37 03 015) 377 € ☐ mounted on left (37 04 015) 377 €

Pedals

Foot fixing pedals – with racing hooks and compensation weights		Home exerciser pedals		Foot cups – width adjustable, with fixing straps & with magnetic lock	
					
Depth		Width			
13 cm		12 cm			
Upcharge / pair		Upcharge / pair		Price / pair	
☐ 37 00 013 91 €		☐ 37 00 014 62 €		☐ Size 1 (37 01 037) 213 €	
				☐ Size 2 (37 02 037) 213 €	
				☐ Size 3 (37 03 037) 213 €	

Foot cups with leg guidance


Foot cups with leg guidance – width adjustable and with fixing straps & with magnetic lock							
	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	21 cm	26 cm
Size 4	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	28 cm	36 cm
							Price / pair
☐ Size 1 (37 01 038)							361 €
☐ Size 2 (37 02 038)							361 €
☐ Size 3 (37 03 038)							361 €
☐ Size 4 (37 04 038)							361 €



Foot cups with dynamic leg guidance

Foot cups with dynamic leg guidance

width adjustable and with fixing straps
& with magnetic lock



	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide Height	Max. leg guide Height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	12 cm	15 cm	8,8 cm	11,8 cm	24 cm	20,5 cm	25,5 cm

Price / pair

- ☐ Size 1 (**37 01039**)
☐ Size 2 (**37 02039**)
☐ Size 3 (**37 03039**)

472 €
472 €
472 €

Configuration possibilities - lighting system

	Art.-No.	Lighting system	for	Upcharge / set
	☐ 37 03021	Upgrade with a lighting system with side runners and reflectors according to StVZO (mounted on the rear wheel)	12"	129 €
	☐ 37 00021	Upgrading of a battery light system attachable incl. reflectors according to StVZO	12" - 26"	169 €
	☐ 37 01021	Upgrading of a battery light system attachable incl. reflectors according to StVZO	16/12"+ 16"+20"	169 €

momo motion. - The tricycle with E-drive

A motor support when using a tricycle is useful for many users and necessary due to the indication. Thus, for some users, the starting situation is often the only obstacle to independent use of the bike. With clinical pictures such as dystrophies, permanent support is required to compensate the relative loss of muscle strength.

In order to make the medical benefits of riding a tricycle accessible to these people, we have developed our **momo motion.**

The corresponding order form for our **momo motion.** as well as further information on the possibilities of motor support can be found on page 94 or on our homepage at **www.schuchmann.de/en/momo-motion**.



92 **follow-me momo.**

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The clever parent-child tandem coupling for momo tricycles.

- To suit all Schuchmann tricycles with free wheeling and / or switchover to free wheeling
- For all 26" and 28" adult's bicycles with and without rear suspension with quick release wheel fixation
- Wheels with full axle or wheel hub gears can be prepared using optionally available adaptors
- The child's seat, luggage carrier and bags can still be fully used on the leading bicycle
- It is possible to attach and remove the tricycle within seconds
- Incl. special axle nuts



Dimensions

	Size 1	Size 2	Size 3
Required size of the leading bicycle	26" - 28"	26" - 28"	26" - 28"
Max. tricycle size	12" (to be used also for 16" + 20")	16/12" + 16" + 20"	24/20" + 24" + 26"
Max. traction weight	40 kg	60 kg	80 kg
Weight	5,1kg	5,1kg	8,7 kg

Basic model		
follow-me momo. size 1 for 12"	follow-me momo. size 2 for 16/12" + 16" + 20"	follow-me momo. size 3 for 24/20" + 24" + 26"
Price / piece	Price / piece	Price / piece
☐ 39 00 100 489 €	☐ 39 01 100 489 €	☐ 39 02 100 720 €

Adapter for rear axle of the leading bicycle				
Hub or chain gears and quick release with a max. frame width of 150 mm	Hub or chain gears and quick release with a max. frame width of 150 mm to 175 mm	Hub or chain gears with full axle and nuts	Hub gears and clickbox	Hub or chain gears and thru axle
				
i One of these four adapters is included. Each additional adapter will be charged at 88 € / set.				
Price / set	Price / set	Price / set	Price / set	Price / set
☐ 39 01 101 incl.	☐ 39 02 101 incl.	☐ 39 01 102 incl.	☐ 39 01 103 incl.	☐ 39 01 104 auf Anfrage

Accessories	
Suspension bracket – for luggage carrier with „Top-Light“, incl. carabiner hooks and assembly instructions (incl. for size 3)	Preparation of the follow-me momo. – adjustment of the follow-me momo. to wheel size of the tricycle (at Schuchmann)
	
Price / piece	Net price (not discountable)
☐ 39 01 109 26 €	☐ 39 01 110 38 €

Note	
Adult bike with rear engine If the adult bike on which the follow-me momo. is to be mounted is equipped with a rear engine, then an adapter other than described above is needed. Please contact us in such case and we will work with you to develop an adapter solution.	momo motion. or momo tricycle. with electric drive in combination with the follow-me momo. Generally we do not recommend to pair a momo motion. or momo tricycle. , which is equipped with an electric drive, via follow-me momo. with an adult bike. This can lead to dangerous situations due to the support that comes from the electric drive. Please contact us if you still require this combination. We will gladly try to find a solution together with you.
i If the cable of a rear motor is routed through the axle, compatibility with the follow-me momo. is not given.	

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The tried-and-tested tricycle with E-drive.

- With Heinzmann Direct Power Motor (36 V, 250 W, impulse torque 60 Nm)
- Aluminium frame with extra-low access point
- Rim brake with separate parking brake
- Angle adjustable handlebar stem
- Handlebar damper for stabilisation when moving straight on
- Lighting system according to StVZO, which is powered by the Li-ion high-power battery (14.25 Ah)
- Incl. mudguards, basket and pennant

new


in size
24/20"

HMV-No.: 22.51.04.0005



Dimensions

		16"	20"	24/20"	24"	26"	26" (XL)
A*	Pedal to upper edge of saddle	42 – 62 cm	49 – 71 cm	54 – 68,5 cm	55 – 77 cm	65 – 87 cm	72 – 94 cm
A**	Pedal to upper edge of saddle	45 – 60 cm	52 – 67 cm	56 – 72 cm	58 – 81 cm	68 – 91 cm	75 – 98 cm
B	Saddle support tube to front tube	39 cm	45 cm	54 cm	52 cm	58 cm	58 cm
C	Wheelbase	80 cm	91 cm	99 cm	104 cm	112 cm	112 cm
D	Access point height	18 cm	20 cm	20 cm	22 cm	23 cm	23 cm
E	Total length	120 cm	140 cm	153 cm	165 cm	178 cm	178 cm
	Total width	65 cm	69 cm	75 cm	75 cm	75 cm	75 cm
	Max. load	60 kg	80 kg	90 kg	100 kg	120 kg	120 kg
	Weight	24,5 kg	27 kg	28,5 kg	29,5 kg	31,5 kg	31,7 kg



* Dimensions with standard saddle post / ** Dimensions with T-saddle post

Basic model

momo motion. 16"	momo motion. 20"	momo motion. 24/20"	momo motion. 24"	momo motion. 26"	momo motion. 26" (XL)
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
☐ 36 02 000 4.350 €	☐ 36 03 000 4.470 €	☐ 36 04 100 4.565 €	☐ 36 04 000 4.565 €	☐ 36 05 000 4.695 €	☐ 36 06 000 4.865 €

Frame colours

Colour selection for 16" + 20"	Colour selection for 24/20" + 24" + 26" + 26" XL
☐ Melon yellow (92 00 006) ☐ Marine blue (92 00 008) ☐ Pure white (92 00 100)	☐ Anthracit metallic (92 00 014) ☐ Pure white (92 00 100)
Price	Price
incl.	incl.
incl.	incl.
incl.	incl.

E - drive

DirectPower E-drive from Heinzmann

The **momo motion.** is equipped with the DirectPower E-drive from Heinzmann. This allows us to offer you a wear-free and silent direct drive. The motor is designed to be brushless and without gears. This minimises maintenance and wear and eliminates motor noise. The recuperative system automatically recharges when driving downhill or braking. This means that the range can be increased by up to 15 %. The existing force of the user is measured by a force sensor integrated in the bottom bracket, and the assistance from the electric drive is adjusted.

Choose the suitable E-drive for the **momo motion.** from the following configurations. You can choose whether the drive should be equipped **with or without start-up assistance**. In addition, you can choose between a **beginner version** and an **advanced version**. The difference here is the **start-up** and the **top speed** of the various modes **Eco**, **Standard** and **Power**, which can be selected via the display. An analysis cable can be used to change the selected version at a later date.

E - drive configuration

DirectPower E-drive from Heinzmann – without start-up assistance

for size	Version	Start-up handling / top speed			Price / piece
		Eco	Standard	Power	
16"	☐ Beginner (36 02 070)	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	incl.
	☐ Advanced (36 02 071)	medium / 6 km/h	strong / 8 km/h	strong / 10 km/h	incl.
20"	☐ Beginner (36 03 070)	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	incl.
	☐ Advanced (36 03 071)	medium / 6 km/h	strong / 8 km/h	strong / 12 km/h	incl.
24" & 24/20"	☐ Beginner (36 04 070)	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	incl.
	☐ Advanced (36 04 071)	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	incl.
26"	☐ Beginner (36 05 070)	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	incl.
	☐ Advanced (36 05 071)	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	incl.

Configuration E - drive

DirectPower E-drive from Heinzmann – with start-up assistance

Select whether the start-up assistance should be activated by pressing the button **(Fig. 1)** or by twist handle **(Fig. 2)**



for size	Version	Start-up assistance	Start-up handling / top speed			Upcharge / piece
			Eco	Standard	Power	
16"	☐ Beginner with twist handle (36 02 072)	4 km/h	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	117 €
	☐ Beginner with button (36 02 073)					117 €
	☐ Advanced with twist handle (36 02 074)	6 km/h	medium / 6 km/h	strong / 8 km/h	strong / 10 km/h	117 €
	☐ Advanced with button (36 02 075)					117 €
20"	☐ Beginner with twist handle (36 03 072)	4 km/h	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	117 €
	☐ Beginner with button (36 03 073)					117 €
	☐ Advanced with twist handle (36 03 074)	6 km/h	medium / 6 km/h	strong / 8 km/h	strong / 12 km/h	117 €
	☐ Advanced with button (36 03 075)					117 €
24" & 24/20"	☐ Beginner with twist handle (36 04 072)	6 km/h	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	117 €
	☐ Beginner with button (36 04 073)					117 €
	☐ Advanced with twist handle (36 04 074)	6 km/h	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	117 €
	☐ Advanced with button (36 04 075)					117 €
26"	☐ Beginner with twist handle (36 05 072)	6 km/h	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	117 €
	☐ Beginner with button (36 05 073)					117 €
	☐ Advanced with twist handle (36 05 074)	6 km/h	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	117 €
	☐ Advanced with button (36 05 075)					117 €

Software

Subsequent conversion/adjustment of the software on the delivered product by our sales force

Net price (not discountable) / Stück

☐ 36 00 090


152 €

Configuration possibilities - drive

 	Art.-No.	Drive	Upcharge / piece
	☐ 3600001	Freewheel brake hub (with back pedal brake)	incl.
	☐ 3601001	3-gear freewheel brake hub (with back pedal brake)	226 €
	☐ 3602001	3-gear freewheel brake hub with reverse option	379 €
	☐ 3603001	7-gear freewheel brake hub (with back pedal brake)	429 €
	☐ 3604001	7-gear-freewheel hub with roll brake (brake lever mounted on the handlebar, without back pedal brake)	540 €

Saddle

Standard saddle (Standard configuration, if nothing else stated)			Gel saddle				Saddle seat with leg hollows				
											
	Width	Length		Rear width	Front width	Length		Rear width	Front width	Rear length	Total length
Size 1	15cm	21cm					Size 1	20,5cm	4 cm	9,5 cm	14,5 cm
Size 2	15cm	24 cm	Size 1	19 cm	4 cm	24 cm	Size 2	24 cm	4 cm	9,5 cm	14,5 cm
Size 3	18cm	26 cm	Size 2	24 cm	7 cm	27 cm					
Price / piece			Upcharge / piece				Upcharge / piece				
☐ Size 1 (3701001)			☐ Size 1 (3701024)				☐ Size 1 (3701003)				
☐ Size 2 (3702001)			☐ Size 2 (3702024)				☐ Size 2 (3702003)				
☐ Size 3 (3703001)											
incl.			57 €				216 €				
incl.			57 €				216 €				
incl.											

Holding brackets				Headrest			
Holding brackets incl. holder – only for use together with the T-saddle support				Headrest – only for use together with a holding bracket			
	max. Height*	Depth					
Size 1	20 cm	12 cm					
Size 2	30 cm	12 cm					
Size 3	37 cm	12 cm					
Size 4	53 cm	12 cm					
Tailored	____ cm (max. 85 cm)	____ cm (max. 30 cm)					
* Max. height: measured from the standard saddle to the top edge of the respective support pad							
					Width	Height	
				Size 1	20 cm	15 cm	
				Size 2	23 cm	18 cm	
				Price / set			
☐ Size 1	(37 01 007)			129 €			
☐ Size 2	(37 02 007)			129 €			
☐ Size 3	(37 03 007)			129 €			
☐ Size 4	(37 04 007)			129 €			
☐ Tailored	(37 09 007)			272 €			
				Price / piece			
				☐ Size 1	(37 01 029)		343 €
				☐ Size 2	(37 02 029)		343 €

Torso supports				Pelvic supports			
Dynamic back support pad padded with strap and bracket (only for use together with a holding bracket)		Width adjustable back support pad padded, incl. bracket (only for use together with a holding bracket)		Dynamic pelvic guide support pad padded with strap and bracket (only for use together with a holding bracket)		Width adjustable pelvic guide support pad – padded, incl. bracket (only for use together with a holding bracket)	
							
	Width		Width		Width		Width
Size 1	20 – 30 cm	Size 1	20 – 28 cm	Size 1	20 – 30 cm	Size 1	20 – 28 cm
Size 2	25 – 35 cm	Size 2	25 – 33 cm	Size 2	25 – 35 cm	Size 2	25 – 33 cm
Size 3	30 – 40 cm	Size 3	30 – 36 cm	Size 3	30 – 40 cm	Size 3	30 – 36 cm
Price / piece		Price / piece		Price / piece		Price / piece	
☐ Size 1	(37 01 008)	☐ Size 1	(37 01 025)	☐ Size 1	(37 01 009)	☐ Size 1	(37 01 026)
	243 €		324 €		243 €		324 €
☐ Size 2	(37 02 008)	☐ Size 2	(37 02 025)	☐ Size 2	(37 02 009)	☐ Size 2	(37 02 026)
	243 €		324 €		243 €		324 €
☐ Size 3	(37 03 008)	☐ Size 3	(37 03 025)	☐ Size 3	(37 03 009)	☐ Size 3	(37 03 026)
	243 €		324 €		243 €		324 €

Harness system							
Chest strap – for width adjustable back support pad				Positioning jacket – for width adjustable back support pad, incl. universal bracket			
							
	Width	Length			Inside width	Total width	Side height
Size 1	5 cm	18,5 cm		Size 1	8 cm	25 cm	10 cm
Size 2	6,5 cm	25,5 cm		Size 2	9 cm	30 cm	12 cm
Size 3	7 cm	30,5 cm		Size 3	12 cm	34 cm	16 cm
Price / piece				Price / piece			
☐ Size 1	(37 01 050)		98 €	☐ Size 1	(37 01 051)		243 €
☐ Size 2	(37 02 050)		98 €	☐ Size 2	(37 02 051)		243 €
☐ Size 3	(37 03 050)		98 €	☐ Size 3	(37 03 051)		243 €

Harness system

T-shape groin harness – for width adjustable pelvic support pad

	Seat width	Front width	Seat length	Total length
Size 1	24 cm	23 cm	15 cm	32 cm
Size 2	26 cm	26 cm	17 cm	34 cm
Size 3	28 cm	28 cm	19 cm	38 cm



Price / piece

- ☐ Size 1 **(37 01 052)**
☐ Size 2 **(37 02 052)**
☐ Size 3 **(37 03 052)**

166 €
166 €
166 €

4-point pelvic harness – for width adjustable pelvic support pad

	Outer part width	Centre section width	Outer part length	Centre section length
Size 1	7 cm	7 cm	12 cm	12 cm
Size 2	9 cm	8,5 cm	13 cm	14 cm
Size 3	11,5 cm	11 cm	15 cm	16 cm



Price / piece

- ☐ Size 1 **(37 01 053)**
☐ Size 2 **(37 02 053)**
☐ Size 3 **(37 03 053)**

170 €
170 €
170 €

Accessories

Universal bracket

for harness used without support pads (width = 27,5 cm)



Steering stop – adjustable up to direction lock (incl. for steering for attendant)



Handlebar extension



for
16"–26"

for
16"–26" Length
10 cm

Price / piece

Price / piece

Price / piece

☐ **37 02 055**

110 €

☐ **37 03 006**

129 €

☐ **37 02 022**

29 €

Handlebar

Universal bracket – black (standard configuration, if nothing else stated)



	Standard for	Width
Size 2	16"	49 cm
Size 3	20"	52 cm
Size 4	24/20" – 26"	58 cm

Price / piece

- ☐ Size 2 **(37 02 011)**
☐ Size 3 **(37 03 011)**
☐ Size 4 **(37 04 011)**

incl.
incl.
incl.

Touring handlebar



	Width	Depth
Size 1	55 cm	17 cm
Size 2	58 cm	17 cm

Upcharge / piece

- ☐ Size 1 **(37 01 012)** 111 €
☐ Size 2 **(37 02 012)** 111 €

Round handlebar



	Width	Depth
Size 2	43 cm	21 cm
Size 3	48 cm	25 cm

Upcharge / piece

- ☐ Size 2 **(37 02 013)** 111 €
☐ Size 3 **(37 03 013)** 111 €

Multifunctional handlebar suitable for 20" – 26"



	Width	Depth
Size 1	61,5 cm	16,5 cm

Upcharge / piece

- ☐ **37 02 014** 162 €

Accessories

Push bar for attendant

height adjustable and removable with integrated twist lock



Price / piece

☐ for alle Größen **(36 00 012)**

179 €

Brake for attendant – incl. motor switch-off and fitting of a disc brake on the front wheel (can only be used in connection with push bar)

Price / piece

☐ **36 07 012**

473 €

Crank

Progressively adjustable crank reducer for 12"-20" there is no widening of the pedal spacing and for 24"+26" there is a widening of the pedal spacing of approx. 18mm on each side 		Crank reducer minimal pedal distance widening, reduction in size by 2.5 or 5 cm respectively 		Special crank for knee contracture for 20" - 26" 	
Price / pair		Price / pair		Price / piece	
☐ for 16"+20"+24/20" (37 00 016)	176 €	☐ for 16"+20"+24/20" (37 02 015)	91 €	☐ mounted on right (37 03 015)	377 €
☐ for 24"+26" (37 04 016)	176 €	☐ for 24"+26" (37 05 015)	91 €	☐ mounted on left (37 04 015)	377 €

Pedals

Foot cups

Foot fixing pedals – with racing hooks and compensation weights 		Home exerciser pedals 		Foot cups – width adjustable, with fixing straps & with magnetic lock 	
Depth		Width			
13 cm		12 cm			
Upcharge / pair		Upcharge / pair		Price / pair	
☐ 37 00 013 91 €		☐ 37 00 014 62 €		☐ Size 1 (37 01 037)	213 €
				☐ Size 2 (37 02 037)	213 €
				☐ Size 3 (37 03 037)	213 €


Foot cups with leg guidance

Foot cups with leg guidance – width adjustable and with fixing straps & with magnetic lock							
	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	21 cm	26 cm
Size 4	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	28 cm	36 cm
							Price / pair
☐ Size 1 (37 01 038)							361 €
☐ Size 2 (37 02 038)							361 €
☐ Size 3 (37 03 038)							361 €
☐ Size 4 (37 04 038)							361 €



Foot cups with dynamic leg guidance

Foot cups with dynamic leg guidance width adjustable and with fixing straps & with magnetic lock 							
	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	12 cm	15 cm	8,8 cm	11,8 cm	24 cm	20,5 cm	25,5 cm
							Price / pair
☐ Size 1 (37 01 039)							472 €
☐ Size 2 (37 02 039)							472 €
☐ Size 3 (37 03 039)							472 €





dynamic pacer.

The dynamic gait trainer.

Sometimes less is more. Even without full equipment, the **dynamic pacer** already provides optimal support for walking, running or romping. For example, equipped as a light version, only with the new hand grips, the new chest system and the MPS.

Further information can be found from page 80.



momo therapy bike.

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____


Commission:
The therapy bike for more balance.

- Aluminium frame with extra-low access point
- Rim brake with separate parking brake
- Angle adjustable handlebar stem
- Steering damper for stabilisation when moving straight on
- Drum brake in front wheel
- Stabiliser frame with dynamic supports and adjustable spring strength
- Lighting system according to StVZO with hub dynamo (from 24" and not possible with an electric drive)
- Incl. luggage carrier

HMV-No.: 22.51.01.0014

Dimensions

		12"	16"	20"	24"	26"
A*	Pedal to upper edge of saddle	30 - 43 cm	42 - 62 cm	49 - 71 cm	58 - 80 cm	70 - 92 cm
A**	Pedal to upper edge of saddle	33 - 40 cm	45 - 60 cm	52 - 67 cm	61 - 84 cm	73 - 96 cm
B	Saddle support tube to front tube	33 cm	39 cm	45 cm	53 cm	58 cm
C	Wheelbase	64 cm	76 cm	89 cm	104 cm	112 cm
D	Total length	99 cm	122 cm	148 cm	168 cm	178 cm
E	Access point height	16 cm	18 cm	20 cm	22 cm	23 cm
F	Total width	43 cm	52 cm	62 cm	71 cm	76 cm
	Max. load	35 kg	50 kg	65 kg	100 kg	100 kg
	Weight	10,5 kg	13 kg	15 kg	18 kg	20 kg



* Dimensions with standard saddle post | ** Dimensions with T-saddle post

Basic model				
momo therapierad. 12"	momo therapierad. 16"	momo therapierad. 20"	momo therapierad. 24"	momo therapierad. 26"
Basic features, without drum brake, additionally with pennant	Basic features, additionally with pennant	Basic features, additionally with pennant	Basic features	Basic features
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
☐ 38 01 000 1.155 €	☐ 38 02 000 1.295 €	☐ 38 03 000 1.410 €	☐ 38 04 000 1.690 €	☐ 38 05 000 1.690 €

Frame colours				
Colour selection for 12"	Colour selection for 16"+20"	Colour selection for 24"+26"		
  	  	 		
Price	Price	Price		
☐ Telemagenta (Pink) (92 00 012) incl.	☐ Melon yellow (92 00 006) incl.	☐ Anthracit metallic (92 00 014) incl.		
☐ Sky blue (92 00 013) incl.	☐ Marine blue (92 00 008) incl.	☐ Pure white (92 00 100) incl.		
☐ Pure white (92 00 100) incl.	☐ Pure white (92 00 100) incl.			


Configuration possibilities - drive			
	Art.-No.	Drive	Upcharge / piece
	☐ 38 00 001	Rigid sprocket (Standard configuration, if nothing else stated)	incl.
	☐ 38 00 101	Freewheel brake hub (with back pedal brake) only for 16"-24"	incl.
	☐ 38 00 002	Extra light transmission (only for above drive types)	incl.
	☐ 38 00 003	3-gear freewheel brake hub (with back pedal brake)	226 €
	☐ 38 00 007	3-gear freewheel brake hub with reverse option only for 20"-26"	324 €
	☐ 38 00 005	7-Gang-Freewheel brake hub (with back pedal brake)	375 €
	☐ 38 00 008	7-gear-freewheel hub with roll brake (brake lever mounted on the handlebar, without back pedal brake) only for 20"-26"	540 €

Saddle											
Standard saddle (standard configuration, if nothing else stated)			Gel saddle				Saddle seat with leg hollows				
											
	Width	Length		Rear width	Front width	Length		Rear width	Front width	Rear length	Total length
Size 1	15cm	21cm									
Size 2	15cm	24cm	Size 1	19cm	4cm	24cm	Size 1	20,5cm	4cm	9,5cm	14,5cm
Size 3	18cm	26cm	Size 2	24cm	7cm	27cm	Size 2	24cm	4cm	9,5cm	14,5cm
Price / piece			Upcharge / piece				Upcharge / piece				
☐ Size 1 (37 01 001)			☐ Size 1 (37 01 024)				☐ Size 1 (37 01 003)				
☐ Size 2 (37 02 001)			☐ Size 2 (37 02 024)				☐ Size 2 (37 02 003)				
☐ Size 3 (37 03 001)											
incl.			57 €				216 €				
incl.			57 €				216 €				
incl.											

Saddle								Saddle supports				
Unicycle saddle – banana shape (caution: min. inside leg length increases by 2 cm)				Moped saddle (caution: min. inside leg length increases by 2 cm)				Standard saddle support (standard configuration, if nothing else stated)		T-saddle support – for horizontal saddle adjustment		
												
Rear width	Middle width	Front width	Length	Rear width	Front width	Rear length	Total length		for		for	
Size 1	12"			Size 1	12"			Size 1	12"			
Size 2	16"			Size 2	16"+20"			Size 2	16"+20"			
Size 3	20"-26"			Size 3	24"+26"			Size 3	24"+26"			
Price / piece								Upcharge / piece				
☐ Size 1 (37 01 004)				incl.				☐ Size 1 (37 01 005)				83 €
☐ Size 2 (37 02 004)				incl.				☐ Size 2 (37 02 005)				83 €
☐ Size 3 (37 03 004)				incl.				☐ Size 3 (37 03 005)				83 €

Holding brackets
Holding brackets incl. holder – only for use together with the T-saddle support

	max. Height*	Depth
Size 1	20 cm	12 cm
Size 2	30 cm	12 cm
Size 3	37 cm	12 cm
Size 4	53 cm	12 cm
Tailored	____ cm (max. 85 cm)	____ cm (max. 30 cm)


Headrest
Headrest – only for use together with a holding bracket

	Width	Height
Size 1	20 cm	15 cm
Size 2	23 cm	18 cm



Max. height: measured from the standard saddle to the top edge of the respective support pad

Price / set			Price / piece		
☐ Size 1	(37 01 007)	129 €	☐ Size 1	(37 01 029)	343 €
☐ Size 2	(37 02 007)	129 €	☐ Size 2	(37 02 029)	343 €
☐ Size 3	(37 03 007)	129 €			
☐ Size 4	(37 04 007)	129 €			
☐ Tailored	(37 09 007)	272 €			

Torso supports
Dynamic back support pad

padded with strap and bracket (only for use together with a holding bracket)



	Width
Size 1	20 - 30 cm
Size 2	25 - 35 cm
Size 3	30 - 40 cm

Price / piece	
☐ Size 1	(37 01 008) 243 €
☐ Size 2	(37 02 008) 243 €
☐ Size 3	(37 03 008) 243 €

Width adjustable back support pad

padded, incl. bracket (only for use together with a holding bracket)



	Width
Size 1	20 - 28 cm
Size 2	25 - 33 cm
Size 3	30 - 36 cm

Price / piece	
☐ Size 1	(37 01 025) 324 €
☐ Size 2	(37 02 025) 324 €
☐ Size 3	(37 03 025) 324 €

Pelvic supports
Dynamic pelvic guide support pad

padded with strap and bracket (only for use together with a holding bracket)



	Width
Size 1	20 - 30 cm
Size 2	25 - 35 cm
Size 3	30 - 40 cm

Price / piece	
☐ Size 1	(37 01 009) 243 €
☐ Size 2	(37 02 009) 243 €
☐ Size 3	(37 03 009) 243 €

Width adjustable pelvic guide support pad

padded, incl. bracket (only for use together with a holding bracket)



	Width
Size 1	20 - 28 cm
Size 2	25 - 33 cm
Size 3	30 - 36 cm

Price / piece	
☐ Size 1	(37 01 026) 324 €
☐ Size 2	(37 02 026) 324 €
☐ Size 3	(37 03 026) 324 €

Harness system
Chest strap – for width adjustable back support pad

	Width	Length
Size 1	5 cm	18,5 cm
Size 2	6,5 cm	25,5 cm
Size 3	7 cm	30,5 cm



Price / piece	
☐ Size 1	(37 01 050) 98 €
☐ Size 2	(37 02 050) 98 €
☐ Size 3	(37 03 050) 98 €

Positioning jacket – for width adjustable back support pad, incl. universal bracket

	Inside width	Total width	Side height	Total height
Size 1	8 cm	25 cm	10 cm	28 cm
Size 2	9 cm	30 cm	12 cm	30 cm
Size 3	12 cm	34 cm	16 cm	35 cm



Price / piece	
☐ Size 1	(37 01 051) 243 €
☐ Size 2	(37 02 051) 243 €
☐ Size 3	(37 03 051) 243 €

Harness system

T-shape groin harness – for width adjustable pelvic support pad						4-point pelvic harness – for width adjustable pelvic support pad					
	Seat width	Front width	Seat length	Total length			Outer part width	Centre section width	Outer part length	Centre section length	
Size 1	24 cm	23 cm	15 cm	32 cm		Size 1	7 cm	7 cm	12 cm	12 cm	
Size 2	26 cm	26 cm	17 cm	34 cm		Size 2	9 cm	8,5 cm	13 cm	14 cm	
Size 3	28 cm	28 cm	19 cm	38 cm	Size 3	11,5 cm	11 cm	15 cm	16 cm		
					Price / piece						Price / piece
☐ Size 1 (37 01 052)					166 €	☐ Size 1 (37 01 053)					170 €
☐ Size 2 (37 02 052)					166 €	☐ Size 2 (37 02 053)					170 €
☐ Size 3 (37 03 052)					166 €	☐ Size 3 (37 03 053)					170 €

Accessories

Universal bracket for harness used without support pads (width = 27.5 cm)		Hand positioner – for all handlebars (if no quantity was specified: standard = 1 piece)		Steering stop – adjustable up to direction lock (incl. for steering for attendant)		Handlebar extension		
		Front width			for			
	Size 0	8,5cm						
	Size 1	10 cm						
	Size 2	12 cm						
		Size 3	14 cm					
Price / piece		Price / piece		Menge	Price / piece		Price / piece	
☐ 37 02 055	110 €	☐ Size 0 (37 00 054)	70 €	_____	☐ Size 1 (37 01 006)	129 €	☐ Size 1 (37 01 022)	29 €
		☐ Size 1 (37 01 054)	70 €	_____	☐ Size 2 (37 03 006)	129 €	☐ Size 2 (37 02 022)	29 €
		☐ Size 2 (37 02 054)	70 €	_____				
		☐ Size 3 (37 03 054)	70 €	_____				

Handlebar

Universal bracket – black (standard configuration, if nothing else stated)			Touring handlebar			Round handlebar (size 1 is not combinable with e-drive)			Multifunctional handlebar suitable for 20" – 26"		
	Standard for	Width		Width	Depth		Width	Depth		Width	Depth
	Size 1	12"									
	Size 2	16"									
	Size 3	20"									
Size 4	24" + 26"	58 cm	Size 2	58 cm	17 cm	Size 3	48 cm	25 cm	Size 1	61,5 cm	16,5 cm
Price / piece			Upcharge / piece			Upcharge / piece			Upcharge / piece		
☐ Size 1 (37 01 011)			☐ Size 1 (37 01 012)			☐ Size 1 (37 01 013)			☐ 37 02 014		
☐ Size 2 (37 02 011)			☐ Size 2 (37 02 012)			☐ Size 2 (37 02 013)					
☐ Size 3 (37 03 011)						☐ Size 3 (37 03 013)					
☐ Size 4 (37 04 011)											

Handlebar brake

Round handlebar with brake lever ring – not combinable with e-drive			Round handlebar with brake function – the brake function is simply activated by pressing down the round handlebars (not combinable with e-drive)				
	Width	Depth			Width	Depth	
Size 1	40 cm	17 cm		Size 1	40 cm	17 cm	
Size 2	43 cm	21 cm		Size 2	43 cm	21 cm	
Size 3	48 cm	25 cm		Size 3	48 cm	25 cm	
Upcharge / piece				Upcharge / piece			
☐ Size 1 (37 01 028) 329 €				☐ Size 1 (37 01 031) 329 €			
☐ Size 2 (37 02 028) 329 €				☐ Size 2 (37 02 031) 329 €			
☐ Size 3 (37 03 028) 329 €				☐ Size 3 (37 03 031) 329 €			

Accessories

Push bar height adjustable and removable, with integrated twist lock 		Parking brake for attendant (can only be used in connection with a push bar) 		Brake for attendant incl. motor switch-off when using the e-drive and fitting of a disc brake on the front wheel (can only be used in conjunction with Push bar)		Steering for attendant via deflection bar, incl. steering stop, the steering damper is omitted (can only be used in connection with a push bar)	
Price / piece		Price / piece		Price / piece		Price / piece	
☐ for 12"	(38 01 009) 179 €	☐ 38 00 011	119 €	☐ 38 02 011	473 €	☐ for 12"	(38 01 010) 432 €
☐ for 16"	(38 02 009) 179 €					☐ for 16"	(38 02 010) 432 €
☐ for 20"-26"	(38 03 009) 179 €					☐ for 20"	(38 03 010) 432 €

Crank

Progressively adjustable crank reducer for 12"-20" there is no widening of the pedal spacing and for 24"+26" there is a widening of the pedal spacing of approx. 18 mm on each side 		Crank reducer – minimal pedal distance widening, reduction in size by 2,5 or 5 cm respectively 		Special crank for knee contracture for 20" – 26" 	
Price / pair		Price / pair		Price / piece	
☐ for 12"-20"	(37 00 016) 176 €	☐ for 16"+20"	(37 02 015) 91 €	☐ mounted on right	(37 03 015) 377 €
☐ for 24"+26"	(37 04 016) 176 €	☐ for 24"+26"	(37 05 015) 91 €	☐ mounted on left	(37 04 015) 377 €

Pedals

Foot fixing pedals – with racing hooks and compensation weights 		Home exerciser pedals 		Foot cups – width adjustable, with fixing straps & with magnetic lock 	
Depth		Width			
13 cm		12 cm			
Upcharge / pair		Upcharge / pair		Price / pair	
☐ 37 00 013	91 €	☐ 37 00 014	62 €	☐ Size 1 (37 01 037)	213 €
				☐ Size 2 (37 02 037)	213 €
				☐ Size 3 (37 03 037)	213 €


Foot cups with leg guidance

Foot cups with leg guidance – width adjustable and with fixing straps & with magnetic lock							
	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	21 cm	26 cm
Size 4	11,5 cm	14 cm	8,1 cm	10,6 cm	23,8 cm	28 cm	36 cm
Price / pair							
☐ Size 1 (37 01 038)							361 €
☐ Size 2 (37 02 038)							361 €
☐ Size 3 (37 03 038)							361 €
☐ Size 4 (37 04 038)							361 €



Foot cups with dynamic leg guidance

Foot cups with dynamic leg guidance
width adjustable and with fixing straps
& with magnetic lock



	Min. front width	Max. front width	Min. rear width	Max. rear width	Length	Min. leg guide height	Max. leg guide height
Size 1	8,7 cm	10,4 cm	5,7 cm	7,4 cm	17,4 cm	15 cm	18,5 cm
Size 2	9,5 cm	12 cm	6,7 cm	9,2 cm	20,1 cm	18 cm	22 cm
Size 3	12 cm	15 cm	8,8 cm	11,8 cm	24 cm	20,5 cm	25,5 cm

Price / pair

☐ Size 1 (37 01039)	472 €
☐ Size 2 (37 02039)	472 €
☐ Size 3 (37 03039)	472 €

Configuration possibilities - lighting system

	Art.-No.	Lighting system	for	Upcharge / set
	☐ 38 01 021	Upgrade with a lighting system with side runners and reflectors according to StVZO	16"+20"	129 €
	☐ 38 00 021	Upgrading of a battery light system attachable incl. reflectors according to StVZO	12" - 26"	169 €
	☐ 38 02 021	Upgrade with a lighting system with hub dynamo and reflectors according to StVZO (not possible in conjunction with an electric drive)	16"+20"	169 €

Accessories

Basket



Quick connection of stabiliser frame for simple disassembly for transport

	Price / piece		Price / set
☐ for 12"+16" (38 01 019)	57 €	☐ for 12" (38 01 026)	134 €
☐ for 20"+24" (38 03 019)	57 €	☐ for 16" (38 02 026)	134 €
☐ for 26" (38 05 019)	57 €	☐ for 20" (38 03 026)	134 €
		☐ for 24" (38 04 026)	134 €
		☐ for 26" (38 05 026)	134 €

Nachrüstsatz

Stabiliser frame for third-party bicycles – with dynamic supports and adjustable spring strength, for attachment to existing Pletscher plate.



If in doubt, clarify the compatibility of the stabiliser frame with your bicycle by making an enquiry before ordering.

	Price / set
☐ for 12" (38 01 101)	688 €
☐ for 16" (38 02 101)	688 €
☐ for 20" (38 03 101)	688 €
☐ for 24" (38 04 101)	753 €
☐ for 26" (38 05 101)	753 €
☐ Assembly at Schuchmann (38 00 025)	Net 165 €

Universal adapter for bikes without Pletscher plate
for mounting the Stabiliser frame for third-party bicycles that do not have a Pletscher plate



	Price / piece
☐ 38 00 100	69 €

E - drive

DirectPower E-drive from Heinzmann

The new **momo therapy bike**. is equipped with the DirectPower E-drive from Heinzmann. This allows us to offer you a wear-free and silent direct drive. The motor is designed to be brushless and without gears. This minimises maintenance and wear and eliminates motor noise. The recuperative system automatically recharges when driving downhill or braking. This means that the range can be increased by up to 15 %. The existing force of the user is measured by a force sensor integrated in the bottom bracket, and the assistance from the electric drive is adjusted.

Choose the suitable E-drive for the **momo therapy bike**. from the following configurations. You can choose whether the drive should be equipped **with or without start-up assistance**. In addition, you can choose between a **beginner version** and an **advanced version**. The difference here is the **start-up** and the **top speed** of the various modes **Eco**, **Standard** and **Power**, which can be selected via the display. An analysis cable can be used to change the selected version at a later date.



The e-drive excludes a rigid gear rim and brings with it an additional weight of approx. 10 kg.

E - drive configuration

DirectPower E-drive from Heinzmann – without start-up assistance

for Größe	Version	Start-up handling / top speed			Price / piece
		Eco	Standard	Power	
16"	☐ Beginner (38 02 085)	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.886 €
	☐ Advanced (38 02 086)	medium / 6 km/h	strong / 8 km/h	strong / 10 km/h	2.886 €
20"	☐ Beginner (38 03 085)	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	2.886 €
	☐ Advanced (38 03 086)	medium / 6 km/h	strong / 8 km/h	strong / 12 km/h	2.886 €
24"	☐ Beginner (38 04 085)	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.886 €
	☐ Advanced (38 04 086)	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.886 €
26"	☐ Beginner (38 05 085)	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	2.886 €
	☐ Advanced (38 05 086)	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	2.886 €

Konfiguration E - drive

DirectPower E-drive from Heinzmann – with start-up assistance

Select whether the start-up assistance should be activated by pressing the button (**Fig. 1**) or by twist handle (**Fig. 2**).



for size	Version	Start-up assistance	Start-up handling / top speed			Price / piece
			Eco	Standard	Power	
16"	☐ Beginner with twist handle (38 02 087)	4 km/h	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	3.004 €
	☐ Beginner with button (38 02 088)					3.004 €
	☐ Advanced with twist handle (38 02 089)	6 km/h	medium / 6 km/h	strong / 8 km/h	strong / 10 km/h	3.004 €
	☐ Advanced with button (38 02 090)					3.004 €
20"	☐ Beginner with twist handle (38 03 087)	4 km/h	gentle / 4 km/h	medium / 4 km/h	strong / 6 km/h	3.004 €
	☐ Beginner with button (38 03 088)					3.004 €
	☐ Advanced with twist handle (38 03 089)	6 km/h	medium / 6 km/h	strong / 8 km/h	strong / 12 km/h	3.004 €
	☐ Advanced with button (38 03 090)					3.004 €
24"	☐ Beginner with twist handle (38 04 087)	6 km/h	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	3.004 €
	☐ Beginner with button (38 04 088)					3.004 €
	☐ Advanced with twist handle (38 04 089)	6 km/h	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	3.004 €
	☐ Advanced with button (38 04 090)					3.004 €
26"	☐ Beginner with twist handle (38 05 087)	6 km/h	gentle / 6 km/h	medium / 8 km/h	strong / 12 km/h	3.004 €
	☐ Beginner with button (38 05 088)					3.004 €
	☐ Advanced with twist handle (38 05 089)	6 km/h	gentle / 10 km/h	medium / 15 km/h	strong / 20 km/h	3.004 €
	☐ Advanced with button (38 05 090)					3.004 €

Software

NSubsequent conversion/adjustment of the software on the delivered product by our sales force

Net price (not discountable) / piece

☐ **38 00 091**

152 €

Video Library.

Would you like to know more? In our video library we have compiled lots of explanatory videos, product information and entertainment for you. Simply select the desired product category to filter the results. Have fun watching the clips.

General.

Sitting + Positioning.

Walking + Travelling.

Walking + Riding.

Standing + Transfer.

Bathroom + Toilet.

malte. / malte outdoor.



The versatile posterior walking aid **malte.** introduces itself!

marcy.



The anterior walking aid **marcy.** introduces itself.

[To the product](#)

dynamic pacer.



The **dynamic pacer.** does a great job with the new 'smart accessories'.



Good positioning for optimal gait training.



Getting to know the MPS Adjustments.

[To the product](#)

momo motion.



The tried-and-tested tricycle with E-drive **momo motion.** introduces itself.

You can find all our videos on our homepage:
schuchmann.de/en/video-library/



Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:


The height adjustable supine standing frame.

- Height and angle adjustment (0–90°) via gas pressure spring for sizes 1+2
- Height and angle adjustment (0–90°) via electric power for sizes 3+4
(incl. safety shut-off and manual operation)
- Frame shaped to work with mobile hoists
- Upper body support adjustable in length, incl. padding elements
to adjust the length of the upper body support
- Incl. chest strap and pelvic harness
- Adapter for leg rests (standard)
- Leg rests (standard), adjustable in length, with up to
5 cm pelvic rotation adjustment
- Split footplates (standard), which can be lowered to the floor
via height adjustment of the center column
- Base frame with 5 (for size 1–3) or 6 (for size 4) castors (75 mm)
with total locking devices, one of them with directional lock
- Frame colour: Telegrey 4 (white grey), attachment parts: black
- Incl. tool for all settings (adjustments can be
made without tools)

HMV-No.: 28.29.02.1057




Dimensions



		Size 1	Size 2	Size 3	Size 4
A	Body size	75 - 115 cm	100 - 135 cm	120 - 160 cm	145 - 185 cm
B	Upper body support	Hip/thorax width	14 - 27 cm	17 - 31 cm	20 - 37 cm
C		Length	30,5 - 43 cm	38 - 48 cm	49 - 61 cm
D	Leg rests	Foot to knee	20,5 - 30,5 cm	27,5 - 37 cm	32,5 - 45 cm
E		Inside leg to knee	12 - 18,5 cm	16,5 - 24 cm	20 - 29,5 cm
		Knee angle	0° - 65°		
		Hip angle	0° - 45°		
		Abduction (each leg)	0° - 30°		
F		Pelvic rotation adjustment	0 - 5 cm		
G	Feet	Dorsiflexion	0° - 20°		
G		Plantarflexion	0° - 15°		
H		Inversion/Eversion	0° - 15° each direction		
I	Height of lying surface (for transfer)	53,5 - 73,5 cm	60 - 85 cm	71,5 - 93,5 cm	82,5 - 104,5 cm
J	Tilt	0° - 90°			
	Transport dimensions (LxWxH)	95x51x75 cm	105x61x85 cm	125x76x97 cm	152x77x110 cm
	Max. load	30 kg	60 kg	80 kg	100 kg
	Weight	30 kg	35 kg	80 kg	90 kg

Settings



The following measurements should be pre-set on the standing frame:

A Body size	_____	F Pelvic width	_____
B Axilla height (Centre lateral pads)	_____	G Hip height (Centre pelvic pads)	_____
C Shoulder height	_____	H Knee height	_____
D Inside leg height	_____	I Neck width	_____
E Chest width	_____	J Tray height	_____

□ 42 00 020

Price

incl.

Basic model

till. size 1 for users with a body size of 75 – 115 cm		till. size 2 for users with a body size of 100 – 135 cm	
with gas spring supported height and angle adjustment	with gas spring supported angle adjustment and power supported height adjustment (mains supply)	with gas spring supported height and angle adjustment	with gas spring supported angle adjustment and power supported height adjustment (mains supply)
Price / piece	Price / piece	Price / piece	Price / piece
☐ 42 01 000 2.550 €	☐ 42 01 100 4.290 €	☐ 42 02 000 2.750 €	☐ 42 02 100 4.500 €

Basic model

till. size 3 for users with a body size of 120 – 160 cm		till. size 4 for users with a body size of 145 – 185 cm	
with power supported height and angle adjustment	with power supported height and angle adjustment and contourable upper body support (divided into shoulder, trunk and pelvis areas – areas can be individually fixed in the desired position with locking levers)	with power supported height and angle adjustment	with power supported height and angle adjustment and contourable upper body support (divided into shoulder, trunk and pelvis areas – areas can be individually fixed in the desired position with locking levers)
Price / piece	Price / piece	Price / piece	Price / piece
☐ 42 03 000 4.520 €	☐ 42 03 100 5.405 €	☐ 42 04 000 4.745 €	☐ 42 04 100 5.735 €

Battery

 Rechargeable battery (Li-Ion 2.25 Ah) – for basic models with power supported height adjustment (Note: can only be retrofitted by trained personnel)	
	Price / piece
☐ 42 01 026	464 €

Covers for size 1+2

Black / Deep Sea Blue (No. 105 + 111)	Black / Pink 2 (No. 105 + 115)
	
Price	Price
☐ 91 00 239 incl.	☐ 91 00 234 incl.

Covers for size 3+4


Black / Petrol (No. 105 + 112)	Black / Dark Grey (No. 105 + 92)
	
Price	Price
☐ 91 00 237 incl.	☐ 91 00 236 incl.

Headrests

Headrest (Standard) padded, incl. cover, height, depth and angle adjustment	Headrest (Shell-shaped) padded, incl. cover, height, depth and angle adjustment	Headrest (Lateral guide) padded, incl. cover, height, depth and angle adjustment
		
Price / piece	Price / piece	Price / piece
☐ for size 1 (42 01 001) 259 € ☐ for size 2 (42 02 001) 259 € ☐ for size 3+4 (42 03 001) 259 €	☐ for size 1 (42 01 003) 353 € ☐ for size 2 (42 02 003) 353 € ☐ for size 3+4 (42 03 003) 353 €	☐ for size 1 (42 01 002) 353 € ☐ for size 2 (42 02 002) 353 € ☐ for size 3+4 (42 03 002) 353 €

Headrests

Headrest (lateral guide, angle adjustable) – padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable



A = Width contact surface
B = Height headrest
C = Width lateral guide

	A	B	C
Size 1	13,5cm	10,5cm	10cm
Size 2	15cm	10,5cm	13,5cm
Size 3	17cm	12,5cm	16,5cm
Price / piece			
☐ Size 1 (42 01 033)			582 €
☐ Size 2 (42 02 033)			582 €
☐ Size 3 (42 03 033)			582 €


Universal adapter for tailored headrests
15x15 mm, height and depth adjustable




Price / piece		Price / piece
☐ 42 01 004		
	199 €	

Lateral pads

Lateral pads for lying surface (standard)
incl. upper arm support, flip away, stepless adjustment for height and width, adjustable for rotation



Pelvic pads for lying surface (standard)
incl. forearm rest, flip away, stepless adjustment for height and width, adjustable for rotation



	Height (A)	Depth (C)	suitable for	Chest width (B)		Height (A)	Depth (C)	suitable for	Chest width (B)
Size 1	10,5cm	9cm	till. Size 1	14 – 27cm	Size 1	12cm	9cm	till. Size 1	14 – 27cm
			till. Size 2	17 – 31cm				till. Size 2	17 – 31cm
Size 2	11,5cm	11cm	till. Size 2	17 – 31cm	Size 2	14cm	10,5cm	till. Size 2	17 – 31cm
			till. Size 3	20 – 37cm				till. Size 3	20 – 37cm
Size 3	12,5cm	16cm	till. Size 3	20 – 37cm	Size 3	17cm	15,5cm	till. Size 3	20 – 37cm
			till. Size 4	29 – 46cm				till. Size 4	29 – 46cm
Size 4	13,5cm	18cm	till. Size 3	20 – 37cm	Size 4	21cm	17,5cm	till. Size 3	20 – 37cm
			till. Size 4	27 – 46cm				till. Size 4	27 – 46cm



The **standard version** of the lateral and pelvic pads is only suitable for the **till. without** **contourable upper body support!**
The **depth adjustable version** of the lateral and pelvic pads is only suitable for the **till. (size 3 & 4) with** **contourable upper body support!**

Price / pair		Price / pair		Price / pair		Price / pair	
☐ Size 1 (42 01 005)	293 €	☐ Size 1, Depth adjustable (42 01 006)	293 €	☐ Size 1 (42 01 007)	293 €	☐ Size 1, Depth adjustable (42 01 008)	293 €
☐ Size 2 (42 02 005)	293 €	☐ Size 2, Depth adjustable (42 02 006)	293 €	☐ Size 2 (42 02 007)	293 €	☐ Size 2, Depth adjustable (42 02 008)	293 €
☐ Size 3 (42 03 005)	293 €	☐ Size 3, Depth adjustable (42 03 006)	293 €	☐ Size 3 (42 03 007)	293 €	☐ Size 3, Depth adjustable (42 03 008)	293 €
☐ Size 4 (42 04 005)	293 €	☐ Size 4, Depth adjustable (42 04 006)	293 €	☐ Size 4 (42 04 007)	293 €	☐ Size 4, Depth adjustable (42 04 008)	293 €

Support pads

Knee support pads with patella recess
individual three-dimensional adjustment, with positioning aid



	Inside width (A)
Size 1	8cm
Size 2	10cm
Size 3	12cm

Price / pair

☐ Size 1 (42 01 035)	519 €
☐ Size 2 (42 02 035)	519 €
☐ Size 3 (42 03 035)	519 €

Flexible knee positioning (FKP) – Individually three-dimensionally adjustable, with length-adjustable contact straps that adapt to the shape of the knee



	Inside width (A)
Size 1	8cm
Size 2	10cm
Size 3	12cm



Price / pair

☐ Size 1 (42 01 036)	585 €
☐ Size 2 (42 02 036)	585 €
☐ Size 3 (42 03 036)	585 €

Padding – for additional support of the knee, to prevent hyperextension (knee overextended)



	Length	Ø
Size 1	8cm	4cm
Size 2	8cm	6cm
Size 3	10cm	8cm
Size 4	12cm	10cm

Price / pair

☐ Size 1 (42 01 011)	161 €
☐ Size 2 (42 02 011)	161 €
☐ Size 3 (42 03 011)	161 €
☐ Size 4 (42 04 011)	161 €

for size 1+2	for size 3+4	Push bar	Safety bolt – for release lever, prevents unintentional release of height and angle adjustments (Note: only for size 1+2)									
<table border="1"> <thead> <tr> <th colspan="2">Price / piece</th> </tr> </thead> <tbody> <tr> <td>☐ for size 1+2 (42 01 012)</td> <td>309 €</td> </tr> <tr> <td>☐ for size 3+4 (42 03 012)</td> <td>309 €</td> </tr> </tbody> </table>			Price / piece		☐ for size 1+2 (42 01 012)	309 €	☐ for size 3+4 (42 03 012)	309 €	<table border="1"> <thead> <tr> <th>Price / piece</th> </tr> </thead> <tbody> <tr> <td>☐ 42 00 013</td> </tr> <tr> <td>75 €</td> </tr> </tbody> </table>	Price / piece	☐ 42 00 013	75 €
Price / piece												
☐ for size 1+2 (42 01 012)	309 €											
☐ for size 3+4 (42 03 012)	309 €											
Price / piece												
☐ 42 00 013												
75 €												

Adapter for leg rests (standard) (standard configuration, if nothing else stated)	Adapter for leg rests (abduction) enables individual and continuous abduction of each leg from 0°–30° (Note: not retrofittable!)	Leg rests (standard) – length adjustable, with 5 cm pelvic rotation adjustment (standard configuration if nothing else stated)	Leg rests (flex) – adjustable in length and angle, with 5 cm pelvic rotation adjustment (Note: not retrofittable!)
Price	Upcharge	Price / pair	Upcharge / pair
☐ for size 1+2 (4201014) incl. ☐ for size 3+4 (4203014) incl.	☐ for size 1+2 (4201015) 794 € ☐ for size 3+4 (4203015) 794 €	☐ for size 1 (4201016) incl. ☐ for size 2 (4202016) incl. ☐ for size 3 (4203016) incl. ☐ for size 4 (4204016) incl.	☐ for size 1 (4201017) 480 € ☐ for size 2 (4202017) 480 € ☐ for size 3 (4203017) 480 € ☐ for size 4 (4204017) 480 €

Split footplates (standard) standard configuration if nothing else stated (Note: not retrofittable!)	Price / pair <table> <tr> <td>☐ for size 1 (42 01018)</td><td>incl.</td></tr> <tr> <td>☐ for size 2 (42 02018)</td><td>incl.</td></tr> <tr> <td>☐ for size 3 (42 03018)</td><td>incl.</td></tr> <tr> <td>☐ for size 4 (42 04018)</td><td>incl.</td></tr> </table>	☐ for size 1 (42 01018)	incl.	☐ for size 2 (42 02018)	incl.	☐ for size 3 (42 03018)	incl.	☐ for size 4 (42 04018)	incl.
☐ for size 1 (42 01018)	incl.								
☐ for size 2 (42 02018)	incl.								
☐ for size 3 (42 03018)	incl.								
☐ for size 4 (42 04018)	incl.								
Split footplates (flex) – individually adjustable, dorsiflex (20°), plantarflex (15°) and inversion/eversion (+/-15°) (Note: not retrofittable!)	Upcharge / pair <table> <tr> <td>☐ for size 1 (42 01019)</td><td>243 €</td></tr> <tr> <td>☐ for size 2 (42 02019)</td><td>243 €</td></tr> <tr> <td>☐ for size 3 (42 03019)</td><td>243 €</td></tr> <tr> <td>☐ for size 4 (42 04019)</td><td>243 €</td></tr> </table>	☐ for size 1 (42 01019)	243 €	☐ for size 2 (42 02019)	243 €	☐ for size 3 (42 03019)	243 €	☐ for size 4 (42 04019)	243 €
☐ for size 1 (42 01019)	243 €								
☐ for size 2 (42 02019)	243 €								
☐ for size 3 (42 03019)	243 €								
☐ for size 4 (42 04019)	243 €								

	Heel edges – incl. belt + magnetic lock (Note: cannot be combined with simple foot strap!)	
		Price / pair
	☐ for size 1 (42 01 021)	176 €
	☐ for size 2 (42 02 021)	176 €
	☐ for size 3 (42 03 021)	176 €
	☐ for size 4 (42 04 021)	176 €

Wood therapy tray (standard) – with edge, removable, depth, height and angle adjustable		Wood therapy tray (with quick adjustment) – with edge, removable, depth, height and angle adjustable via table linkage and additionally via quick adjustment		Padded tray cover – covered with black wipe clean material, height: 1 cm	
Price / piece		Price / piece		Price / piece	
☐ for size 1 (42 01 024)	750 €	☐ for size 1 (42 01 034)	905 €	☐ for size 1 (42 01 025)	157 €
☐ for size 2 (42 02 024)	750 €	☐ for size 2 (42 02 034)	905 €	☐ for size 2 (42 02 025)	157 €
☐ for size 3 (42 03 024)	750 €	☐ for size 3 (42 03 034)	905 €	☐ for size 3 (42 03 025)	157 €
☐ for size 4 (42 04 024)	750 €	☐ for size 4 (42 04 034)	905 €	☐ for size 4 (42 04 025)	157 €

Harness

Chest strap – with pad and adjustable clip on both sides, colour: black (standard configuration if nothing else stated),



	Size 1	Size 2	Size 3
Pad width	4 cm	6 cm	6 cm
Pad length per side	18 cm	25 cm	30 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2,5 cm	2,5 cm	2,5 cm

Price / piece

- ☐ Size 1 **(42 01 029)** incl.
☐ Size 2 **(42 02 029)** incl.
☐ Size 3 **(42 03 029)** incl.

Wide Chest strap – padded, colour: black



	Size 1	Size 2	Size 3
Pad width	8 cm	10 cm	12 cm
PolsterLength	10 cm	25 cm	35 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2,5 cm	2,5 cm	2,5 cm


Upcharge / piece

- ☐ Size 1 **(42 01 030)** 56 €
☐ Size 2 **(42 02 030)** 56 €
☐ Size 3 **(42 03 030)** 56 €

Harness

Positioning jacket – made of elastic, breathable and moisture repellent material, shoulder straps and pelvis straps (25 mm) with clips, colour: black

	A	B	C	D
Size 1	14 cm	26,5 cm	14 cm	33 cm
Size 2	15 cm	29,5 cm	16 cm	36 cm
Size 3	17 cm	35,5 cm	19 cm	41 cm
Size 4	18 cm	38,5 cm	21 cm	43,5 cm



Extension for positioning jacket – for the belt attachment of the positioning jacket

recommended from a neck width of

for size 1	13 cm
for size 2	17 cm
for size 3	21 cm
for size 4	30 cm



Upcharge / piece

Price / set

- ☐ Size 1 **(42 01 027)** 227 €
☐ Size 2 **(42 02 027)** 227 €
☐ Size 3 **(42 03 027)** 227 €
☐ Size 4 **(42 04 027)** 227 €

- ☐ for size 1+2 **(42 01 028)** 54 €
☐ for size 3+4 **(42 03 028)** 54 €

Harness

Pelvic harness – with pad and adjustable clip on both sides (standard configuration if nothing else stated), colour: black



	Size 1	Size 2	Size 3
for pelvic width	16 – 25 cm	25 – 40 cm	38 – 46 cm
Pad width	4 cm	6 cm	6 cm
Pad length per side	18 cm	25 cm	30 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2,5 cm	2,5 cm	2,5 cm

Price / piece

- ☐ Size 1 **(42 01 031)** incl.
☐ Size 2 **(42 02 031)** incl.
☐ Size 3 **(42 03 031)** incl.

Pelvic harness – by „Pelvi.Loc“, colour: black



	Size 1 (XS)	Size 2 (S)	Size 3 (M)	Size 4 (L)
for Pelvic width	16 – 25 cm	25 – 34 cm	32 – 40 cm	38 – 46 cm
Pad width	6,5 cm	7 cm	7,5 cm	8,5 cm
Pad length per side	16 cm	19 cm	22 cm	24 cm
Strap length per side	25 cm	25 cm	45 cm	45 cm
Strap width	2,5 cm	2,5 cm	2,5 cm	2,5 cm

Upcharge / piece

- ☐ Size 1 (XS) **(42 01 032)** 293 €
☐ Size 2 (S) **(42 02 032)** 293 €
☐ Size 3 (M) **(42 03 032)** 293 €
☐ Size 4 (L) **(42 04 032)** 293 €

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The sporty and agile standing frame.


- **Basic frame (standard)** with height and depth adjustable footplate and four height adjustable support castors (75 mm), rear with springs and lock
- **Basic frame (with abduction)** with individual adjustable split footplate, two height adjustable support castors (75 mm) in the front and one support castor (40 mm) in the back
- Incl. chest support pad, rear support pad and three-dimensionally adjustable knee support pads with patella recess
- Mechanical angle adjustment of the central column from 0 - 15°
- Incl. drum brakes, easy operation via locking lever
- incl. tray support and storage net (for size 3+4)
- Adjustable camber



HMV-No.: 28.29.01.3011

Please send us the following measurements for optimum presetting of the standing frame

Settings



The following measurements should be pre-set on the standing frame:

A Body size _____

B Sternum height _____

C Elbow height _____

D Crotch height _____

H Knee height _____

E Chest width _____

F Pelvic width _____


G Knee width _____

☐ **49 00 020**

Price

incl.

Dimensions



		Size 1		Size 2			Size 3		Size 4		
	Body size	75 – 110 cm		100 – 130 cm			120 – 155 cm		150 – 175 cm		
	Wheel size	24"	28"	30"	32"	36"	32"	36"	36"	40"	42"
A	Grip height: footplate – hand rim at top	40 – 48 cm	50 – 58 cm	57 – 65 cm	64 – 69 cm	72 – 80 cm	60 – 68 cm	71 – 79 cm	82 – 90 cm	87 – 95 cm	87 – 95 cm
	Recommended from elbow height	50 cm	60 cm	70 cm	80 cm	92 cm	82 cm	92 cm	92 cm	97 cm	99 cm
B	Width of hand rim at top	0°	-	-	-	-	65 cm	65 cm	67 cm	67 cm	67 cm
		3°	-	-	54 cm	53 cm	61 cm	60 cm	63 cm	62 cm	61 cm
		6°	-	-	51 cm	49 cm	57 cm	55 cm	59 cm	57 cm	56 cm
		9°	45 cm	44 cm	48 cm	45 cm	54 cm	51 cm	55 cm	53 cm	51 cm
		12°	42 cm	41 cm	45 cm	41 cm	-	-	-	-	-
		15°	40 cm	-	-	-	-	-	-	-	-
C	Width of hand rim at the bottom (min. – max.)	68 – 74 cm	70 – 73 cm	71 – 83 cm	71 – 83 cm	71 – 84 cm	73 – 83 cm	73 – 85 cm	76 – 86 cm	76 – 86 cm	76 – 90 cm
D	Height from footplate to centre of knee support pads	19 – 30 cm		28 – 36 cm			34 – 44 cm		42 – 51 cm		
E	Height from footplate to support of buttocks support pad	26 – 42 cm		43 – 58 cm			53 – 69 cm		65 – 80 cm		
	Recommended for inside leg length	30 – 47 cm		48 – 63 cm			58 – 74 cm		70 – 84 cm		
F	Abduction split footplate	0° – 6°		0° – 6°			0° – 30°		0° – 30°		
F	Abduction split footplate (with abducted basic model)	0° – 40°		0° – 40°			-		-		
G	Pelvis dept*	11 – 24 cm		11 – 24 cm			15 – 28,5 cm		15 – 28,5 cm		
	Pelvis dept** (adjustment range of pelvis position via crank max. 6 cm)	11 – 20 cm (+6 cm)		11 – 20 cm (+6 cm)			15 – 27 cm (+6 cm)		15 – 27 cm (+6 cm)		
H	Sternum height (height from footplate to centre of chest support pad)	54 – 80 cm		68 – 101 cm			92 – 123 cm		111 – 140 cm		
I	Tray height from footplate	43 – 66 cm		61 – 79 cm			72 – 96 cm		92 – 110 cm		
	Max. load	30 kg		45 kg			65 kg		80 kg		
	Weight	28 kg		31 kg			38 kg		42 kg		

* with bracket for buttocks support pad (standard) ** with bracket for buttocks support pad (flip away)

Basic model

todd. size 1 for users with a body size 75 – 110 cm		todd. size 2 for users with a body size 100 – 130 cm		todd. size 3 for users with a body size 120 – 155 cm		todd. size 4 for users with a body size 150 – 175 cm	
Price / piece		Price / piece		Price / piece		Price / piece	
☐ Standard	(49 01 000) 2.795 €	☐ Standard	(49 02 000) 2.930 €	☐ Standard	(49 03 000) 3.520 €	☐ Standard	(49 04 000) 3.975 €
☐ Abductable	(49 01 100) 3.135 €	☐ Abductable	(49 02 100) 3.270 €				

Frame colours

Dormant red	Melon yellow	Turquoise blue	Lasurit blue metallic	Anthracit metallic	Pure white
					
Price	Price	Price	Price	Price	Price
☐ 92 00 003 incl.	☐ 92 00 006 incl.	☐ 92 00 016 incl.	☐ 92 00 010 incl.	☐ 92 00 014 incl.	☐ 92 00 100 incl.

Coated covers

Flexi Black		Flexi light grey	
	Price		Price
☐ 91 00 043	incl.	☐ 91 00 050	incl.


Support pads

Chest system – adjustable in height			Chest support pad with lateral guide can be adjusted in height, width and depth			Chest support pad with lateral guide can be adjusted in height, width and depth, additionally angle adjustable			Back strap – with magnetic lock for chest support pad with lateral guide		
											
	Inside width	Height		Chest width			Chest width			Pad width	Pad height
Size 1	27 cm	13 cm	Size 1	16 – 24 cm		Size 1	16 – 24 cm		Size 1	17 cm	8,5 cm
Size 2	36 cm	16 cm	Size 2	18 – 26 cm		Size 2	18 – 26 cm		Size 2	28 cm	10 cm
Size 3	43 cm	19 cm	Size 3	24 – 34 cm		Size 3	24 – 34 cm		Size 3	39 cm	11,5 cm
			Size 4	30 – 42 cm		Size 4	30 – 42 cm				
Price / piece			Upcharge / piece			Upcharge / piece			Price / piece		
☐ Size 1 (49 01 001)		incl.	☐ Size 1 (49 01 002)		342 €	☐ Size 1 (49 01 003)		342 €	☐ Size 1 (49 01 024)		145 €
☐ Size 2 (49 02 001)		incl.	☐ Size 2 (49 02 002)		342 €	☐ Size 2 (49 02 003)		342 €	☐ Size 2 (49 02 024)		145 €
☐ Size 3 (49 03 001)		incl.	☐ Size 3 (49 03 002)		342 €	☐ Size 3 (49 03 003)		342 €	☐ Size 3 (49 03 024)		145 €
			☐ Size 4 (49 04 002)		342 €	☐ Size 4 (49 04 003)		342 €			

Accessories

Back support pad – can be adjusted in height and angle			Headrest – can be adjusted in height, depth and angle, only for use together with the back support pad		
	Width	Height			
for size 1+2	20 cm	23 cm			
for size 3+4	25 cm	28 cm			
	Price / piece		Price / piece		
☐ for size 1 (49 01 025)		404 €	☐ for size 1 (49 01 026)		311 €
☐ for size 2 (49 02 025)		404 €	☐ for size 2 (49 02 026)		311 €
☐ for size 3+4 (49 03 025)		404 €	☐ for size 3 (49 03 026)		311 €
			☐ for size 4 (49 04 026)		311 €

Configuration options - knee pads

Knee support pads with patella recess individual three-dimensional adjustment			Flexible knee positioning (FKP) Individually three-dimensionally adjustable, with length-adjustable contact straps that adapt to the shape of the knee			Knee pad support – 2 different options are available to accommodate the adjacent knee pads		
								
	Inside width			Inside width			Rigid	
Size 1	8 cm		Size 1	8 cm			angle, depth and height adjustable	
Size 2	10 cm		Size 2	10 cm			Ball head	
Size 3	12 cm		Size 3	12 cm			angle-, depth-, height- and additionally tilt angle adjustable	
							(Rigid = Standard configuration unless otherwise specified)	
Please select one of the adjacent knee pads.								
Price / pair			Upcharge / pair			Upcharge / pair		
☐ Size 1 (49 01 031)		incl.	☐ Size 1 (49 01 032)		109 €	☐ Rigid (49 00 030)		incl.
☐ Size 2 (49 02 031)		incl.	☐ Size 2 (49 02 032)		109 €	☐ Ball head (49 01 030)		143 €
☐ Size 3 (49 03 031)		incl.	☐ Size 3 (49 03 032)		109 €			

FKP.
flexible knee positioning

Configuration possibilities - Bracket for buttock and pelvic positioning

Configuration 1:



Bracket for buttocks support pad (standard)

depth adjustable via quick adjustment

	Depth adjustment Pelvic pad up to buttock pad	Standard configuration
Size 1	11 – 24 cm	for size 1+2
Size 2	15 – 28,5 cm	for size 3+4
		Price / piece
☐ Size 1 (4901005)		incl.
☐ Size 2 (4903005)		incl.

Configuration 2:




Bracket for buttocks support pad (flip away)

for step-by-step adjustment in the depth via hand crank and quick adjustment, fold-down mechanism for easier transfer (approx. 20 ° (10 cm) downwards)

	Depth adjustment Pelvic pad up to buttock pad	Depth adjustment via Crank	Standard configuration
Size 1	11 – 20 cm	6 cm	for size 1+2
Size 2	15 – 27 cm	6 cm	for size 3+4
			Upcharge / set
☐ Size 1 (4901004)			370 €
☐ Size 2 (4903004)			370 €

Configuration 3:



Pelvic frame (with swivel arm) – variable depth adjustment, additional depth adjustment by hand crank for gradual and graded hip extension of the user, incl. holder for rear support pad and back support pad

 Please consider the combination options below and select a pad configuration on page 120.

	Width outside	Width inside	Depth adjustment Pelvic pad up to buttock pad	Depth adjustment via Crank
Size 1	37,5 cm	31 cm	11 – 22 cm	6 cm
Size 2	47,5 cm	41 cm	14 – 28 cm	6 cm
Size 3	57,5 cm	51 cm	18 – 35 cm	6 cm
				Upcharge / set
☐ Size 1 (4901029)				490 €
☐ Size 2 (4902029)				490 €
☐ Size 3 (4903029)				490 €

Combination options

The pelvic frame (with swivel arm) is available in three different sizes. Please take into account that not every camber setting is possible with every size of pelvic frame, as a collision with the drive wheels may occur above a certain camber setting. The following overview shows which size of pelvic frame can be combined with which basic frame.


	Camber	todd. (size 1)		todd. (size 2)			todd. (size 3)		todd. (size 4)		
Wheel size		24"	28"	30"	32"	36"	32"	36"	36"	40"	42"
	0°						● □ ▲	● □ ▲	● □ ▲	● □ ▲	● □ ▲
	3°			● □	● □	● □	● □	● □	● □	● □	● □
	6°			●	●	●	● □	● □	● □	● □	● □
	9°	●	●	●	●	●	● □	●	● □	● □	●
	12°			●							
	15°										

● = Pelvic frame (size 1), □ = Pelvic frame (size 2), ▲ = Pelvic frame (size 3)

Configuration possibilities - Pelvic positioning

Option 1:

Pelvis support pad with lateral guide (mounted at the front of the pelvis frame) + buttocks support pad (mounted at the back of the pelvis frame)



Pelvis support pad with lateral guide

width and depth adjustable



	Pelvic width	Depth adjustment
Size 1	16 – 24 cm	5 cm
Size 2	18 – 26 cm	5 cm
Size 3	24 – 34 cm	5 cm
Size 4	30 – 42 cm	5 cm
		Price / piece
☐ Size 1 (49 01 006)		381 €
☐ Size 2 (49 02 106)		381 €
☐ Size 3 (49 02 006)		381 €
☐ Size 4 (49 03 006)		381 €


Buttocks support pad



	Width x Height
Size 1	16x10 cm
Size 2	20x10 cm
Size 3	25x14 cm
Price / piece	
☐ Size 1 (49 01 007)	incl.
☐ Size 2 (49 02 007)	incl.
☐ Size 3 (49 03 007)	incl.

Option 2:

Pelvic positioner (mounted at the front of the pelvis frame) + pelvis support pad with lateral guide (mounted at the back of the pelvis frame)



Pelvic positioner – depth adjustable



	Width x Height	Depth adjustment
Size 1	16 x 10 cm	5 cm
Size 2	20 x 10 cm	5 cm
Size 3	25 x 14 cm	5 cm
		Price / piece
☐ Size 1 (49 01 008)		145 €
☐ Size 2 (49 02 008)		145 €
☐ Size 3 (49 03 008)		145 €

Pelvis support pad with lateral guide


width and depth adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Upcharge / piece	
☐ Size 1 (4901010)	235 €
☐ Size 2 (4902110)	235 €
☐ Size 3 (4902010)	235 €
☐ Size 4 (4903010)	235 €

Option 2:

Pelvic positioner (mounted at the front of the pelvis frame) + pelvis support pad with lateral guide (mounted at the back of the pelvis frame)



Pelvic positioner

depth adjustable



	Width x Height	Depth adjustment
Size 1	16 x 10 cm	5 cm
Size 2	20 x 10 cm	5 cm
Size 3	25 x 14 cm	5 cm
		Price / piece
☐ Size 1 (49 01 008)		145 €
☐ Size 2 (49 02 008)		145 €
☐ Size 3 (49 03 008)		145 €

Combined spine and pelvic support pads – width, depth and height adjustable



	Pelvic width	Depth adjustment
Size 1	14 – 24 cm	5 cm
Size 2	22 – 32 cm	5 cm
Size 3	32 – 42 cm	5 cm
		Price / pair
☐ Size 1 (49 01 009)		381 €
☐ Size 2 (49 02 009)		381 €
☐ Size 3 (49 03 009)		381 €

Buttocks support pad



	Width x Height
Size 1	16x10 cm
Size 2	20x10 cm
Size 3	25x14 cm
Price / piece	
☐ Size 1 (49 01 007)	incl.
☐ Size 2 (49 02 007)	incl.
☐ Size 3 (49 03 007)	incl.

Accessories



Therapy tray – made of wood, removable, with edge, incl. height adjustable attachment for chest support pad

	Width	Depth
for size 3+4	70 cm	58 cm
Price / piece		
☐ for size 3+4 (4903028)		370 €

for size 1+2



for size 3+4



Therapy tray with edge – made of Perspex, removable, not usable in combination with basket

	Width	Depth
for size 1+2	52 cm	40 cm
for size 3+4	62 cm	50 cm
Price / piece		
☐ for size 1+2 (4901013)		370 €
☐ for size 3+4 (4903027)		370 €

Accessories



Basket – removable, not usable in combination with therapy tray

Price / piece	
☐ for size 1+2 (4900014)	162 €



Mounting position of the brake – the brake lever can be mounted under the tray on the left or right, as desired (standard configuration = under the tray on the right)

Price	
☐ under tray on left (4901015)	incl.
☐ under tray on right (4902015)	incl.

Rear wheels

Drive wheels – with drum brake, push rims and black PU tyres



Drive wheels with closed hand rim with drum brake and PU tyres



for 24" and 28" = black PU tyres
for 30" to 40" = grey PU tyres

Price / pair for original equipment			Price / pair for retrofitting		
24" for size 1	☐ 4900016	415 €	☐ 4900016ET	930 €	
28" for size 1	☐ 4901016	415 €	☐ 4901016ET	930 €	
30" for size 2	☐ 4902016	496 €	☐ 4902016ET	1.409 €	
32" for size 2+3	☐ 4903016	496 €	☐ 4903016ET	1.409 €	
36" for size 2-4	☐ 4904016	496 €	☐ 4904016ET	1.409 €	
40" for size 4	☐ 4905016	623 €	☐ 4905016ET	1.659 €	
42" for size 4	☐ 4906016	623 €	☐ 4906016ET	1.659 €	

Rear wheels



Spoke guard – transparent

Price / pair	
☐ 24" for size 1 (4900018)	214 €
☐ 28" for size 1 (4901018)	214 €
☐ 30" for size 2 (4902018)	214 €
☐ 32" for size 2+3 (4903018)	311 €
☐ 36" for size 3 (4904018)	311 €
☐ 40" for size 4 (4905018)	381 €



Wheel camber adjustment to set the hand rim width

Price	
☐ 0° for size 3+4 (4900019)	incl.
☐ 3° for size 2-4 (4901019)	incl.
☐ 6° for size 2-4 (4902019)	incl.
☐ 9° for size 1-4 (4903019)	incl.
☐ 12° for size 1+2 (4904019)	incl.
☐ 15° for size 1 (4905019)	incl.

Accessories

Footplate – height and depth adjustable



Split footplate – individual height, depth, width and angle adjustable

already included in basic model „with abduction“



Heel edges with ankle straps – angle, width and depth adjustable, for mounting on the footplate(s) with magnetic locks



Price / piece		Upcharge / pair		Price / pair	
☐ for size 1 (4901020)	incl.	☐ for size 1 (4901021)	392 €	☐ Size 1 (4901027)	182 €
☐ for size 2 (4902020)	incl.	☐ for size 2 (4902021)	392 €	☐ Size 2 (4902027)	182 €
☐ for size 3 (4903020)	incl.	☐ for size 3+4 (4903021)	392 €		
☐ for size 4 (4904020)	incl.				

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The standing frame – from prone to supine.


- Basic frame with height and depth adjustable footplate
- Angular adjustment of the central column via gas pressure spring from (-)15° – 30°
- Incl. chest support pad with lateral guide and buttocks support pad
- Height and depth adjustable pelvis frame, pelvis position adjustable via a crank
- Incl. three-dimensionally adjustable knee support pads with patella cutout
- Incl. therapy tray
- 4 castors with locking devices (75 mm) and total lock



HMV-No.: 28.29.01.1059

Please send us the following measurements for optimum presetting of the standing frame

Settings



The following measurements should be pre-set on the standing frame:

- | | |
|-------------------------------|-----------------------------|
| A Body size _____ | E Chest width _____ |
| B Sternum height _____ | F Pelvic width _____ |
| C Elbow height _____ | G Knee width _____ |
| D Crotch height _____ | H Knee height _____ |


☐ **48 00 020**

Price

incl.

Dimensions

		Size 1	Size 2	Size 3
	Body size	75 - 110 cm	100 - 140 cm	120 - 160 cm
A	Height from footplate to tray	48.5 - 80.5 cm	71 - 95 cm	80 - 112 cm
B	Pelvis depth	12 - 22 cm	16 - 28 cm	20 - 35 cm
	Depth adjustment of pelvis frame (via crank)	6 cm	6 cm	6 cm
C	Total length	84 cm	90 cm	100 cm
D	Total width	57 cm	65 cm	71 cm
E	Abduction divided footplate	0° - 27°	0° - 27°	0° - 27°
F	Pelvis height	45 - 60 cm	57 - 80 cm	65 - 90 cm
G	Sternum height	55 - 85 cm	75 - 110 cm	95 - 125 cm
	Max. load	40 kg	60 kg	80 kg
	Weight	27 kg	29 kg	32 kg



Basic model

tim. Size 1 For users with a body size 75 - 110 cm		tim. Size 2 For users with a body size 100 - 140 cm		tim. Size 3 For users with a body size 120 - 160 cm	
Price / piece		Price / piece		Price / piece	
☐ 48 01 100	3.300 €	☐ 48 02 100	3.300 €	☐ 48 03 100	3.570 €

Frame colours

Colour selection for size 1	Colour selection for size 2	Colour selection for size 3	Coated with an individual colour (on request) Please provide the requested colour here: RAL _____
			
Price	Price	Price	Upcharge
☐ Melon yellow (92 00 006) incl.	☐ Dormant Red (92 00 003) incl.	☐ Lasurit blue metallic (92 00 010) incl.	☐ 48 00 090 210 €

Coated covers

Flexi Black	Flexi light grey
	
Price	Price
☐ 91 00 043 incl.	☐ 91 00 050 incl.


Lateral pads

Chest support pad with lateral guide can be adjusted in height, width and depth		Chest support pad with lateral guide can be adjusted in height, width and depth, additionally angle adjustable		Chin pad – height adjustable		Back strap – with magnetic lock	
							
Chest width		Chest width		Inside width		Pad width	Pad height
Size 1	16 - 24 cm	Size 1	16 - 24 cm	Size 1	6 cm	Size 1	17 cm
Size 2	18 - 26 cm	Size 2	18 - 26 cm	Size 2	8 cm	Size 2	28 cm
Size 3	24 - 34 cm	Size 3	24 - 34 cm	Size 3	10 cm	Size 3	39 cm
Size 4	30 - 42 cm	Size 4	30 - 42 cm				
Price / piece		Upcharge / piece		Price / piece		Price / piece	
☐ Size 1 (48 01 101) incl.		☐ Size 1 (48 01 102) 109 €		☐ Size 1 (48 01 012) 243 €		☐ Size 1 (48 01 013) 145 €	
☐ Size 2 (48 02 201) incl.		☐ Size 2 (48 02 202) 109 €		☐ Size 2 (48 02 012) 243 €		☐ Size 2 (48 02 013) 145 €	
☐ Size 3 (48 02 101) incl.		☐ Size 3 (48 02 102) 109 €		☐ Size 3 (48 03 012) 243 €		☐ Size 3 (48 03 013) 145 €	
☐ Size 4 (48 03 101) incl.		☐ Size 4 (48 03 102) 109 €					

Combinations - pelvic frame

Option 1:

Pelvis support pad with lateral guide (mounted at the front of the pelvic frame) + buttocks support pad (mounted at the back of the pelvic frame)



Pelvis support pad with lateral guide width adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Price / piece	
☐ Size 1 (48 01 005)	381 €
☐ Size 2 (48 02 105)	381 €
☐ Size 3 (48 02 005)	381 €
☐ Size 4 (48 03 005)	381 €


Buttocks support pad



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
☐ Size 1 (48 01 003)	incl.
☐ Size 2 (48 02 003)	incl.
☐ Size 3 (48 03 003)	incl.

Option 2:

Pelvic positioner (mounted at the front of the pelvic frame) + pelvis support pad with lateral guide (mounted at the back of the pelvic frame)



Pelvic positioner



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
☐ Size 1 (48 01 004)	145 €
☐ Size 2 (48 02 004)	145 €
☐ Size 3 (48 03 004)	145 €


Pelvis support pad with lateral guide width adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Upcharge / piece	
☐ Size 1 (48 01 006)	235 €
☐ Size 2 (48 02 106)	235 €
☐ Size 3 (48 02 006)	235 €
☐ Size 4 (48 03 006)	235 €

Option 3:

Pelvic positioner and combined spine and pelvic support pads (mounted at the front of the pelvic frame) + buttocks support pad (mounted at the back of the pelvic frame)



Pelvic positioner



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
☐ Size 1 (48 01 004)	145 €
☐ Size 2 (48 02 004)	145 €
☐ Size 3 (48 03 004)	145 €

Combined spine and pelvic support pads – width, depth and height adjustable



	Pelvic width
for size 1	14 – 24 cm
for size 2	22 – 32 cm
for size 3	32 – 42 cm
Price / pair	
☐ for size 1 (48 01 007)	381 €
☐ for size 2 (48 02 007)	381 €
☐ for size 3 (48 03 007)	381 €

Buttocks support pad



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price / piece	
☐ Size 1 (48 01 003)	incl.
☐ Size 2 (48 02 003)	incl.
☐ Size 3 (48 03 003)	incl.

Accessories

Back support pad – mounted at the back of the pelvis frame, height and angle adjustable

	Width	Height
for size 1	20 cm	23 cm
for size 2	20 cm	23 cm
for size 3	25 cm	28 cm
Price / piece		
☐ for size 1 (48 01 010)		404 €
☐ for size 2 (48 02 010)		404 €
☐ for size 3 (48 03 010)		404 €



Headrest – height, depth and angle adjustable, only for use together with the back support pad

	Price / piece
☐ for size 1 (48 01 111)	311 €
☐ for size 2 (48 02 111)	311 €
☐ for size 3 (48 03 111)	311 €



Configuration options - knee pads

Knee support pads with patella recess
individual three-dimensional adjustment



	Inside width		Inside width
Size 1	8 cm	Size 1	8 cm
Size 2	10 cm	Size 2	10 cm
Size 3	12 cm	Size 3	12 cm

Flexible knee positioning (FKP)
Individually three-dimensionally adjustable, with length-adjustable contact straps that adapt to the shape of the knee



Knee pad support – 2 different options are available to accommodate the adjacent knee pads

Rigid
angle, depth and height adjustable

Ball head
angle-, depth-, height- and additionally tilt angle adjustable

(Rigid = Standard configuration unless otherwise specified)

Please select one of the adjacent knee pads.

Price / pair		Upcharge / pair		Upcharge / pair	
☐ Size 1 (48 01 027)	incl.	☐ Size 1 (48 01 028)	109 €	☐ Rigid (48 00 026)	incl.
☐ Size 2 (48 02 027)	incl.	☐ Size 2 (48 02 028)	109 €	☐ Ball head (48 01 026)	143 €
☐ Size 3 (48 03 027)	incl.	☐ Size 3 (48 03 028)	109 €		

Accessories

Therapy tray



Therapy tray – with adjustable armrests



Therapy tray – with adjustable armrests and bowl (incl. cover)



Price / piece		Upcharge / piece		Upcharge / piece	
☐ for size 1 (48 01 114)	incl.	☐ for size 1 (48 01 116)	208 €	☐ for size 1 (48 01 117)	449 €
☐ for size 2 (48 02 114)	incl.	☐ for size 2 (48 02 116)	208 €	☐ for size 2 (48 02 117)	449 €
☐ for size 3 (48 03 114)	incl.	☐ for size 3 (48 03 116)	208 €	☐ for size 3 (48 03 117)	449 €

Accessories

Pad for therapy tray without armrests covered with black Flexi cover, height: 1 cm 	Pad for therapy tray with armrests covered with black Flexi cover, height: 1 cm	Angular adjustment of the center column by gas spring from (-) 15° - 30° (standard configuration = with foot pedal, if nothing else stated) 
Price / piece	Price / piece	Price / piece
☐ for size 1 (48 01 018) 162 € ☐ for size 2 (48 02 018) 162 € ☐ for size 3 (48 03 018) 162 €	☐ for size 1 (48 01 019) 208 € ☐ for size 2 (48 02 019) 208 € ☐ for size 3 (48 03 019) 208 €	☐ by foot pedal on the base frame (48 00 025) incl. ☐ by manual release under the tray (48 01 025) incl.

Accessories

Footplate – height and depth adjustable	Individual footplates – height, width and angle adjustable (abduction of legs)	Heel edges with ankle straps – angle, width and depth adjustable, for mounting on the footplate(s) with magnetic locks						
								
		<table><tr><td></td><td>Front width</td></tr><tr><td>Size 1</td><td>8,5 – 11 cm</td></tr><tr><td>Size 2</td><td>10 – 12,5 cm</td></tr></table>		Front width	Size 1	8,5 – 11 cm	Size 2	10 – 12,5 cm
	Front width							
Size 1	8,5 – 11 cm							
Size 2	10 – 12,5 cm							
Price / piece	Upcharge / pair	Price / pair						
☐ for size 1 (48 01 120) incl.	☐ for size 1 (48 01 121) 392 €	☐ Size 1 (48 01 024) 182 €						
☐ for size 2 (48 02 120) incl.	☐ for size 2 (48 02 121) 392 €	☐ Size 2 (48 02 024) 182 €						
☐ for size 3 (48 03 120) incl.	☐ for size 3 (48 03 121) 392 €							



activity chair.

The Hi/Lo base frame.

Further information can be found from page 44.



Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:


The transfer and mobility device.

- Abductable frame
- Padded trunk support
- Electric height adjustment
- Push bar
- 2 leg straps
- Battery and charger
- 4 dual castors (2 with brakes)




Dimensions

			tram.
A	Min. total width		70 cm
B	Max. total width		102 cm
C	Max. trunk support width		56 cm
D	Min. trunk support height		110 cm
E	Max. trunk support height		147 cm
F	Total length		116 cm
G	Push bar height		110 cm
H	Floor to top of base leg	low base	11 cm
		standard base	17,5 cm
	Girth		56 – 152 cm
	Turning diameter		127 cm
	Min. user armpit height		77 cm
	Max. user height		193 cm
	Weight		32 kg
	Max. load		158 kg



Basic model			
tram.		tram. with integrated scale	
		Price / piece	Price / piece
☐ Low base	(47 00 200) (K322+K304)	4.745 €	☐ Low base (47 00 300) (K321+K304) 5.675 €
☐ Standard base	(47 00 000) (K302+K304)	4.745 €	☐ Standard base (47 00 100) (K301+K304) 5.675 €

Accessories			
Arm prompts – incl. bracket 		Arm platforms – incl. bracket 	
	Padded surface Adjustment range (elbow to fist)		Pelvic support 
Size 1	B 11 x T 17 cm	20 – 31 cm	
Size 2	B 13 x T 21 cm	26 – 40 cm	
Price / pair		Price / pair	Price / piece
☐ Size 1 (47 01 001) (K305)	644 €	☐ 47 00 007 (K397)	672 €
☐ Size 2 (47 02 001) (K315)	644 €		
			Range of circumference Size 1 20 – 40 cm Size 2 40 – 70 cm Size 3 50 – 90 cm
			☐ Size 1 (47 01 002) (K306) 322 € ☐ Size 2 (47 02 002) (K316) 346 € ☐ Size 3 (47 03 002) (K326) 370 €

Accessories															
Thigh straps	Seat belt	Hip positioner – adjustable	Pad – for hip positioner												
		 													
		<table><tr><td></td><td>Size 1</td><td>Size 2</td></tr><tr><td>Crotch height (a)</td><td>6 cm</td><td>9 cm</td></tr><tr><td>Width (b)</td><td>3 cm</td><td>4,5 cm</td></tr><tr><td>Length (c)</td><td>23 cm</td><td>28 cm</td></tr></table>		Size 1	Size 2	Crotch height (a)	6 cm	9 cm	Width (b)	3 cm	4,5 cm	Length (c)	23 cm	28 cm	
	Size 1	Size 2													
Crotch height (a)	6 cm	9 cm													
Width (b)	3 cm	4,5 cm													
Length (c)	23 cm	28 cm													
Price / pair	Price / piece	Price / piece													
☐ Narrow (13 cm) (47 01 003) (K308) incl.	☐ Narrow (13 cm) (47 00 003) (K318) 95 €	☐ Size 1 (47 01 007) (K667)	601 €												
☐ Wide (18 cm) (47 01 106) (K338) incl.	☐ Wide (18 cm) (47 00 106) (K328) 95 €	☐ Size 2 (47 02 007) (K668)	668 €												
		☐ for size 1 (47 01 008) (K541)	217 €												
		☐ for size 2 (47 02 008) (K542)	217 €												

Accessories			
Switch pole – allows either the caregiver or client to operate the electric lift  		Additional battery 	
		Swivel lock  	
Price / piece		Price	
☐ 47 00 008 (K666)	324 €	☐ 47 00 004 (K307)	273 €
		☐ 1 Piece (47 05 005) (K323) for low version 108 € ☐ 1 Piece (47 00 005) (K309) for tall version 108 € ☐ 1 Pair (47 06 005) (K324) for low version 216 € ☐ 1 Pair (47 01 005) (K319) for tall version 216 €	

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____


Commission:
The comfortable toilet and shower chair.

- Anatomically-shaped seat pan made of PE, colour steel blue
- Subframe made of stainless steel
- Seat angle adjustment via gas pressure spring:
130 N (size 0), 200 N (sizes 1–3) and 2x 130 N (size 4)
- Mechanical height adjustment
- 4 dual castors (100 mm) with locking devices

HMV-No.: 33.40.04.1014

Dimensions

		Size 0	Size 1	Size 2	Size 3	Size 4
A	Seat width (rear)	26 cm	31 cm	36 cm	41 cm	46 cm
B	Seat depth	31 cm	36 cm	41 cm	46 cm	48 cm
C	Seat shell width (front)	31 cm	36 cm	41 cm	46 cm	51 cm
D	Back height	49 cm	51 cm	56 cm	61 cm	65 cm
E	Tilt-in-space	(-) 5° - 30°	(-) 5° - 30°	(-) 5° - 30°	(-) 5° - 30°	(-) 5° - 30°
F	Can be positioned over a toilet	47 - 57 cm	45 - 55 cm	42 - 52 cm	42 - 52 cm	42 - 52 cm
G	Total width	56 cm	61 cm	63 cm	71,5 cm	74 cm
H	Total length	65 cm	65 cm	65 cm	70 cm	70 cm
I	Max. lower leg length	38 cm	38 cm	48 cm	48 cm	48 cm
J	Front edge of hygiene opening to back edge of subframe	47 cm	49 cm	51 cm	54 cm	57 cm
K	Front opening width	10 cm	10 cm	10 cm	10 cm	10 cm
L	Rear opening width	12 cm	12 cm	13,5 cm	13,5 cm	13,5 cm
M*	Rear opening width	-	-	18,5 cm	18,5 cm	18,5 cm
	Max. load	60 kg	80 kg	100 kg	100 kg	150 kg
	Weight	28 kg	30 kg	33 kg	35 kg	38 kg


*** only applies to sizes 2b, 3b and 4b**

Basic model				
boris. size 0	boris. size 1	boris. size 2	boris. size 3	boris. size 4
Price / piece	Price / piece	Price / piece	Price / piece	Price / piece
☐ 52 00 100 2.540 €	☐ 52 01 100 2.565 €	☐ Size 2a (52 02 100/A) 2.565 € ☐ Size 2b (52 02 200/A) 2.640 €	☐ Size 3a (52 03 100/A) 2.640 € ☐ Size 3b (52 03 200/A) 2.715 €	☐ Size 4a (52 04 100/A) 2.955 € ☐ Size 4b (52 04 200/A) 3.040 €

Accessories			
PU-Headrest – blue, adjustable in height and depth with stainless steel bracket 	PU insert pad – two-piece for opening form a 	PU insert pad – two-piece for opening form b 	Tray – with padded arm rest segment and holder 
Price / piece	Price / piece	Price / piece	Price / piece
☐ 52 00 001 259 €	☐ for size 0 (52 00 101) 367 € ☐ for size 1 (52 01 101) 367 € ☐ for size 2a (52 02 101) 367 € ☐ for size 3a (52 03 101) 367 € ☐ for size 4a (52 04 101) 367 €	☐ for size 2b (52 02 201) 411 € ☐ for size 3b (52 03 201) 411 € ☐ for size 4b (52 04 201) 411 €	☐ for size 0 (52 00 012) 417 € ☐ for size 1 (52 01 002) 417 € ☐ for size 2 (52 02 002) 417 € ☐ for size 3 (52 03 002) 417 € ☐ for size 4 (52 04 002) 417 €

Accessories					
Toilet bucket – incl. lid and bracket as assembly set		Abductor wedge – with integrated splash guard (Note: a subsequent assembly at size 4 has to take place in our factory)		Footplate – height, knee and foot angle adjustment, also folding	
					
			Width	Height from seat	Adjustable lower leg length
for size 0+1		7 cm	14 cm	for size 0+1	30 – 38 cm
for size 2-4		8 cm	17.5 cm	for size 2-4	35 – 48 cm
Price / set		Price / piece		Price / piece	
☐ for size 0+1+3 (52 00 102) 188 €		☐ for size 0+1 (52 01 103) 201 €		☐ for size 0+1 (52 01 004) 360 €	
☐ for size 2+4 (52 02 102) 188 €		☐ for size 2+4 (52 02 103) 201 €		☐ for size 2-4 (52 02 004) 360 €	
		☐ for size 3 (52 03 103) 201 €			
				PU-Calf plate – height adjustable (only for use in combination with the footplate)	
					
				Height	15 cm
				Width	38 cm
				Price / piece	
				☐ 52 01 009 227 €	

Accessories		Strap holder		
Positioning jacket – made of spacer fabric incl. strap holder for shoulder straps, colour: blue 	Pelvic harness – made of spacer fabric, colour: blue 	Strap holder – for shoulder straps when using other harnesses apart from the adjacent positioning jacket 		
Price / set	Price / piece	Art.-No.	Nutzung	max. SchulterWidth
☐ Size 1 (52 01 010) 278 € ☐ Size 2 (52 02 010) 278 € ☐ Size 3 (52 03 010) 278 €	☐ Size 1 (52 01 011) 117 € ☐ Size 2 (52 02 011) 117 € ☐ Size 3 (52 03 011) 117 €	☐ 52 01 005	with Headrest	30 cm
		☐ 52 01 008	without Headrest	30 cm
				Price / piece
				117 €

Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:
The universal toilet seat system.

- Usable as an attachment on the toilet, with subframe freestanding or as a shower chair
- Incl. pelvic harness and removable armrests
- Adjustable seat depth and back angle

HMV-No.: 33.40.01.4005
(use as a toilet attachment)

HMV-No.: 33.40.04.1008
(use as a toilet chair)

Basic model

	hts. size 1	hts. size 2	hts. size 3
	Price / piece	Price / piece	Price / piece
	☐ 57 01 000 (Z111) 1.070 €	☐ 57 02 000 (Z121) 1.070 €	☐ 57 03 000 (Z131) 1.165 €

Subframes
Mobile subframe hts. – incl. footplate and 75 mm castors (with locking device and directional lock), adjustable in height, with seat angle adjustment (-)5° – 25°

Mobile subframe hts. – incl. footplate and 75 mm castors (with locking device and directional lock), adjustable in height, with seat angle adjustment (-)15° – 15°

Fixed subframe hts.
incl. footplate, adjustable in height


	Price / piece		Price / piece		Price / piece
☐ for size 1 (57 01 025) (Z166)	1.430 €	☐ for size 1 (57 01 002) (Z166+Z169)	1.430 €	☐ for size 1 (57 01 026) (Z184)	1.080 €
☐ for size 2 (57 02 025) (Z167)	1.430 €	☐ for size 2 (57 02 002) (Z167+Z169)	1.430 €	☐ for size 2 (57 02 026) (Z185)	1.080 €
☐ for size 3 (57 03 025) (Z168)	1.430 €	☐ for size 3 (57 03 002) (Z168+Z169)	1.430 €	☐ for size 3 (57 03 026) (Z186)	1.080 €

Subframe

Transfer and care subframe hts. – made of stainless steel, powder-coated, incl. footplate and 100-mm double castors, height-adjustable. With the transfer and care subframe HTS, we provide a further subframe for the Rifton HTS shower and toilet chair in order to better fulfil the transfer and care requirements of the user and their assistant. With the help of a pedal on both sides of the subframe the height can be adjusted. The seat can simply be tucked forward and backwards.

	for size 1+2	for size 3
Width	67 cm	71 cm
Length	84 cm	87 cm
Tilt-in-space	(-)5° - 25°	(-)5° - 25°
Seat height	51 - 78 cm	51 - 78 cm
Lower leg length	27 - 40 cm	40 - 53 cm
Price / piece		
☐ for size 1+2 (57 01 030)	2.840 €	
☐ for size 3 (57 03 030)	3.070 €	



Subframe

Bath subframe hts. – not possible for size 3, only for use without footplate

Width (bottom / top)	41 / 49 cm
Length	58 cm
Height (front)	25 cm
Price / piece	
☐ for size 1+2 (57 01 003) (2138)	390 €



Accessories

Attachment bar – to mount on the toilet (standard configuration, if a subframe has not been selected)





☐ for size 1+2 (57 01 004) (2115)	185 €
☐ for size 3 (57 03 004) (2135)	185 €

Price / piece

Dimensions

			Size 1	Size 2	Size 3
A	Width	Mobile subframes	53 cm	53 cm	57 cm
		Fixed subframes	51 cm	51 cm	53 cm
		Transfer and care subframe	67 cm	67 cm	71 cm
B	Seat height	Mobile subframes	39 - 55 cm	39 - 55 cm	44 - 58 cm
		Fixed subframes	33 - 43 cm	36 - 44 cm	38 - 48 cm
		Transfer and care subframe	51 - 78 cm	51 - 78 cm	51 - 78 cm
C	Seat width		20* / 25* / 30 cm	20* / 25* / 30* / 36 cm	28* / 33* / 38* / 43 cm
D	Seat depth		20 / 23 / 25 / 28 cm	29 / 32 / 34 / 37 cm	37 / 40 / 42 / 44 cm
E	Armrest height		13,5 cm	17 cm	20 cm
F	Back width		29 cm	34 cm	41 cm
G	Lower leg length toilet seat		22 - 34 cm	27 - 39 cm	34 - 44 cm
H	Back height		42 cm	50 cm	62 cm
I	Max. position above toilet height	Mobile subframes	52 cm	52 cm	52,5 cm
		Fixed subframes	38 cm	38 cm	42 cm
J	Length	Mobile subframes	74 cm	74 cm	81 cm
		Fixed subframes	71 cm	71 cm	74 cm
		Transfer and care subframe	84 cm	84 cm	87 cm
K	Front edge of hygiene opening to back edge of subframe		48 cm	48 cm	48 cm
	Back angle		(-)6°, 6°, 18°	(-)6°, 6°, 18°	(-)6°, 6°, 18°
	Tilt-in-space		(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
			(-)15° - 15°	(-)15° - 15°	(-)15° - 15°
	Max. load		34 kg	68 kg	113 kg

* adjustable with optional thigh guide



Armlehnen

	Arm rests for size 1 – detachable (standard configuration for size 1, if nothing else stated)		Arm rests for size 1-3 – detachable (standard configuration for size 2, if nothing else stated)		Arm rests for size 2+3 – detachable (standard configuration for size 3, if nothing else stated)	
	Height	13.5 cm	Height	17 cm	Height	20 cm
	Price / pair		Price / pair		Price / pair	
	☐ 57 01 001 (Z101) incl.		☐ 57 02 001 (Z102) incl.		☐ 57 03 001 (Z103) incl.	

Polster

Back pad and seat pad (closed at back)		Back pad and seat pad (open at back)		Back pad and seat pad (completely closed)		Seat pad (completely closed)	
							
Hole size in seat pad		Hole size in seat pad					
for size 1	8 x 20 cm	for size 1	10 x 23 cm				
for size 2	9 x 25 cm	for size 2	11 x 30 cm				
for size 3	11 x 30 cm	for size 3	15 x 36 cm				
Price / set		Price / set		Price / set		Price / piece	
☐ for size 1 (57 01 013) (Z154) 371 €		☐ for size 1 (57 01 012) (Z153) 371 €		☐ for size 1 (57 01 028) (Z187) 445 €		☐ for size 1 (57 01 029) (Z191) 251 €	
☐ for size 2 (57 02 013) (Z164) 371 €		☐ for size 2 (57 02 012) (Z163) 371 €		☐ for size 2 (57 02 028) (Z188) 445 €		☐ for size 2 (57 02 029) (Z192) 251 €	
☐ for size 3 (57 03 013) (Z174) 371 €		☐ for size 3 (57 03 012) (Z173) 371 €		☐ for size 3 (57 03 028) (Z189) 445 €		☐ for size 3 (57 03 029) (Z193) 251 €	

Accessories

Lateral pads – width and height adjustable				Lateral pads flip away – width and height adjustable			
	Size 1	Size 2			Size 1	Size 2	
Height	10 cm	13 cm		Height	10 cm	13 cm	
Depth	13 cm	17 cm		Depth	13 cm	17 cm	
to suit	Size 1+2	Size 2+3		to suit	Size 1+2	Size 2+3	
Chest width for size 1	15 – 28 cm			Chest width for size 1	15 – 28 cm		
Chest width for size 2	19 – 32 cm			Chest width for size 2	19 – 32 cm		
Chest width for size 3	28 – 42 cm			Chest width for size 3	28 – 42 cm		
			Price / pair				Price / pair
☐ Size 1	(57 01 007) (Z171)		330 €	☐ Size 1	(57 01 032) (Z1058)		460 €
☐ Size 2	(57 02 007) (Z172)		330 €	☐ Size 2	(57 02 032) (Z1068)		460 €
☐ Size 1 with additional chest strap	(57 01 008)		390 €	☐ Size 1 with additional chest strap	(57 01 033) (Z1078)		520 €
☐ Size 2 with additional chest strap	(57 02 008)		390 €	☐ Size 2 with additional chest strap	(57 02 033) (Z1088)		520 €

Accessories

Thigh guides width adjustable						Tray / arm rest		Shaped headrest – height, angle and depth adjustable	
	Pad height	Pad width	Seat width range						
			Size 1	Size 2	Size 3				
57 00 010	7,5 cm	18 cm	20/25 cm	25/30 cm	33/38 cm				
57 01 010	7,5 cm	28 cm	–	25/30 cm	33/38 cm				
57 02 010	7,5 cm	28 cm	–	20/25 cm	28/33 cm				
Price / pair						Price / piece		Price / piece	
☐ for size 1-3 (57 00 010) (Z109)					177 €	☐ for size 1 (57 01 011) (Z145)		☐ 57 00 006 (Z104)	
☐ for size 2+3 (57 01 010) (Z149)					199 €	☐ for size 2 (57 02 011) (Z155)			
☐ for size 2+3 (57 02 010) (Z198)					220 €	☐ for size 3 (57 03 011) (Z165)			

Accessories

Standard deflector – can only be used with the open seat padding and in conjunction with an abduction wedge 	Deflector with splash guard – can only be used with the open seat padding 	Abductor wedge 	Toilet bucket 
Price / piece	Price / piece	Price / piece	Price / piece
☐ for size 1+2 (57 01 014) (Z118) 77 € ☐ for size 3 (57 03 014) (Z128) 77 €	☐ for size 1+2 (57 01 015) (Z119) 159 € ☐ for size 3 (57 03 015) (Z129) 159 €	☐ for size 1+2 (57 01 016) (Z124) 181 € ☐ for size 3 (57 03 016) (Z144) 181 €	☐ for size 1 (57 01 017) (Z117) 129 € ☐ for size 2 (57 02 017) (Z127) 129 € ☐ for size 3 (57 03 017) (Z137) 129 €

Accessories

Footplate – when used as a toilet attachment 	Calf guide – only for use with the footplate 	Ankle cuffs – only for use with the footplate 
Price / piece	Price / piece	Price / pair
☐ for size 1 (57 01 018) (Z146) 495 € ☐ for size 2 (57 02 018) (Z147) 495 € ☐ for size 3 (57 03 018) (Z148) 495 €	☐ for size 1+2 (57 01 019) (Z158) 103 € ☐ for size 3 (57 03 019) (Z159) 103 €	☐ for size 1+2 (57 01 020) (Z156) 118 € ☐ for size 3 (57 03 020) (Z157) 118 €

Accessories

Chest-shoulder strap 	Push handles – for mobile subframes 	Castors 100mm – for mobile subframes 
Price / piece	Price / pair	Upcharge / set
☐ for size 1 (57 01 009) (Z116) 260 € ☐ for size 2 (57 02 009) (Z126) 260 € ☐ for size 3 (57 03 009) (Z136) 260 €	☐ 57 00 030 (Z199) 223 €	☐ 57 00 031 (Z181) 148 €

Accessories

Travel base hts. – with bag, max. toilet height = 47 cm, not possible for size 3, only for use without footplate (only for temporary use when traveling) 
Price / piece
☐ for size 1+2 (57 01 022) (Z139) 731 €

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The bathing and shower cradle.

- UV-resistant plastic frame with removable cover
- Seat and back section with adjustable angle
- Back angle adjustable with practical single hand operation
- Two heights possible via adjustable feet

HMV-No.: 04.40.04.0007



Basic model

wave. size 1 - 3 Standard



	Price / piece
☐ Size 1 (58 01 000) (Z211)	900 €
☐ Size 2 (58 02 000) (Z221)	955 €
☐ Size 3 (58 03 000) (Z231)	1.020 €

wave. size 1 - 3 with angle adjustable calf rest



	Price / piece
☐ Size 1 (58 01 100) (Z211+Z216)	1000 €
☐ Size 2 (58 02 100) (Z221+Z226)	1.055 €
☐ Size 3 (58 03 100) (Z231+Z236)	1.125 €

Basic model incl. accessory pack

wave size 1 - 3 standard, incl. headrest, pelvis strap and leg straps





	Price / piece
☐ Size 1 (58 01 050) (Z211+Z214+Z213+Z215)	945 €
☐ Size 2 (58 02 050) (Z221+Z224+Z223+Z225)	1.104 €
☐ Size 3 (58 03 050) (Z231+Z234+Z233+Z235)	1.167 €

wave size 1 - 3 with angle adjustable calf rest, incl. headrest, pelvis strap and leg straps



	Price / piece
☐ Size 1 (58 01 150) (Z211+Z214+Z213+Z215+Z216)	1.039 €
☐ Size 2 (58 02 150) (Z221+Z224+Z223+Z225+Z226)	1.199 €
☐ Size 3 (58 03 150) (Z231+Z234+Z233+Z235+Z236)	1.263 €

Bezugfarben			
Blue		Pink	
Price / piece		Price / piece	
☐ 9100 300 (BLU)	incl.	☐ 9100 301 (PNK)	incl.

Dimensions					
		Size 1	Size 2	Size 3	
A	Seat depth	27 cm	33 cm	43 cm	
B	Seat width	38 cm	38 cm	44 cm	
C	Back height	58 cm	67 cm	81 cm	
D	Seat height	with plastic frame	5 cm / 18 cm	5 cm / 18 cm	
		with bath cradle subframe	32 - 47 cm	32 - 47 cm	
		with shower cradle subframe	70 cm	70 cm	
		with transfer subframe	53 - 75 cm	53 - 75 cm	
E	Width	46 cm	46 cm	53 cm	
	Length of bathing cradle, completely flat, without calf pad	89 cm	104 cm	128 cm	
	Length of bathing cradle, completely flat, with calf pad	108 cm	124 cm	155 cm	
	Length, completely flat	75 cm	75 cm	75 cm	
	Height of bath subframe	20 cm	20 cm	20 cm	
	Length of bath subframe	99 cm	99 cm	99 cm	
	Width of bath subframe	44 cm	44 cm	44 cm	
	Height of shower subframe	56 cm	56 cm	56 cm	
	Width of shower subframe	69 cm	69 cm	69 cm	
	Length of shower subframe	86 cm	86 cm	86 cm	
	Max. load	34 kg	68 kg	114 kg	
	Back and seat angle	0°, 15°, 30°, 45°, 60°, 75°	0°, 15°, 30°, 45°, 60°, 75°	0°, 15°, 30°, 45°, 60°, 75°	

Untergestelle			
Bath subframe wave. to increase the bathing height 	Shower cradle subframe wave. mobile with 4x100 mm dual castors with brakes 	Transfer subframe wave. – mobile with 4 castors (2 with brakes), pivots to the right (by means of conversion also pivotable to the left), foldable, incl. wall hook to hang it on the wall. (we recommend the use in combination with a bathing and shower cradle with calf rest) 	
Price / piece	Price / piece	Price / piece	
☐ 58 00 500 (Z217) 635 €	☐ 58 00 510 (Z218) 1.350 €	☐ 58 00 520 (Z227 + Z265) 1.565 €	

Accessories			
Headrest / torso support 	Chest strap / pelvis strap 	Chest strap with lateral positioning 	Leg straps 
Price / piece	Price / piece	Price / piece	Price / piece
☐ for size 1+2 (58 01 001) (Z214) 142 €	☐ for size 1 (58 01 002) (Z213) 91 €	☐ for size 1 (58 01 003) (Z212) 113 €	☐ for size 1 (58 01 004) (Z215) 91 €
☐ for size 3 (58 03 001) (Z234) 142 €	☐ for size 2 (58 02 002) (Z223) 91 €	☐ for size 2 (58 02 003) (Z222) 113 €	☐ for size 2 (58 02 004) (Z225) 91 €
	☐ for size 3 (58 03 002) (Z233) 91 €	☐ for size 3 (58 03 003) (Z232) 113 €	☐ for size 3 (58 03 004) (Z235) 91 €

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

The versatile bathing and shower cradle.


- Black, powder-coated frame with removable, washable cover
- Incl. back angle, knee angle and head angle adjustment
- Incl. chest strap

HMV-No.: 04.40.04.0006



Dimensions

		Size 1	Size 2	Size 3	Size 4
A	Seat width	30 cm	35 cm	35 cm	40 cm
B	Seat depth	30 cm	34 cm	38 cm	44 cm
C	Back height	45 cm	61 cm	71 cm	82 cm
D	Standing space length	69 cm	74 cm	88 cm	93 cm
E	Lower leg length	19 cm	20 cm	28 cm	28 cm
F	Total width	36 cm	41 cm	41 cm	46 cm
	Back angle	90° - 180°	90° - 180°	90° - 180°	90° - 180°
	Knee angle	0° - 90°	0° - 90°	0° - 90°	0° - 90°
	Headrest angle	0° - 90°	0° - 90°	0° - 90°	0° - 90°
	Seat height	0° / 5 cm	0° / 5 cm	0° / 5 cm	0° / 5 cm
	Max. load	35 kg	45 kg	55 kg	75 kg
	Weight	4 kg	4,1 kg	4,3 kg	4,8 kg
* in combination with the cover with pelvic mould depression					



Basic model

basti. size 1	basti. size 2	basti. size 3	basti. size 4 back angle adjustment via stainless steel ratchet joint
Price / piece	Price / piece	Price / piece	Price / piece
☐ 56 01 000 945 €	☐ 56 02 000 1.005 €	☐ 56 03 000 1.060 €	☐ 56 04 000 1.200 €

Accessories

Cover with hollow for pelvis and adjustable tube edge



Upcharge / piece

- ☐ for size 1 **(56 01 001)**
- ☐ for size 2 **(56 02 001)**
- ☐ for size 3 **(56 03 001)**
- ☐ for size 4 **(56 04 001)**

153 €
153 €
182 €
182 €

Head guide pad – adjustable with velcro




Price / piece

- ☐ for size 1 **(56 01 002)**
- ☐ for size 2 **(56 02 002)**
- ☐ for size 3 **(56 03 002)**
- ☐ for size 4 **(56 04 002)**

250 €
250 €
250 €
250 €

Accessories

Additional retainer strap – adjustable with velcro for the chest or pelvis



Price / piece

- ☐ for size 1 **(56 01 003)**
- ☐ for size 2+3 **(56 03 003)**
- ☐ for size 4 **(56 04 003)**

121 €
121 €
121 €

Shower cradle subframe basti. – mobile with 4 castors with brakes





Price / piece

- ☐ **56 00 200**

830 €

Dimensions - subframe

A	Length	92 cm
B	Width	57 cm
C	Seat height with subframe	61* / 66 cm
	* in combination with the cover with pelvic mould depression	



Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:


The lightweight and compact bathing seat.

- Aluminium frame with removable cover
- Back section with adjustable angle
- Incl. pelvic harness



Dimensions

		bibi.	bibi. with long back
A	Seat depth	30 cm	30 cm
B	Seat width	36 cm	36 cm
C	Back height	45 cm	55 cm
D	Total length	60 cm	60 cm
	Back angle	90° - 180°	90° - 180°
	Max. load	30 kg	30 kg
	Weight	2,3 kg	2,3 kg



Basic model

bibi.		bibi. with long back	
	Price / piece		Price / piece
☐ 5501000	480 €	☐ 5501001	510 €

Accessories

Additional pelvic harness



	Price / piece
☐ 5501002	130 €

Name embroidery.

Invoice address:

Dispatch address:

☐ See invoice address

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____

Commission:

We hereby offer a way of customizing our aids **smilla**, **sofie**, **madita**, **madita-fun**, and **madita-fun mini**. You can provide the respective back pad with the name of the user.

Simply assign the name embroidery to the product and tell us the desired name. Then you have to only send this form together with the order form of the particular product and you will receive the ordered product with the desired name on the back pad.

- The height of the embroidery is 2.5cm (2cm for **madita-fun mini**.) – width is therefore given by the length of the name
- The embroidery is placed on the right side of the back pad (in the middle at **madita-fun mini**.)
- For the embroidery we use a silver-white yarn
- You can see the font that is used for the embroidery in the accompanying picture

Note: an embroidery of our Flexi fabric is not possible



Configuration name embroidery

Indication of the name please
in BLOCK LETTERS
(max. 12 letters)

--	--	--	--	--	--	--	--	--	--	--	--	--

Product assignment

☐ **smilla.**
☐ **sofie.**

☐ **madita.**
☐ **madita-fun.**

☐ **madita-fun mini.**

142 **Fabric collection.**
Invoice address:
Dispatch address:
☐ See invoice address

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Orderer: _____

Date: _____

Signature: _____


Commission:
Our Fabric collection.


Article description	Art.-No.	No.	Material composition (Upper)	goods width	Quantity	Price / rm	Price from 10 rm	Price from 30 rm
Woven fabric								
☐ Leaf-Green	(9001 123)	123	100% Polyester / Trevira CS	1,35 m	_____ m	36 €	30 €	25 €
☐ Slate grey	(9001 124)	124	100% Polyester / Trevira CS	1,35 m	_____ m			
☐ Flamigo-Rosé	(9001 125)	125	100% Polyester	1,35 m	_____ m			
☐ Pazifik-Mint	(9001 126)	126	100% Polyester	1,35 m	_____ m			
Knitted fabric								
☐ Red	(9001 091)	91	100% Polyester	1,50 m	_____ m	24 €	20 €	16 €
☐ Dark Grey	(9001 092)	92	100% Polyester	1,50 m	_____ m			
☐ Black	(9001 105)	105	100% Polyester	1,50 m	_____ m			
☐ Light blue	(9001 106)	106	100% Polyester	1,50 m	_____ m			
☐ Uni Blue	(9001 088)	88*	100% Polyester	1,50 m	_____ m			
☐ Uni Apple green	(9001 089)	89*	100% Polyester	1,50 m	_____ m			
☐ Uni Dark green	(9001 093)	93*	100% Polyester	1,50 m	_____ m			
☐ Deep Sea Blue	(9001 111)	111	100% Polyester / Trevira CS	1,55 m	_____ m	36 €	30 €	25 €
☐ Petrol	(9001 112)	112	100% Polyester / Trevira CS	1,55 m	_____ m			
☐ Pink-Telemagenta	(9001 114)	114	100% Polyester / Trevira CS	1,55 m	_____ m			
☐ Pink 2	(9001 115)	115	100% Polyester / Trevira CS	1,55 m	_____ m			
☐ Blue, with gray boxes	(9001 097)	97*	100% Polyester / Trevira CS	1,55 m	_____ m			
☐ Deep Orange	(9001 107)	107*	100% Polyester / Trevira CS	1,55 m	_____ m			
☐ Azur	(9001 110)	110*	100% Polyester / Trevira CS	1,55 m	_____ m			
☐ Rosa	(9001 113)	113*	100% Polyester / Trevira CS	1,55 m	_____ m			
☐ Jeans	(9001 116)	116*	100% Polyester / Trevira CS	1,55 m	_____ m			
Fabric Flexi, backed								
☐ Black	(9001 043)	43	56% Polyamid; 44% Polyurethan	1,40 m	_____ m	46 €	38 €	31 €
☐ Dark Blue	(9001 044)	44	63% Polyamid; 37% Polyurethan	1,40 m	_____ m			
☐ Light grey	(9001 050)	50	56% Polyamid; 44% Polyurethan	1,40 m	_____ m			
☐ Red	(9001 045)	45*	63% Polyamid; 37% Polyurethan	1,40 m	_____ m			
Fabric Flexi, unbacked								
☐ Red	(9001 042)	42*	63% Polyamid; 37% Polyurethan	1,40 m	_____ m	41 €	33 €	28 €
Mesh fabric								
☐ Blue	(9001 046)	46**	100% Polyester	1,50 m	_____ m	46 €	38 €	31 €
Spacer fabric								
☐ Black, 3mm	(9001 047)	47	100% Polyester	1,55 m	_____ m	36 €	30 €	25 €
☐ Blue, 3mm	(9001 048)	48**	100% Polyester	1,55 m	_____ m			
Spacer fabric (with closed surface)								
☐ Black-Black	(9001 101)	101	89% Polyester; 11% Elasthan	1,30 m	_____ m	36 €	30 €	25 €
☐ Black-Grey	(9001 109)	109	89% Polyester; 11% Elasthan	1,30 m	_____ m			
☐ Black-Blue	(9001 102)	102*	89% Polyester; 11% Elasthan	1,30 m	_____ m			

Upgrade your products with our fabrics, as they combine excellent quality with superb design. We provide you with a collection of the latest designs, different materials and the following convincing benefits:

- Simple and easy processing
- **High abrasion resistance and durability**
- Flame retardant in accordance with DIN EN 1021-1+2
- Skin-compatible in accordance with DIN EN ISO 10993-5
- Fulfills the requirements laid down in the MPG (Medical Products Directive) for your seat shells and special constructions

Woven fabric.


9001123 (No. 123)

9001124 (No. 124)

9001125 (No. 125)

9001126 (No. 126)

Knitted fabric.


9001091 (No. 91)

9001092 (No. 92)

9001105 (No. 105)

9001106 (No. 106)

9001111 (No. 111)*

9001112 (No. 112)*

9001114 (No. 114)*

9001115 (No. 115)*



Fabric Flexi, backed.


9001043 (No. 43)


9001044 (No. 44)

9001050 (No. 50)

Mesh fabric. / Spacer fabric.


9001046 (No. 46)**

9001047 (No. 47)

9001048 (No. 48)**

9001101 (No. 101)

9001109 (No. 109)

Discontinued models.


9001088 (No. 88)*

9001089 (No. 89)*

9001093 (No. 93)*

9001097 (No. 97)*

9001107 (No. 107)*

9001110 (No. 110)*

9001113 (No. 113)*

9001116 (No. 116)*

9001045 (No. 45)*

9001102 (No. 102)*

Note:

We reserve the right to carry out design, colour and material modifications according to the latest technology and within the scope of standard market production development. The images of the products and colours merely serve as orientation aids. Due to printing processes, deviations from the original cannot be excluded. All data without guarantee.

*Discontinued model

**This product is not DIN EN 1021-1+2 compliant.

Cleaning + disinfection.

Introduction to using our products.

Please note that **people with special needs, reading difficulties and impaired vision** must ask a parent, supervisor or technical advisor to read out or explain this technical documentation so that safe usage of our products can be guaranteed.

We would like to point out that only **spare parts and accessories** from our company may be used and installed so as not to endanger the safety of the users. In individual cases, the combination of our products with those of other manufacturers may require clarification with your specialist retailer or with us as the manufacturer.

Settings and adjustments to the product or accessories may only be made by people who have been given the necessary instructions by a medical product advisor. Please ensure that none of the user's extremities are in the respective area when making adjustments of any kind to minimise the risk of injury

The respective product should be reduced to its smallest dimensions **for transportation**. For example, reduce the seat, the armrests, the back etc. to the smallest dimensions. Remove the parts of the product that can be easily removed, e.g. the wheels of the **todd.** standing frame, the handlebars of the **momo.** etc. Fold up the foot support on products with this feature. Please see the User Manual belonging to the respective product for more detailed information.

Instructions for use on the respective products are available on www.schuchmann.de/en/media-library.

If applicable, the **product-specific cleaning guidelines** in the respective instructions for use are to be taken into consideration.

Cleaning guidelines for the different materials.

The care and cleaning of our products should be carried out at regular intervals to ensure that they remain functional over their useful life. The choice of suitable cleaning agents is crucial.

Wood: Wooden parts should be cleaned with a damp cloth or similar. When using detergents, please consult the manufacturer's instructions regarding detergent concentration. Caustic household detergents or corrosive substances are not suitable. Never leave natural wood in a damp or wet condition over long periods, as it will take on water and can change in shape (swelling). Therefore, only wipe and clean damp (i.e. with wrung-out cloths). Test a small, non-visible area for compatibility.

Powder-coated parts: Aluminium parts should be cleaned with a sponge or damp cloth using normal antibacterial household detergents.

Stainless steel: Stainless steel parts should be cleaned using a sponge or damp cloth with a mild soap. Bleach must not be used.

Covers / Fabric: Covers and fabric should be washed at 40° C in the washing machine using normal detergent. Do not tumble dry.

Plastics: Plastics should be cleaned using a sponge or damp cloth with mild soap. Bleach must not be used.

PU (polyurethane): PU parts should only be cleaned using a sponge or damp cloth with a mild soap. Bleach must not be used.

Should you have any
queries, you are welcome
to contact us directly!



Disinfection measures by material type.

Various products can be used for **surface disinfection** of **metal**, **wood** and **plastic parts**. The disinfectants used should have a wide range of efficacy (bactericidal, levurocidal and virucidal) and be approved for medical devices.

Suitable disinfectants can be found in the KRINKO/BfArM recommendation.

Liquid disinfectants are available as a ready-to-use solution which is sprayed on and evenly applied with a soft cloth. Alternatively, cloths soaked with disinfectant can be used, with to which the products are wiped all over. In both cases, a complete wetting be respected.

The exposure times may vary and can be found in the manufacturer's instructions for the agents used. However, products based on alcohol should not be used.

Disinfection systems achieve very good results. There the contaminated devices are completely and disinfected by cold fogging. Methods of disinfection by icing are not permitted.

If necessary, the **product-specific disinfection instructions** in the respective instructions for use are to be taken into consideration.

Hygiene measures for all our products.

1. Who should carry out the hygiene measures?

Users or accompanying persons themselves, service staff from specialist dealers and other authorised institutions.

2. How regularly?

- Cleaning: Regularly and as required (in case of soiling)
- Cleaning & disinfection: Regularly and as required (in case of possible contamination and before maintenance and servicing)
- Prior to cleaning and disinfection: Damaged sub-frames or components must first be removed or repaired.
- Processing by trained personnel (intensive cleaning and disinfection): Before passing on to other users.

3. With what should the work be carried out?

We recommend wearing protective gloves (e.g. made of nitrile or butyl rubber)! Avoid the penetration of liquids into Bowden cables, frame tubes or similar. Remove visible soiling with a damp sponge cloth and, if necessary, a cleaning agent. Do not use scouring agents or aids!

4. How should the work be conducted exactly?

For cleaning, all accessible surfaces are cleaned with a damp sponge cloth and, if necessary, a suitable cleaning medium. Then dry the surfaces with a cotton or disposable cloth. Hard-to-access places can also be dried using a hairdryer on its lowest setting. The use of high-pressure cleaners is not permitted, as their use can lead to damage to the products.

For disinfection the disinfectant should be used according to the manufacturer's instructions.

Ensure **complete wetting**. If necessary, remove excess preparation quantities with unstained disposable tissue. It is particularly important to observe the exposure times of the manufacturers!





schuchmann.de